Thioglycollate Medium w Vitamin K1 & Hemin

Liquid medium for the isolation and cultivation of fastidious or slow growing anaerobic organisms.

DESCRIPTION

Thioglycollate Medium enriched with Vitamin K1 and Hemin is a general purpose medium used in qualitative procedures for the cultivation of fastidious as well as nonfastidious microorganisms, including aerobic and anaerobic bacteria, from a variety of clinical and nonclinical materials.

TYPICAL FORMULA	(g/l)
Pancreatic Digest of Casein	17.0
Papaic Digest of Soybean Meal	3.0
Dextrose	6.0
Sodium Chloride	2.5
Sodium Thioglycollate	0.5
Sodium Sulfite	0.1
L-Cystine	0.25
Vitamin K1	0.001
Hemin	0.005
Agar	0.7

Final pH 7.0 ± 0.2 at 25°C

METHOD PRINCIPLE

Pancreatic digest of casein and papaic digest of soybean meal provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Dextrose is a source of energy. Sodium chloride maintains the osmotic balance of the medium. Sodium thioglycollate, sodium sulfite and L-cystine are included to reduce the redox potential of the medium and create an anaerobic atmosphere. Vitamin K1 is a growth requirement for some strains of *Prevotella melaninogenica* and is reported to enhance the growth of some strains of *Bacteroides* species and Gram-positive nonsporeformers. Hemin is the source of the X factor which stimulates the growth of many microorganisms. The small amount of agar assists in the maintenance of a low redox potential by stabilizing the medium against convection currents, thereby maintaining anaerobiosis in the lower depths of the medium.

PREPARATION

Medium in tubes/bottles For anaerobic incubation, liquid media may be reduced immediately prior to use by boiling in water bath, with caps loosened, and cooling with tightened caps, to room temperature before inoculation.

TEST PROCEDURE

Inoculate the specimen as soon as it arrives in the laboratory. Specimens known or suspected to contain obligate anaerobes should be inoculated near the bottom of the tube.

Incubate tubes with tight caps aerobically at 35 ± 2°C or other appropriate temperature. Growth of strictly aerobic bacteria can be improved by slightly loosening the cap. Broth cultures should be held at least 1 week before discarding as negative.

INTERPRETING RESULTS

Turbidity of the medium compared to uninoculated control indicates microbial growth.

Strict aerobes tend to grow in a thin layer at the surface of the broth; obligate anaerobes will grow only in that portion of the broth below the upper oxidized layer.

APPEARANCE

Slightly opalescent, light amber.

STORAGE

Store bottles and tubes at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Medium in bottles: 2 years.

Medium in tubes: 1 year.

QUALITY CONTROL

The medium is inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity: ≤ 100 CFU.

Incubation conditions: $35 \pm 2^\circ\text{C}$ for 48 h in an aerobic atmosphere with the caps screwed down tightly (except for *S. aureus* that is incubated with loosened cap).

QC Table.

Microorganism		Growth
<i>Peptostreptococcus anaerobius</i>	ATCC® 27337	Good
<i>Clostridium sporogenes</i>	ATCC® 19404	Good
<i>Clostridium perfringens</i>	ATCC® 13124	Good
<i>Staphylococcus aureus</i>	ATCC® 6538	Good

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. Dowell, V.R., Jr., and T.M. Hawkins (1987) Laboratory methods in anaerobic bacteriology. CDC laboratory manual. HHS Publication No. (CDC) 7-8272. Centers for Disease Control, Atlanta.
2. Dowell, V.R., Jr., G.L. Lombard, F.S. Thompson, and A.Y. Arm eld (1977) Media for isolation, characterization, and identi cation of obligately anaerobic bacteria. CDC laboratory manual. Center for Disease Control, Atlanta.
3. Finegold, S.M., V.L. Sutter, H.R. Attebery, and J.E. Rosenblatt (1974) Isolation of anaerobic bacteria, p. 365-375. In E.H. Lennette, G.H. Spaulding, and J.P. Truant (ed.), Manual of clinical microbiology, 2nd ed. American Society for Microbiology, Washington, D.C.
4. Gibbons, R.J. and J.B. MacDonald (1960) Hemin and vitamin K compounds as required factors for the cultivation of certain strains of *Bacteroides melaninogenicus*. J. Bacteriol. 80:164-170.
5. Brewer J.H. (1940) Clear liquid medium for the aerobic "cultivation" of anaerobes. J. Am. Med. Assoc. 115:598-600.

PRESENTATION

		Contents	Ref.
Thioglycollate Medium w Vitamin K1 & Hemin	Tubes	10 x 10 ml tubes	20171
Thioglycollate Medium w Vitamin K1 & Hemin	Bottles	6 x 100 ml bottles	442500

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care	 Keep away from sunlight
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse	

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BIBLIOGRAFIA

- Dowell, V.R., Jr., and T.M. Hawkins (1987) Laboratory methods in anaerobic bacteriology. CDC laboratory manual. HHS Publication No. (CDC) 7-8272. Centers for Disease Control, Atlanta.
- Dowell, V.R., Jr., G.L. Lombard, F.S. Thompson, and A.Y. Arm eld (1977) Media for isolation, characterization, and identi cation of obligately anaerobic bacteria. CDC laboratory manual. Center for Disease Control, Atlanta.
- Finegold, S.M., V.L. Sutter, H.R. Attebery, and J.E. Rosenblatt (1974) Isolation of anaerobic bacteria, p. 365-375. In E.H. Lennette, G.H. Spaulding, and J.P. Truant (ed.), Manual of clinical microbiology, 2nd ed. American Society for Microbiology, Washington, D.C.
- Gibbons, R.J. and J.B. MacDonald (1960) Hemin and vitamin K compounds as required factors for the cultivation of certain strains of *Bacteroides melaninogenicus*. J. Bacteriol. 80:164-170.
- Brewer J.H. (1940) Clear liquid medium for the aerobic "cultivation" of anaerobes. J. Am. Med. Assoc. 115:598-600.

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		Contenuto	Ref.
Thioglycollate Medium w Vitamin K1 & Hemin	Provette	Provette 10 x 10 ml	20171
Thioglycollate Medium w Vitamin K1 & Hemin	Flaconi	Flaconi 6 x 100 ml	442500

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