## Nutrient Broth

Liquid medium for the cultivation of nonfastidious microorganisms.

### DESCRIPTION

Nutrient Broth is a liquid medium used for the cultivation of a wide variety of organisms from clinical specimens and other materials.

This medium can be enriched with other ingredients such as blood, serum, sugars, etc., for special purposes.

### TYPICAL FORMULA

	(g/l)
Beef Extract	1.0
Peptone	5.0
Yeast Extract	2.0
Sodium Chloride	5.0
Final pH 6.8 ± 0.2 at 25°C	

### METHOD PRINCIPLE

Beef extract and peptone provide amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Yeast extract is a source of vitamins, particularly of B-group. Sodium chloride maintains the osmotic balance of the medium.

### PREPARATION

Dehydrated medium Suspend 13 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Sterilize in autoclave at 121°C for 15 minutes.

### TEST PROCEDURE

Inoculate broth with test sample. Incubate at 35 ± 2°C for 18-24 hours or longer if necessary.

### INTERPRETING RESULTS

Turbidity indicates microbial growth.

### APPEARANCE

Dehydrated medium: free-flowing, homogeneous, white to light beige.

Prepared medium: clear to slightly opalescent, light amber.

### STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and tubes at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

### SHELF LIFE

Dehydrated medium: 4 years.

Medium in tubes/bottles: 2 years.

**QUALITY CONTROL**

The medium is inoculated with the microbial strains indicated in the QC table.

Inoculum for productivity:  $\leq 100$  CFU

Incubation conditions: aerobically at  $35 \pm 2^\circ\text{C}$  for 18-24 hours.

**QC Table.**

Microorganism		Growth
<i>Escherichia coli</i>	ATCC® 25922	Good
<i>Staphylococcus aureus</i>	ATCC® 25923	Good

**WARNING AND PRECAUTIONS**

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

**DISPOSAL OF WASTE**

Disposal of waste must be carried out according to national and local regulations in force.

**BIBLIOGRAPHY**

1. Association of Official Analytical Chemists (1995) Official methods of analysis of AOAC International, 16<sup>th</sup> ed.
2. Marshall, R.T. (ed.) (1993) Standard methods for the microbiological examination of dairy products, 16<sup>th</sup> ed.
3. American Public Health Association (1923) Standard methods of water analysis, 5<sup>th</sup> ed.

PRESENTATION		Contents	Ref.
Nutrient Broth	Tubes	20 x 10 ml tubes	24103
Nutrient Broth	Tubes	50 x 5 ml tubes	27503
Nutrient Broth	Bottles	6 x 100 ml bottles	402000
Nutrient Broth	Bottles	6 x 500 ml bottles	470050
Nutrient Broth	Dehydrated medium	500 g of powder	610037
Nutrient Broth	Dehydrated medium	100 g of powder	620037
Nutrient Broth	Dehydrated medium	5 kg of powder	6100375

**TABLE OF SYMBOLS**

<b>LOT</b> Batch code	<b>IVD</b> <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
<b>REF</b> Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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