VTM

Transport medium for Viruses, *Chlamydia*, *Mycoplasma* and *Ureaplasma*.

Instructions For Use
ENGLISH

DESCRIPTION

VTM is a liquid medium used for transport of clinical specimens from the collection site to the testing laboratory.

This medium allows survival of even fragile organisms for relatively long periods of time at room temperature.

Some of the viruses and bacteria that can be maintained in VTM are the following:

- SARS-CoV-2 (the virus responsible for COVID-19)
- Herpes Simplex virus
- Varicella-Zoster virus
- Influenza Type A
- Respiratory Syncytial virus (RSV)
- Cytomegalovirus (CMV)
- *Chlamydia trachomatis*
- *Mycoplasma hominis*
- *Ureaplasma urealyticum*
- *Mycoplasma pneumoniae*

COMPOSITION

Sucrose

Hanks Balanced Salt Solution (HBSS)

Bovine Serum Albumin (BSA)

Buffered Solution

Gelatin

Amino acids

Antimicrobial Agents

Phenol Red

Final pH 7.3 ± 0.2 at 25°C

VTM does not contain guanidium thiocyanate.

METHOD PRINCIPLE

Sucrose acts as cryoprotectant ensuring organism viability during freeze-thaw. HBSS provides essential inorganic ions and along with buffered solution maintains pH and osmotic balance. BSA, gelatin and amino acids help to stabilise virus particles and to sustain viability for bacteria. Antibiotics and antimycotics are included to inhibit growth of bacteria and fungi. Phenol red is the pH indicator ensuring medium integrity at the time of specimen collection.

DIRECTIONS FOR USE

1. Collect swab specimens according to standard technique using a swab with a synthetic tip and a plastic shaft.
2. After collection, immediately place the swab specimen into a VTM tube.
3. Identify the tube containing the specimen.
4. Send promptly to the laboratory.

Specimens may be dispatched at ambient (10-25°C) or refrigerated temperature to arrive at the laboratory for processing within 48 hours.

If a delay in testing or shipping is expected, specimens should be frozen at -70°C or below (repeated freezing and thawing of specimens should be avoided as may reduce the recovery of viable organisms).

Notes:

- VTM is intended for transport of specimens only and must not be taken internally.
- VTM must not be used for premoistening or prewetting the swab prior to collecting the sample or for rinsing or irrigating the sampling sites.
- SWABS are NOT PROVIDED with this device.

- When choosing the collection swab, it is important to know that certain swab components may reduce recovery of some microorganisms or interfere with molecular detection methods. For example, calcium alginate swabs and wooden shaft swabs should not be used.
- All specimens collected for laboratory investigations should be regarded as potentially infectious.
- Condition, timing, and site of sampling as well as the type and volume of specimen are significant variables in obtaining reliable culture results. Refer to appropriate guidelines for detailed explanation of test procedure and sample collection.

APPEARANCE

Clear, light orange-red.

STORAGE

Store at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

The medium is tested using appropriate quality control strains and specifications as outlined in the Certificate of Analysis (CoA) available for each lot on Liofilchem's website.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

PRESENTATION	Format	Packaging	Ref.
VTM	16 x 100 mm Tube with internal shaped conical bottom and screw cap	100 x 3mL tubes	26490

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

This document is available from the online Support Center:
liofilchem.com/ifu-sds

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1. World Health Organization (2020). <https://www.who.int/emergencies/diseases/novel-coronavirus-2019/technical-guidance>
2. Center for Disease Control and Prevention (2020). Interim Guidelines for Collecting, Handling, and Testing Clinical Specimens from Persons for Coronavirus Disease 2019 (COVID-19). <https://www.cdc.gov/coronavirus/2019-ncov/lab/guidelines-clinical-specimens.html>
3. Public Health England (2017). Chlamydia trachomatis infection – testing by Nucleic Acid Amplification Tests (NAATs). UK Standards for Microbiology Investigations. V 37 Issue 4. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
4. Clinical and Laboratory Standards Institute (2014). Quality Control of Microbiological Transport Systems; Approved Standard - Second Edition. CLSI document M40-A2.
5. Leland, D.S. (1992). Concepts of clinical diagnostic virology, p. 3-43, In E.H. Lennette (ed.), Laboratory Diagnosis of Viral Infections, Second Edition. Marcel Dekker, Inc., New York.
6. Johnson, F.B. (1990). Transport of Viral Specimens. Clinical Microbiology Reviews 3(2):120-131.

TABLE OF SYMBOLS / TABELLA DEI SIMBOLI / TABLA DE SÍMBOLOS / TABELA DE SÍMBOLOS / SYMBOL TAFEL / TABLE DES SYMBOLES / TABEL MET SYMBOLEN / TABULKA SYMBOLŮ

LOT	Batch code / Codice del lotto / Número de lote / Código do lote / Batch-Code / Code du lot / Batchcode / Číslo šarže
REF	Catalogue number / Numero di catalogo / Referencia de catálogo / Número de catálogo / Katalognummer / Numéro de catalogue / Catalogus nummer / Katalogové číslo
IVD	<i>In Vitro</i> Diagnostic Medical Device / Dispositivo Medico Diagnostico <i>in Vitro</i> / IVD Dispositivo médico de diagnóstico <i>in Vitro</i> / IVD Dispositivo Médico de Diagnóstico <i>in Vitro</i> / <i>In-Vitro</i> -Diagnostik Medizinprodukt / Dispositif médical de Diagnostic <i>In Vitro</i> / <i>In Vitro</i> Diagnostisch medisch hulpmiddel / <i>In Vitro</i> Diagnostisch medisch hulpmiddel / <i>In vitro</i> Diagnostický zdravotnický prostředek
	Manufacturer / Fabbrikante / Fabricante / Fabricante / Hersteller / Fabricant / Fabrikant / Výrobce
	Use by / Utilizzare entro / Fecha de caducidad / Utilizar até / Verwendung durch / Utiliser par / Gebruiken voor / Spotřebujte do
	Fragile, handle with care / Fragile, maneggiare con cura / Frágil, maneje con cuidado / Frágil, manusear com cuidado / Zerbrechlich! Mit Vorsicht behandeln / Křehké, zacházejte opatrně
	Temperature limitation / Limiti di temperatura / Límite de temperatura / Limitação de temperatura / Temperaturbegrenzung / Fragile manipuler avec soin / Breekbaar - voorzichtig behandelen / Teplotní rozmezí
	Contains sufficient for <n> tests / Contenuto sufficiente per <n> saggi / Contiene suficiente para <n> pruebas / Contém o suficiente para <n> testes / Enthält ausreichend für <n> Tests / Contient suffisamment pour <n> tests / Bevat voldoende voor <n> testen / Obsahuje dostatečné množství pro <n> testů
	Caution, consult Instruction For Use / Attenzione, consultare le istruzioni per l'uso / Atención, consulte las instrucciones de uso / Consulte as instruções de utilização / Vorsicht, lesen Sie die Gebrauchsanweisung / Attention, consultez le mode d'emploi / Let op, raadpleeg de gebruiksaanwijzing / Pozor, čtěte Návod k použití
	Do not reuse / Non riutilizzare / No reutilizar / Não reutilizar / Nicht wiederverwenden / Ne pas réutiliser / Niet hergebruiken / Nepoužívat opakovaně

**LIOFILCHEM® s.r.l.**Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.com