Eugon LT SUP

Liquid culture medium for detection of microorganisms in cosmetics.

Instructions For Use
ENGLISH

DESCRIPTION

Eugon LT SUP is a highly nutritious medium used for the detection and enumeration of microorganisms in cosmetic products.

This medium complies with the requirements of the following standards: ISO 16212, ISO 17516, ISO 18415, ISO 18416, ISO 21149, ISO 21150, ISO 22717 and ISO 22718.

Eugon LT SUP is formulated to allow a better dispersion of the test sample while neutralizing preservatives that may be present in cosmetics.

TYPICAL FORMULA*	(g/l)
Pancreatic Digest of Casein	15.0
Papaic Digest of Soybean Meal	5.0
L-Cystine	0.7
Sodium Chloride	4.0
Sodium Sulfite	0.2
Glucose	5.5
Egg Lecithin	1.0
Polysorbate 80	15.0
Sodium Lauryl Sulfate	1.56
Final pH 7.0 ± 0.2 at 25°C	

*Formula may be adjusted and/or supplemented as required to meet performance specifications;
Grams per litre of purified water.

METHOD PRINCIPLE

Pancreatic digest of casein and papaic digest of soybean meal provide nitrogen, carbon, minerals and vitamins for microbial growth. L-cystine is a sulfur-containing amino acid that along with sodium sulfite generates a low redox potential (Eh), permitting the growth of some anaerobic organisms. Sodium chloride maintains the osmotic balance of the medium. Glucose is incorporated as a source of energy. Egg lecithin and polysorbate 80 are neutralizers effective against many preservatives such as phenyl derivatives, aldehydes and quaternary ammonium salts. Sodium lauryl sulfate is a dispersing agent.

TEST PROCEDURE

Following the procedure described by ISO 21149

For the detection of microorganisms by enrichment

- Transfer the sample to Eugon LT SUP, to obtain a 1/10 dilution (1 g or 1 mL in 9 ml / 10 g or 10 mL in 90 mL). If necessary, perform additional decimal dilutions using the same diluent.
- Incubate the suspension at 32.5 ± 2.5°C for at least 20 hours. The user is responsible for choosing the appropriate incubation temperature for the intended use, according to the reference ISO standard.
- Transfer 0.1 to 0.5 ml of the incubated suspension on the surface of an agar medium specific for the target microorganisms.

NOTE Do not invert the inoculated plate (or wait for the absorption of the incubated suspension by the agar before inverting).

For the enumeration of microorganisms

- Dilute the sample in Eugon LT SUP, usually by 1/10.
- NOTE The specimen must not remain in the broth for more than 45 minutes before inoculation onto the solid culture medium.
- Perform successive dilutions in the broth, if necessary.
 - Subculture onto the non-selective agar medium intended for the enumeration of the target microorganisms, such as TSA, by pour plate, surface spread or membrane filtration method.

INTERPRETING RESULTS

Refer to the instructions for use of the agar medium used.

STORAGE

Store at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Medium in bottles: 2 years.

Medium in tubes: 1 year.

QUALITY CONTROL

Appearance: Light amber; limpid, may have a slight precipitate.

Expected cultural response:

Control strain		Inoculum	Incubation	Specification
<i>Aspergillus brasiliensis</i>	WDCM 00053 (ATCC 16404; NCPF 2275)	$\leq 100 \text{ CFU}$	24-72 h / 32.5 \pm 2.5°C	Good growth
<i>Bacillus subtilis</i>	WDCM 00003 (ATCC 6633; NCTC 10400)			
<i>Candida albicans</i>	WDCM 00054 (ATCC 10231; NCPF 3179)			
<i>Escherichia coli</i>	WDCM 00012 (ATCC 8739; NCTC 12923)			
<i>Pseudomonas aeruginosa</i>	WDCM 00026 (ATCC 9027; NCTC 12924)			
<i>Staphylococcus aureus</i>	WDCM 00032 (ATCC 6538; NCTC 10788)			
<i>Staphylococcus aureus</i> + AAP	WDCM 00032 (ATCC 6538; NCTC 10788)			
<i>Pseudomonas aeruginosa</i> + AAP	WDCM 00026 (ATCC 9027; NCTC 12924)		20 h / 32.5 \pm 2.5°C	Good growth ($P_R \geq 0.5$ on TSA)

AAP: Antimicrobial Activity Product.

A productivity ratio (P_R) of 0.5 is equivalent to a recovery rate of 50%.

Please refer to the actual batch related Certificate of Analysis (CoA).

WARNING AND PRECAUTIONS

For professional use only. Operators must be trained and have certain experience. Please read the instructions carefully before using this product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.

Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

The product is available in the configurations listed on the next page. There may be additional product ref. numbers as well. For an updated listing of available products, visit lioofilchem.com

Product	Format	Packaging	Ref.
Eugon LT SUP	Tube	100 x 9 ml	26492
	Plastic Tube	100 x 9 ml	26492P
	Bottle	6 x 90 ml (screw cap)	402750
	Bottle	6 x 500 ml (perforable cap)	463300
	Bottle	6 x 1 000 ml (perforable cap)	463310

This IFU document and the SDS are available from the online Support Center:

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liofilchem.com/ifu-sds

BIBLIOGRAPHY / BIBLIOGRAFIA

1. ISO 16212 Cosmetics – Microbiology – Enumeration of yeast and mould.
2. ISO 17516 Cosmetics – Microbiology – Microbiological limits.
3. ISO 18415 Cosmetics – Microbiology – Detection of specified and non-specified microorganisms.
4. ISO 18416 Cosmetics – Microbiology – Detection of *Candida albicans*.
5. ISO 21149 Cosmetics – Microbiology – Enumeration and detection of aerobic mesophilic bacteria.
6. ISO 21150 Cosmetics – Microbiology – Detection of *Escherichia coli*.
7. ISO 22717 Cosmetics – Microbiology – Detection of *Pseudomonas aeruginosa*.
8. ISO 22718 Cosmetics – Microbiology – Detection of *Staphylococcus aureus*.
9. Vera, H. D. (1947) The ability of peptones to support surface growth of lactobacilli. J. Bacteriol. 54:14.
10. Guisno, R., I.W. Gibby & M.J. Foter (1946) A neutralizing medium for evaluation of the germicidal potency of the quaternary ammonium salts. Amer J Pharm 118:320-323.

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