# A.F. Genital System

**ENGLISH**

## System for detection, count and susceptibility testing of pathogenic urogenital microorganisms.

### DESCRIPTION

A.F. Genital System is a 24-well system containing desiccated biochemical and antibiotic substrates for detection, presumptive identification and susceptibility testing of microorganisms from urogenital specimens. The system also provides a semi-quantitative assessment of the presence of urogenital mycoplasmas, i.e. *Mycoplasma hominis* and *Ureaplasma urealyticum*. The system is inoculated with the suspension of the clinical sample and incubated at  $36 \pm 1^\circ\text{C}$ . The test results are read after 24 hours of incubation and interpreted by assessing the change in color of the various wells and performing a microscopic examination.

**Note:** The species *Ureaplasma urealyticum* has been divided into two new species: *Ureaplasma parvum* and *Ureaplasma urealyticum*. They are considered together as *Ureaplasma* spp. when using this device and indicated as Uu.

### CONTENT OF THE KIT (kit for 20 tests)

- 20 A.F. Genital System
- Instructions Sheet also available at [www.liofilchem.com/ifu-sds](http://www.liofilchem.com/ifu-sds)
- 20 Vials of Physiological Solution (7.0 mL/vial)
- Test Results Form

### ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

- A.F. Genital System Reagent (ref. 80258)
- Mycoplasma Transport Broth (ref. 20158)
- Strepto B Latex Kit (ref. 96153)
- Sampling swabs
- Cover slides and slides
- Microscope

### CONFIGURATION

Well	<b>Count and identification of mycoplasmas / ureaplasmas</b>
1-Uu $10^3$	<i>Ureaplasma</i> spp. (titre = $10^3$ CFU/mL)
2-Uu $10^4$	<i>Ureaplasma</i> spp. (titre = $10^4$ CFU/mL)
3-Uu $\geq 10^5$	<i>Ureaplasma</i> spp. (titre $\geq 10^5$ CFU/mL)
4-Mh $10^4$	<i>Mycoplasma hominis</i> (titre = $10^4$ CFU/mL)
5-Mh $\geq 10^5$	<i>Mycoplasma hominis</i> (titre $\geq 10^5$ CFU/mL)
Well	<b>Detection of <i>T. vaginalis</i> and <i>Candida</i> spp.</b>
6-TR/YE	<i>Trichomonas vaginalis</i> / <i>Candida</i> spp.
Well	<b>Susceptibility test of the mycoplasmas / ureaplasmas</b>
7-TE	Tetracycline - 8 $\mu\text{g/mL}$
8-PEF	Pefloxacin - 16 $\mu\text{g/mL}$
9-OFX	Ofloxacin - 4 $\mu\text{g/mL}$
10-DO	Doxycycline - 8 $\mu\text{g/mL}$
11-E	Erythromycin - 16 $\mu\text{g/mL}$
12-CLA	Clarithromycin - 16 $\mu\text{g/mL}$
13-MN	Minocycline - 8 $\mu\text{g/mL}$
14-JOS	Josamycin - 8 $\mu\text{g/mL}$
15-CD	Clindamycin - 8 $\mu\text{g/mL}$
Well	<b>Detection and presumptive identification of other microorganisms</b>
16-ESC	<i>Escherichia coli</i>
17-PRO	<i>Proteus</i> spp. / <i>Providencia</i> spp.
18-PSE	<i>Pseudomonas</i> spp.
19-GAR	<i>Gardnerella vaginalis</i>
20-STF	<i>Staphylococcus aureus</i>
21-STR	<i>Enterococcus faecalis</i>
22-NES	<i>Neisseria gonorrhoeae</i>
23-STG	<i>Streptococcus agalactiae</i> (Group B)
24-CAN	<i>Candida</i> spp.

## PRINCIPLE OF THE METHOD

A.F. Genital System allows the detection, semi-quantitative count, presumptive identification and susceptibility test of *Mycoplasma hominis* and *Ureaplasma* spp., and the detection and presumptive identification of the microorganisms most frequently isolated from vaginal and urethral swabs and seminal fluid: *Trichomonas vaginalis*, *Escherichia coli*, *Proteus* spp./*Providencia* spp., *Pseudomonas* spp., *Gardnerella vaginalis*, *Staphylococcus aureus*, *Enterococcus faecalis*, *Neisseria gonorrhoeae*, *Streptococcus agalactiae* (Group B) and *Candida* spp.

- The semi-quantitative count of *Ureaplasma* spp. is indicated by a color change from yellow to red in wells **1-Uu 10<sup>3</sup>**, **2-Uu 10<sup>4</sup>** and **3-Uu  $\geq 10^5$** .
- The semi-quantitative count of *Mycoplasma hominis* is evidenced by a color change from yellow to red in wells **4-Mh 10<sup>4</sup>** and **5-Mh  $\geq 10^5$** .
- The presence of *Trichomonas vaginalis* and *Candida* spp. is determined by microscope (40x) examination of a drop of culture liquid taken from well **6-TR/YE** to assess the presence of mobile ciliated trophozoites for identification of *Trichomonas vaginalis* and/or the presence of chlamydozoites and hyphae for identification of *Candida* spp.
- The red color change of the medium in the wells containing an antibiotic (**7-TE**, **8-PEF**, **9-OFX**, **10-DO**, **11-E**, **12-CLA**, **13-MN**, **14-JOS**, **15-CD**) indicates the presence of bacterial growth and hence resistance to the antibiotic concentration being tested. The yellow color of the medium indicates the absence of bacterial growth and hence susceptibility to the antibiotic concentration being tested.
- The presence of *Escherichia coli* is detected by a color change from grey-red to blue in the well **16-ESC**. Confirm the presence of *Escherichia coli* by adding two drops of Kovacs Reagent into the well **17-PRO**; wait for the formation of a pink-red ring.
- The presence of *Proteus* spp./*Providencia* spp. is shown by a color change from yellow to brown-black in the well **17-PRO**.
- The presence of *Pseudomonas* spp. is shown by a color change from yellow to turbid green in the well **18-PSE**. Confirm by oxidase test (A.F. Genital System Reagent, ref. 80258).
- The presence of *Gardnerella vaginalis* is shown by a color change from red to yellow-orange in the well **19-GAR**. Confirm by isolation on selective media (*Gardnerella* Agar, ref. 11054) and biochemical tests (*Integral System Gardnerella*, ref. 71724).
- The presence of *Staphylococcus aureus* is shown by the appearance of a black ring on the bottom of well **20-STF**. Confirm by coagulase test (ref. 88030).
- The presence of *Enterococcus faecalis* is shown by a color change from yellow to black in the well **21-STR**.
- The presence of *Neisseria gonorrhoeae* is indicated by a positive oxidase test from the well **22-NES** (A.F. Genital System Reagent, ref. 80258). Confirm the result by isolation on selective media (*Thayer Martin Agar*, ref. 11040), and by performing biochemical and/or serological tests.
- The presence of *Streptococcus agalactiae* (Group B) is shown by a color change from yellow to green in the well **23-STG**. Confirm the presence of *Streptococcus agalactiae* with the following methods:
  - a) – Agglutination test *Strepto B Latex Kit* (ref. 96153)
  - b) – Cultural isolation with traditional methods
- The presence of *Candida* spp. is shown by a color change from green to yellow in the well **24-CAN**. Watch for the presence of chlamydozoites and hyphae at the microscope.

## COLLECTION AND STORAGE OF THE SAMPLES

Obtain the vaginal or urethral secretion with synthetic fibre swabs.

Obtain the seminal fluid according to the technique envisaged for the microbiological culture examination.

The samples must be sent for inoculation in the A.F. Genital System immediately after they are taken.

Do not store in the refrigerator for even limited periods of time, since low temperatures may damage the vitality of some particularly sensitive microorganisms (e.g. *Trichomonas vaginalis*, *Neisseria gonorrhoeae* etc.) and alter the final result.

## TEST PROCEDURE

### PREPARATION OF THE CLINICAL SAMPLE

Carefully open a vial of Physiological Solution\* using the ampule cracker provided with the kit.

#### 1. VAGINAL SWAB – URETHRAL SWAB

Immerse the swab (after obtaining the clinical material) in the vial of physiological solution and wait 5 minutes. Carefully squeeze the swab against the vial wall so that the clinical material is dispersed homogeneously in the Physiological Solution.

**Note:** Keep the swab, dipped in a nutritive broth, until the conclusion of the test.

#### 2. SEMINAL FLUID

Add 0.2 mL of sample to the vial of Physiological Solution\*, shake and wait 5 minutes before inoculation.

#### 3. URINE

Centrifuge about 10 mL of urinary sample. Collect 0.1 mL of urinary sediment and add to the vial of Physiological Solution contained in the kit; shake and wait 5 minutes before inoculating the system.

#### 4. CLINICAL MATERIAL FROM Mycoplasma Transport Broth (ref. 20158)

Add 1 mL Mycoplasma Transport Broth, containing the sample under examination, to the vial of Physiological Solution, shake gently and wait 5 minutes before inoculation.

**Note:** Conserve the specimen in the vial of Mycoplasma Transport Broth until the conclusion of the test.

\* Physiological Solution (g\L): Sodium chloride 9 g; Distilled water 1000 mL; pH 6.8 ± 0.2

### INOCULATION OF THE SYSTEM

1. Take a system from its wrapper and bring it to room temperature.
2. Write down the name of the patient, the date of the start of the examination, and the type of clinical material.
3. Transfer 0.2 mL of suspension of the clinical sample into each well of the system.
4. Overlay with 1 drop vaseline oil for microbiological use all the wells including 19-GAR, except for the wells 6, 16, 17, 18, 20, 21, 22 and 23.
5. Cover the system with the lid provided and incubate at 36 ± 1°C for 24 hours, and up to 48 hours.
6. Observe and confirm the presence of mycoplasmas within 48 hours incubation.

### INTERPRETATION OF THE RESULTS

#### Count and Identification of *mycoplasmas* / *ureaplasmas*

Watch for the color change in the wells **1-Uu 10<sup>3</sup>** to **5-Mh ≥ 10<sup>5</sup>** and interpret the results using Table 1.

#### Detection of *Trichomonas vaginalis* and *Candida* spp.

Take a drop of liquid from the well, deposit it on a glass slide, place a cover slip on top and examine at the microscope (40x) for the presence of *T. vaginalis* and *Candida* spp. Interpret as indicated in Table 1.

#### Antimicrobial Susceptibility Testing of *mycoplasmas* / *ureaplasmas*

Watch for the color change in the wells **7-TE** to **15-CD** and interpret the results using Table 1.

#### Identification of *Escherichia coli*

Watch for the color change in the well **16-ESC** and interpret the results using Table 1. To confirm the presence of *E. coli*, add two drops of Kovacs Reagent to the well and wait for the formation of a pink-red ring.

#### Identification of *Neisseria gonorrhoeae*

Take a drop of liquid from the well **22-NES**, deposit it on an Oxidase Test Stick and observe the immediate development (in about 10 seconds) of a blue color (positive oxidase test). Confirm the presence of *N. gonorrhoeae* through isolation on selective media: streak the vaginal or urethral swab, previously kept in a nutritive broth (see point 1 in the TEST PROCEDURE), directly on the surface of a selective plate (Thayer Martin Agar ref. 11040). Incubate the plate at 36 ± 1°C under microaerobic atmosphere, and perform biochemical and serological tests for final identification.

#### Identification of *Streptococcus agalactiae*

Watch for the color change in the well **23-STG** and interpret the results using Table 1. Confirm the presence of *S. agalactiae* with Strepto B Latex Kit Agglutination test (ref. 96153) and by cultural isolation with traditional methods.

#### Presumptive Identification of *other microorganisms*

Watch for the color change in the wells **17-PRO** to **24-CAN** ed and interpret the results using Table 1.

Note the results on the TEST RESULTS FORM (copy as many forms as necessary).

Table 1.

Well	COUNT AND IDENTIFICATION OF MYCOPLASMAS / UREAPLASMAS	Well colour	
		Positive reaction*	Negative reaction
1-Uu 10 <sup>3</sup>	<i>Ureaplasma</i> spp. (titre = 10 <sup>3</sup> CFU/mL) <sup>(1)</sup>	red	yellow
2-Uu 10 <sup>4</sup>	<i>Ureaplasma</i> spp. (titre = 10 <sup>4</sup> CFU/mL) <sup>(2)</sup>	red	yellow
3-Uu ≥ 10 <sup>5</sup>	<i>Ureaplasma</i> spp. (titre ≥ 10 <sup>5</sup> CFU/mL) <sup>(3)</sup>	red	yellow
4-Mh 10 <sup>4</sup>	<i>Mycoplasma hominis</i> (titre = 10 <sup>4</sup> CFU/mL)	red	yellow
5-Mh ≥ 10 <sup>5</sup>	<i>Mycoplasma hominis</i> (titre ≥ 10 <sup>5</sup> CFU/mL)	red	yellow

CFU: Colony Forming Units

\* An orangey colouring should be considered as a positive test.

Well	T. VAGINALIS AND CANDIDA spp.	Microscope examination (40x)
6-TR/YE	<i>Trichomonas vaginalis</i> / <i>Candida</i> spp.	<i>T. vaginalis</i> : mobile ciliated protozoites <i>Candida</i> spp.: chlamydospores and hyphae

Well	SUSCEPTIBILITY TESTING OF THE MYCOPLASMAS / UREAPLASMAS	Well colour	
		S	R**
7-TE	Tetracycline - 8 µg/mL	yellow	red
8-PEF	Pefloxacin - 16 µg/mL	yellow	red
9-OFX	Ofloxacin - 4 µg/mL	yellow	red
10-DO	Doxycycline - 8 µg/mL	yellow	red
11-E	Erythromycin - 16 µg/mL	yellow	red
12-CLA	Clarithromycin - 16 µg/mL	yellow	red
13-MN	Minocycline - 8 µg/mL	yellow	red
14-JOS	Josamycin - 8 µg/mL	yellow	red
15-CD	Clindamycin - 8 µg/mL	yellow	red

S = Sensitive

R = Resistant

\*\* An orangey colouring should be considered as presence of growth and hence resistance to the tested antibiotic concentration.

Well	PRESUMPTIVE IDENTIFICATION OF OTHER MICROORGANISMS	Well colour	
		Positive reaction	Negative reaction
16-ESC	<i>Escherichia coli</i>	blue	grey-red
17-PRO	<i>Proteus</i> spp. / <i>Providencia</i> spp.	brown-black	yellow
18-PSE	<i>Pseudomonas</i> spp.	turbid green	yellow-blue
19-GAR	<i>Gardnerella vaginalis</i>	yellow-orange	red
20-STF	<i>Staphylococcus aureus</i>	black ring	yellow
21-STR	<i>Enterococcus faecalis</i>	black	yellow
22-NES	<i>Neisseria gonorrhoeae</i>	blue <sup>(4)</sup>	colourless <sup>(4)</sup>
23-STG	<i>Streptococcus agalactiae</i> (Group B)	green	yellow
24-CAN	<i>Candida</i> spp.	turbid yellow	green

#### LEGEND

(1) = Equivalent to 5-20 colonies grown on Mycoplasma Agar A7

(2) = Equivalent to 20-50 colonies grown on Mycoplasma Agar A7

(3) = Over 50 colonies grown on Mycoplasma Agar A7

(4) = Oxidase Test Result

**QUALITY CONTROL**

Every batch of A.F. Genital System is subjected to quality control using the following reference microorganisms:

<i>Mycoplasma hominis</i>	ATCC 23114	<i>Gardnerella vaginalis</i>	ATCC 14018
<i>Ureaplasma urealiticum</i>	ATCC 27618	<i>Staphylococcus aureus</i>	ATCC 25923
<i>Candida albicans</i>	ATCC 10231	<i>Streptococcus faecalis</i>	ATCC 19433
<i>Escherichia coli</i>	ATCC 25922	<i>Neisseria gonorrhoeae</i>	ATCC 19424
<i>Proteus mirabilis</i>	ATCC 25933	<i>Streptococcus agalactiae</i>	ATCC 13813
<i>Pseudomonas aeruginosa</i>	ATCC 27853	<i>Trichomonas vaginalis</i>	ATCC 30001

**FACTORS THAT MAY INVALIDATE THE RESULTS**

Poor standardisation of the inoculum; clinical material unsuitable; use of expired systems and/or reagents; non compliance with temperatures and times of incubation.

**LIMITS AND WARNINGS**

- *Neisseria gonorrhoeae*, *Streptococcus agalactiae* and *Gardnerella vaginalis* require appropriate culture media, and additional biochemical and/or serological tests for complete identification.
- Some bacteria that are present in quantities of  $\geq 10^{6-7}$  CFU/mL and produce urease, may cause all the wells in the panel to change colour. The presence of these can be verified by reisolating on agar medium (e.g. chocolate agar) from the original broth (see TEST PROCEDURE).

**PRECAUTIONS**

The product A.F. Genital System is not classified as hazardous under current legislation, however see the safety data sheet for a correct use. A.F. Genital System is a disposable device to be used only for diagnostic use *in vitro*. It must be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

**STORAGE**

Store at 2-8°C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

**DISPOSAL OF USED MATERIAL**

After use, A.F. Genital System and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

**TABLE OF SYMBOLS**

<b>IVD</b> <i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
<b>REF</b> Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	<b>LOT</b> Batch code

Product	Packaging	Ref.
A.F. Genital System	20 tests	74156
A.F. Genital System	4 tests	79156

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## BIBLIOGRAPHY / BIBLIOGRAFIA

1. Horner P, Donders G, Cusini M, *et al.* Should we be testing for urogenital *Mycoplasma hominis*, *Ureaplasma parvum* and *Ureaplasma urealyticum* in men and women? - a position statement from the European STI Guidelines Editorial Board. *Eur Acad Dermatol Venereol.* 2018;32(11):1845-1851.
2. Capoccia R, Greub G, Baud D. *Ureaplasma urealyticum*, *Mycoplasma hominis* and adverse pregnancy outcomes. *Curr Opin Infect Dis.* 2013;26(3):231-40.
3. Waites KB, Bebear CM, Robertson JA, Talkington DF, and Kenny GE. Laboratory diagnosis of mycoplasmal infections. Cumulative techniques and procedures in clinical microbiology (Cumitech 34). American Society for Microbiology Press, Washington, D.C. 2001; 1-30.
4. Ullmann U, Schubert S, Krausse R. Comparative in-vitro activity of levofloxacin, other fluoroquinolones, doxycycline and erythromycin against *Ureaplasma urealyticum* and *Mycoplasma hominis*. *J Antimicrob Chemother.* 1999;43 Suppl C:33-6.
5. Taylor-Robinson D, Bébéar C. Antibiotic susceptibilities of mycoplasmas and treatment of mycoplasmal infections. *J Antimicrob Chemother* 1997;40(5):622-30.
6. Grattard F, Soleihac B, De Barbeyrac B *et al.* Epidemiologic and molecular investigations of genital mycoplasmas from women and neonates at delivery. *Pediatr Infect Dis J.* 1995;14(10):853-8.
7. Bébéar C, de Barbeyrac B, Dewilde A *et al.* Multicenter study of the in vitro sensitivity of genital mycoplasmas to antibiotics. *Pathol Biol.* 1993;41(4):289-93.
8. Kenny GE, Cartwright FD. Effect of pH, inoculum size, and incubation time on the susceptibility of *Ureaplasma urealyticum* to erythromycin in vitro. *Clin Infect Dis.* 1993;17 Suppl 1:S215-8.
9. Bauriaud R, Séror C, Lareng MB, Lefèvre JC. In vitro sensitivity to antibiotics of genital mycoplasmas isolated in Toulouse. Study of new molecules (macrolides and quinolones). *Pathol Biol.* 1992;40(5):479-82.
10. Renaudin H, Bébéar C. Comparative in vitro activity of azithromycin, clarithromycin, erythromycin and lomefloxacin against *Mycoplasma pneumoniae*, *Mycoplasma hominis* and *Ureaplasma urealyticum*. *Eur J Clin Microbiol Infect Dis.* 1990;9(11):838-41.
11. Bebear C, de Barbeyrac B, Bernet C, Renaudin H. Methods of studying mycoplasma infections. *Ann Biol Clin.* 1989;47(7): 415-20.

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