# ComASP™ Colistin

System for colistin susceptibility testing with the broth microdilution method.

## DESCRIPTION

Clinicians are re-thinking to colistin as therapeutic option for the treatment of severe infections caused by multidrug-resistant microorganisms, such as *Pseudomonas aeruginosa*, *Acinetobacter baumannii* and carbapenem-resistant Enterobacteriaceae. However, colistin-resistant strains are worldwide disseminated leading to the need to determine the value of minimum inhibitory concentration (MIC) before the drug is prescribed. Although several commercial tests based on different techniques have been developed, broth microdilution (BMD) is considered the best method for performing colistin susceptibility testing so far.

ComASP Colistin is a 4-test panel containing the dried up antibiotic in 7 two-fold dilutions (0.25 - 16 µg/ml). The system is used to perform the broth microdilution (BMD) method for the antimicrobial susceptibility testing of colistin (polymyxin E) as recommended by international standards (i.e. CLSI, EUCAST, ISO) but in a simpler and less time-consuming way.

## KIT CONTENT

- 4 Systems (panels) of ComASP Colistin (panels individually packed in foil with silica gel desiccant)
- 16 Tubes of Mueller Hinton II Broth (3.6 ml)
- Sealing Film
- Instructions Sheet and Test Results Form

## ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

- McFarland 0.5 Barium Sulphate Standard (ref.80400)
- Solution reservoir for multichannel pipette (ref. 96761)
- Physiological Solution (ref. 20095)
- Tips for multichannel pipette (ref.96758)
- Multichannel pipette 30-300 µl (ref.96759)

## CONFIGURATION

Test	Colistin Concentration (µg/ml)						
<b>A</b> Growth	0.25	0.5	1	2	4	8	16
<b>B</b> Growth	0.25	0.5	1	2	4	8	16
<b>C</b> Growth	0.25	0.5	1	2	4	8	16
<b>D</b> Growth	0.25	0.5	1	2	4	8	16

Growth indicates growth control: No antimicrobial agent in the well.

## PRINCIPLE OF THE METHOD

A panel allows the susceptibility testing of 4 different isolates. All wells in a single row (A, B, C or D) are rehydrated with a standardized microbial suspension. After incubation in thermostat the result is read and interpreted.

## COLLECTION AND STORAGE OF THE SAMPLE

ComASP Colistin is not for use directly with clinical or other specimens. The microorganism to be tested must first be isolated on a suitable non-selective culture medium. In case of mixed culture, selected colonies should be purified by subculturing.

## TEST PROCEDURE

1. Take a panel from its envelope and leave it at room temperature for 10 min,  
DO NOT DISCARD THE ENVELOPE until all 4 tests have been carried out.
2. Prepare a suspension of the test organism using either the direct colony suspension or growth method.
3. Standardize the suspension to the density of a McFarland 0.5 standard.
4. Optimally within 15 min of preparation, dilute the adjusted suspension 1:20 in saline; this will be the **Solution A**.
5. Add 0.4 ml of Solution A to a tube of MH II Broth\* provided in the kit to obtain the **Solution B**.
6. Dispense 100 µl of Solution B into each well in a row.
7. Cover the panel with the lid provided and incubate at  $36 \pm 2^\circ\text{C}$  for 16-20 hours in ambient air.

\* Mueller Hinton II Broth (g/l): Beef Extract 3.0 g; Acid Hydrolysate of Casein 17.5 g; Starch 1.5 g; Distilled Water 1000 ml;  
pH  $7.3 \pm 0.1$   
(adjusted with appropriate salts to provide 20-25 mg/l of calcium and 10-12.5 mg/l of magnesium)

## READING THE RESULTS

At the end of the incubation period observe the growth in the wells and establish the MIC, i.e. the lowest concentration of antibiotic that inhibits visible growth.

The result can be read visually or automatically. Reading by the naked eye can be improved by use of bright indirect lighting against a dark background.

Growth appears as turbidity or as a button at the bottom of the well (compare with the amount of growth in the growth-control well).

The growth-control well should be evaluated first. Make sure it is positive for growth. If not, check the viability of the colonies picked and repeat the test using a new row in the same panel or a new panel and a microbial culture of recent growth.

Note the results on the Test Results Form included in the kit (copy as many form as necessary).

## INTERPRETATION OF THE RESULTS

The MIC obtained should be interpreted according to current EUCAST or CLSI interpretative criteria.

**NOTE:** Each single panel can be used up to four times. If less than 4 tests have been performed, use the film provided in the kit to seal the inoculated rows so to prevent any leakage of contaminated fluids. Then, return the panel into its own desiccant envelope and into the refrigerator (see STORAGE).

## USER QUALITY CONTROL

Quality control of ComASP Colistin is performed using the following reference strains:

1. *Escherichia coli* ATCC® 25922; NCTC 12241
2. *Pseudomonas aeruginosa* ATCC® 27853; NCTC 12903
3. *Escherichia coli* NCTC 13846

Strain 1 and 2 are both susceptible to colistin and the acceptable MIC range is 0.25–2 µg/ml for *E. coli* and 0.5–4 µg/ml for *P. aeruginosa*. For *E. coli* NCTC 13846 (mcr-1 positive), which is instead resistant to colistin, the MIC target value is 4 µg/ml and should only on occasion be 2 or 8 µg/ml.

## FACTORS THAT MAY INVALIDATE THE RESULTS

Contaminated culture; poor standardization of the inoculum; clinical material unsuitable; use of expired panels or expired supplementary reagents; non compliance with temperatures and times of incubation.

## PRECAUTIONS

The product ComASP Colistin does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. ComASP Colistin is a disposable device to be used only for diagnostic use *in vitro*. The product must be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

## STORAGE

Store ComASP Colistin at 2-8°C in the original packaging. Once an envelope is opened the panel should be used within 7 days. Keep away from direct sunlight and direct heat. Do not use the panels beyond the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

## DISPOSAL OF USED MATERIAL

After use, ComASP Colistin and material that has come into contact with the sample must be decontaminated and disposed of in accordance with guidelines used in the laboratory for decontamination and disposal of potentially infected material.



## REFERENCES

- Carretto E. *et al.* Clinical validation of the SensiTest™ Colistin, a broth microdilution based method to evaluate colistin MICs. J Clin Microbiol. 2018;17. pii: JCM.01523-17. <http://jcm.asm.org/content/early/2018/01/11/JCM.01523-17>
- Carretto E. *et al.* Detection of mcr-4 positive Salmonella enterica serovar Typhimurium in clinical isolates of human origin, Italy, October to November 2016. Euro Surveill. 2018;23(2):pii=17-00821. <https://doi.org/10.2807/1560-7917.ES.2018.23.2.17-00821>
- CLSI. Performance Standards for Antimicrobial Susceptibility Testing; 28th ed. CLSI Supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
- CLSI. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; 11th ed. CLSI standard M07. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
- The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 8.0, 2018. <http://www.eucast.org>.
- The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. Version 8.0, 2018. <http://www.eucast.org>.
- ISO 20776-1:2006. Clinical laboratory testing and in vitro diagnostic test systems -- Susceptibility testing of infectious agents and evaluation of performance of antimicrobial susceptibility test devices -- Part 1: Reference method for testing the in vitro activity of antimicrobial agents against rapidly growing aerobic bacteria involved in infectious diseases.

## PRESENTATION

Product	µg/ml	Packaging	Ref.
ComASP™ Colistin	0.25 - 16	4x4 tests	75001

## TABLE OF SYMBOLS

<b>IVD</b> <i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limits
<b>REF</b> Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	<b>LOT</b> Batch code

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
## BIBLIOGRAFIA

1. Carretto E. *et al.* Clinical validation of the SensiTest™ Colistin, a broth microdilution based method to evaluate colistin MICs. J Clin Microbiol. 2018;17. pii: JCM.01523-17. <http://jcm.asm.org/content/early/2018/01/11/JCM.01523-17>
2. Carretto E. *et al.* Detection of mcr-4 positive Salmonella enterica serovar Typhimurium in clinical isolates of human origin, Italy, October to November 2016. Euro Surveill. 2018;23(2):pii=17-00821. <https://doi.org/10.2807/1560-7917.ES.2018.23.2.17-00821>
3. CLSI. Performance Standards for Antimicrobial Susceptibility Testing; 28th ed. CLSI Supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
4. CLSI. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; 11th ed. CLSI standard M07. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
5. The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 8.0, 2018. <http://www.eucast.org>.
6. The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. Version 8.0, 2018. <http://www.eucast.org>.
7. ISO 20776-1:2006. Clinical laboratory testing and in vitro diagnostic test systems -- Susceptibility testing of infectious agents and evaluation of performance of antimicrobial susceptibility test devices -- Part 1: Reference method for testing the in vitro activity of antimicrobial agents against rapidly growing aerobic bacteria involved in infectious diseases.

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

## REFERÊNCIAS

1. Carretto E. *et al.* Clinical validation of the SensiTest™ Colistin, a broth microdilution based method to evaluate colistin MICs. J Clin Microbiol. 2018;17. pii: JCM.01523-17. <http://jcm.asm.org/content/early/2018/01/11/JCM.01523-17>
2. Carretto E. *et al.* Detection of *mcr-4* positive Salmonella enterica serovar Typhimurium in clinical isolates of human origin, Italy, October to November 2016. Euro Surveill. 2018;23(2):pii=17-00821. <https://doi.org/10.2807/1560-7917.ES.2018.23.2.17-00821>
3. CLSI. Performance Standards for Antimicrobial Susceptibility Testing; 28th ed. CLSI Supplement M100S. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
4. CLSI. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; 11th ed. CLSI standard M07. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
5. The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 8.0, 2018. <http://www.eucast.org>.
6. The European Committee on Antimicrobial Susceptibility Testing. Routine and extended internal quality control for MIC determination and disk diffusion as recommended by EUCAST. Version 8.0, 2018. <http://www.eucast.org>.
7. ISO 20776-1:2006. Clinical laboratory testing and *in vitro* diagnostic test systems -- Susceptibility testing of infectious agents and evaluation of performance of antimicrobial susceptibility test devices -- Part 1: Reference method for testing the *in vitro* activity of antimicrobial agents against rapidly growing aerobic bacteria involved in infectious diseases.

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