ESC swab (ISO 18593, ISO 17604)

Buffered Peptone Water

DESCRIPTION


ESC (Easy Surface Checking) swab - Buffered Peptone Water is a device used for the microbiological monitoring of surfaces by swab sampling method. The medium is recommended for the detection and enumeration of *Salmonella* spp, Enterobacteriaceae and *Listeria* spp within foodstuffs.

CONTENTS OF THE PACKAGES

The product consists of a transparent plastic tube containing liquid medium and a sterile swab attached to the cap.

Each package contains:

- 100 tubes with swab and liquid medium;
- 1 instructions sheet.

METHOD PRINCIPLE

Buffered Peptone Water contains peptone and sodium chloride in a phosphate buffer. Its formula allows recovery of even sub-lethally injured cells from food and associated samples prior to selective enrichment and isolation of target organism. No neutralizing agents are included. For areas where residual disinfectants are presumably present, it is recommended that an alternative collection medium be used such as Neutralizing Rinse Solution (ref. 85601), D/E Neutralizing Broth (ref. 85605), or Letheen Broth (ref. 85607).

FORMULA (per liter of water)

Peptone 10 g, sodium chloride 5 g, disodium hydrogen phosphate 3.57 g, potassium dihydrogen phosphate 1.5 g, pH 7.0 ± 0.2 at 25°C.

TEST PROCEDURE

1. Take a tube from the package and unscrew the cap with an attached swab.
2. Press the tip of the swab against the wall of the tube to remove excess liquid.
3. Streak the swab on a surface horizontally and vertically while rotating swab. A Sampling Template 10x10 (ref. 96762) can be used to delimit and sample an area of 100 cm².
4. Return the swab to the tube.
5. Screw the cap and record date and the sampling point.
6. Transport the device to the laboratory and examine according to methods in use. Incubation during the pre-enrichment step can be directly performed in the tube itself.

QUALITY CONTROL

Each batch of ESC swab - Buffered Peptone Water is submitted to quality control using the following strains:
S. Typhimurium WDCM 00031, S. Enteritidis WDCM 00030, E. coli WDCM 00012, L. monocytogenes 4b WDCM 00021.

STORAGE

10-25°C away from light, until the expiry date on the label. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

ESC Swab - Buffered Peptone Water does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product must be used only by properly trained operators.

BIBLIOGRAPHY

- ISO 18593:2018. Microbiology of the food chain – Horizontal method for surface sampling.
- ISO 17604:2015. Microbiology of the food chain – Carcass sampling for microbiological analysis.
- ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
- ISO 6887-1:2017. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination. Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
- ISO 21528-1:2017. Microbiology of the food chain - Horizontal method for the detection and enumeration of Enterobacteriaceae – Part 1: Detection of Enterobacteriaceae.
- ISO 6579-1:2017. Microbiology of the food chain – Horizontal method for the detection, enumeration and serotyping of *Salmonella* – Part 1: Detection of *Salmonella* spp.

PRESENTATION	Ref.	Contents
ESC swab - Buffered Peptone Water	85603	100 x 10 ml tubes
ESC swab - Buffered Peptone Water	85604	100 x 5 ml tubes

TABLE OF SYMBOLS

LOT Batch code	Fragile, handle with care	Manufacturer	Use by
REF Catalogue number	Do not reuse	Contains sufficient for <n> tests	Caution, consult Instruction For Use

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
BIBLIOGRAFIA

- ISO 18593:2018. Microbiology of the food chain – Horizontal method for surface sampling.
- ISO 17604:2015. Microbiology of the food chain – Carcass sampling for microbiological analysis.
- ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
- ISO 6887-1:2017. Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination. Part 1: General rules for the preparation of the initial suspension and decimal dilutions.
- ISO 21528-1:2017. Microbiology of the food chain - Horizontal method for the detection and enumeration of *Enterobacteriaceae* – Part 1: Detection of *Enterobacteriaceae*.
- ISO 6579-1:2017. Microbiology of the food chain – Horizontal method for the detection, enumeration and serotyping of *Salmonella* – Part 1: Detection of *Salmonella* spp.

PRESENTAZIONE	Ref	Contenuto
ESC swab - Buffered Peptone Water	85603	Provette 100 x 10 ml
ESC swab - Buffered Peptone Water	85604	Provette 100 x 5 ml

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