ESC swab (ISO 18593, ISO 17604)

D/E Neutralizing Broth

DESCRIPTION

ESC (Easy Surface Checking) swab - D/E Neutralizing Broth is a device used for the microbiological monitoring of surfaces by swab sampling method. Neutralizing agents included in the medium allow to assess the effectiveness of sanitization procedures in areas where residual disinfectants are presumably present. D/E Neutralizing Broth inactivates higher concentrations of residual antimicrobials when compared with other standard neutralizing formulations.

CONTENTS OF THE PACKAGES

The product consists of a transparent plastic tube containing the collection broth and a sterile swab attached to the cap.

Each package contains:

- 100 tubes with swab and liquid medium;
- 1 instructions sheet.

METHOD PRINCIPLE

The medium containing thioglycollate, sodium thiosulphate, sodium bisulphite, lecithin and polysorbate 80 has the ability to neutralize a broad spectrum of antiseptic and disinfectant chemicals including quaternary ammonium compounds, phenolics, iodine, chlorine preparations, mercurials, formaldehyde and glutaraldehyde. Bromocresol purple is incorporated as pH indicator to demonstrate the production of acid from the fermentation of glucose.

FORMULA (per liter of medium)

Enzymatic digest of casein 5 g, yeast extract 2.5 g, dextrose 10 g, sodium thioglycollate 1 g, sodium thiosulphate 6 g, sodium bisulfite 2.5 g, lecithin 7.0 g, polysorbate 80 5 g, bromocresol purple 0.02 g, final pH 7.6 ± 0.2 at 25°C.

QUALITY CONTROL

Each batch of ESC swab - D/E Neutralizing Broth is submitted to quality control using the following strains: *E. coli* ATCC® 25922, *S. aureus* ATCC® 6538, *S. Typhimurium* ATCC® 14028, *B. subtilis* ATCC® 6633, *P. aeruginosa* ATCC® 27853.

TEST PROCEDURE

1. Take a tube from the package and unscrew the cap with an attached swab.
2. Press the tip of the swab against the wall of the tube to remove excess liquid.
3. Streak the swab on a surface horizontally and vertically. A Sampling Template 10x10 (ref. 96762) can be used to delimit and sample an area of 100 cm².
4. Return the swab to the tube.
5. Screw the cap and record date and the sampling point.
6. Transport the device to the laboratory and examine according to methods in use.

STORAGE

10-25°C away from light, until the expiry date on the label. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

ESC swab - D/E Neutralizing Broth does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use only and must be used by properly trained operators.


BIBLIOGRAPHY

- ISO 18593:2018. Microbiology of the food chain - Horizontal method for surface sampling.
- ISO 17604:2015. Microbiology of food the food chain - Carcass sampling for microbiological analysis.
- Dey B.P. and Engley Jr F.B. (1995) Comparison of Dey and Engley (D/E) Neutralising medium to Lethen medium and Standards Methods Medium for recovery of *Staphylococcus aureus* from sanitised surfaces. J. Ind. Microbiol. 14:21-25.

PRESENTATION

PRESENTATION	Ref.	Contents
ESC swab - D/E Neutralizing Broth	85605	100 x 10 ml tubes
ESC swab - D/E Neutralizing Broth	85606	100 x 5 ml tubes

TABLE OF SYMBOLS

LOT Batch code	 Fragile, handle with care	 Manufacturer	 Use by
REF Catalogue number	 Do not reuse	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use

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