E.coli Test

Rapid test for detection of β -galactosidase activity, indole production and MUG fluorescence.

INTENDED PURPOSE

Rapid test for *Escherichia coli* presumptive identification. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

E.coli Test is used for *E. coli* detection through the qualitative determination of activity of three enzymes: β -galactosidase, tryptophanase and β -glucuronidase.

KIT CONTENT

- 30 Tubes of E.coli Test
- 1 Tube of Kovac's Reagent (3.0 ml)

METHOD PRINCIPLE

The enzyme β -galactosidase, which allows bacteria to metabolize lactose by splitting it into glucose and galactose, hydrolyzes the colourless substrate ONPG (ortho-nitrophenyl- β -D-galactopyranoside) with the formation of o-nitrophenol, a yellow compound. Bacteria that produce the enzyme tryptophanase are able to degrade the amino acid tryptophan into pyruvic acid, ammonium and indole. The indole produced reacts with p-dimethylaminobenzaldehyde of the Kovac's (indole) reagent to form a red-violet compound. The enzyme β -glucuronidase cleaves the colourless substrate MUG (4-methylumbelliferyl- β -D-glucuronide) to produce 4-methylumbelliferone, which fluoresces under long-wave ultra-violet light (approximately 366 nm).

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: incubator, inoculating loop, pipettes, physiological solution (0.85% saline), culture media, Wood's lamp, quality control organisms.

REAGENTS

- A tube of E.coli Test contains desiccated biochemical substrates and nutritive components such as ONPG, MUG and tryptone, which includes tryptophane.
- Kovac's Reagent consists of 5% p-dimethylaminobenzaldehyde dissolved in a solution of 25% hydrochloric acid and 75% isobutyl alcohol.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Take the number of E.coli Test tubes needed from the fridge and allow them to reach room temperature.
2. Add 0.25 ml of physiological solution to each tube.
3. With a sterile loop pick up one or more than one well isolated colony from a fresh culture plate (18-24 hours old) and emulsify in the tube.
4. Incubate tubes at $35 \pm 2^\circ\text{C}$ for up to 24 hours.
5. Examine tubes at intervals of 1, 4 and 24 hours for broth colour change and fluorescence under long-wave UV light.
6. Add 2 drops of Kovac's Reagent to perform the Indole Test.

Note: Positive and negative controls should be run simultaneously with the organism to be tested (see QUALITY CONTROL).

INTERPRETING RESULTS

A broth colour change to yellow (ONPG positive), blue fluorescence (MUG positive) and development of a red ring (Indole positive) at the surface of the broth within 30 seconds after addition of Kovac's Reagent, indicate the presence of *E. coli*.

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

Appearance of E.coli Test: yellowish powder, yellow solution once reconstituted

Appearance of Kovac's Reagent: yellow-green liquid, clear

Control strain		Characteristic reactions
<i>Klebsiella pneumoniae</i>	ATCC® 13383	Positive ONPG reaction: Development of a yellow colour Negative MUG reaction: No fluorescence Negative indole reaction: Yellow/brown ring
<i>Escherichia coli</i>	ATCC® 25922	Positive ONPG reaction: Development of a yellow colour Positive MUG reaction: Fluorescence Positive indole reaction: Formation of a red ring
<i>Proteus mirabilis</i>	ATCC® 25923	Negative ONPG reaction: No colour change Negative MUG reaction: No fluorescence Negative indole reaction: Yellow/brown ring

PERFORMANCE CHARACTERISTICS

Performance testing of E.coli Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

This test is only part of the overall scheme for the identification of *E. coli*. Further biochemical and/ or serological testing is required for definitive identification.

Due to nutritional variation, some strains may result in poor growth or fail to grow on this medium providing false negative results.

Inoculation of test organism should be done only from lactose-containing medium.

The reaction speed depends on the size of the inoculum.

Allow at least 30 seconds for the color to develop before considering the indole test negative. Pigmented organisms may yield ambiguous results.

Indole produced by certain anaerobes, such as *Clostridium* species, can be rapidly degraded resulting in false negative reactions.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**

- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
E.coli Test	30 tests	88006

Revision History

Revision	Release Date	Change Summary
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

This IFU document and the SDS are available from the online Support Center:

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References / Riferimenti

1. Public Health England. (2018). ONPG (β -Galactosidase) test. UK Standards for Microbiology Investigations. TP 24 Issue 4. <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
2. Lenette EH (1995) Manual of Clinical Microbiology.
3. Blazevic DJ, Ederer GM (1975) Principles of biochemical tests in diagnostic microbiology. 63-67. New York, John Wiley & Sons Lapage SP, Jayaraman MS (1964). Beta-galactosidase and lactose fermentation in the identification of enterobacteria including salmonellae. J Clin Path. 17:117.
4. Lederberg J (1950) The β -galactosidase of escherichia coli, strain K-12. J Bact. 60:381.

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