Hippurate Test

Rapid test tube for the detection of Hippurate hydrolysis.

INTENDED PURPOSE

Kit to perform hippurate test on bacterial isolate. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Hippurate Test is used to detect the hydrolysis of sodium hippurate by certain bacteria. This test allows for example to differentiate Lancefield Group B Streptococci (GBS) from other β -hemolytic streptococci.

KIT CONTENT

- 1 Vial of Ninhydrin Reagent (5 ml) + dropper bottle
- 1 Cartridge with 30 discs of Hippurate Test
- 30 Empty Tubes

METHOD PRINCIPLE

Group B streptococci (*Streptococcus agalactiae*) and some other bacteria such as *Gardnerella vaginalis*, *Legionella pneumophila*, and *Campylobacter jejuni*, produce the enzyme hippurate hydrolase (hippuricase) which hydrolyzes hippuric acid to sodium benzoate and glycine. The hydrolysis of hippuric acid is detected with ninhydrin, a strong oxidizing agent that deaminates α -amino groups of glycine with the release of ammonia and carbon dioxide. Then, ammonia reacts with residual ninhydrin to form a blue to purple colored complex.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loop, pipettes, culture media, physiological solution, quality control organisms.

REAGENTS

- Ninhydrin Reagent: 7% ninhydrin dissolved in an organic solvent.
- Hippurate Test discs: paper discs containing approximately 8 mg of sodium hippurate.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Add 0.2 ml of sterile physiological solution to a test tube.
2. Inoculate with 1-3 colonies of a fresh isolate and mix well.
3. Using forceps, drop a Hippurate Test disc into the microbial suspension.
4. Incubate tubes at $35 \pm 2^\circ\text{C}$ for 2 hours.
5. Add 2 drops of Ninhydrin Reagent to the test tube, shake gently, re-incubate for 10 minutes. Do not exceed 30 minutes as false reactions may occur.

INTERPRETING RESULTS

Observe for a color change of the test suspension: A positive hippurate reaction is indicated by a color change of the test suspension to blue-purple. No color change is considered as a negative result.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

2 years.

QUALITY CONTROL**Expected Response:**

Appearance of Ninhydrin Reagent: Yellow liquid, clear.

Control strain		Characteristic reactions
<i>Streptococcus agalactiae</i>	ATCC® 13813 NCTC 8181	Positive reaction: Color change to blue
<i>Streptococcus pyogenes</i>	ATCC® 19615 NCTC 12696	Negative reaction: No color change

Please refer to the actual batch related Certificate of Analysis (CoA).

PERFORMANCE CHARACTERISTICS

Performance testing of Hippurate Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

The Hippurate rapid test can be used as a presumptive diagnostic aid for identification of microorganisms capable of hydrolyzing hippurate, such as *Streptococcus agalactiae*. However further tests are needed to complete the diagnostic results.

Some strains of group D streptococci may give a weak positive result. Group B and D streptococci can be distinguished by subculture on Bile Aesculin Agar. Only group D streptococci hydrolyze esculin yielding a black color of the medium.

Incubation of the test culture for more than 30 min after the addition of Ninhydrin Reagent can give false-positive results.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Hippurate Test	30 tests	88007

Revision History

Revision	Release Date	Change Summary
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

This IFU document and the SDS are available from the online Support Center:

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