Listeria Mono Test

Rapid test for differentiation of *Listeria monocytogenes* from other *Listeria* species.

INTENDED PURPOSE

Rapid test for *Listeria monocytogenes* identification from clinical and non-clinical specimens. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Listeria Mono Test is for rapid screening and differentiation of *Listeria* species based on qualitative determination of D-alanyl-aminopeptidase enzyme activity.

KIT CONTENT

- 20 Tubes of Listeria Mono Test
- 1 Cartridge of 20 ABN discs
- 1 Tube of Listeria Mono Reagent (3.0 ml)

METHOD PRINCIPLE

The genus *Listeria* includes the pathogen *Listeria monocytogenes* which causes human listeriosis. *Listeria* species, with the exception of *L. monocytogenes*, produce the enzyme D-alanyl-aminopeptidase (DALAase), which hydrolyzes the substrate aminoacyl- β -naphthylamide (ABN) to give a yellow colour after combining with the reagent p-dimethylamino-benzaldehyde.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: incubator, inoculating loop, pipettes, culture media, physiological solution (0.85% NaCl), quality control organisms.

REAGENTS

- A tube of Listeria Mono Test contains desiccated biochemical substrates and nutritive components
- Each ABN disc is impregnated with a solution of aminoacyl- β -naphthylamide
- Listeria Mono Reagent contains p-dimethylamino-benzaldehyde

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

Before following below protocol, ensure that the culture to be tested is a Gram-positive, catalase-positive, oxidase-negative organism.

1. To minimize condensation on the discs, allow the cartridge to reach room temperature before use.
2. Add 0.2 ml (2 drops) of sterile physiological solution to a tube of Listeria Mono Test.
3. Inoculate heavily with the test organism from a fresh overnight pure culture.
4. Incubate the inoculated tube at $35 \pm 2^\circ\text{C}$ for 12 to 24 hours.
5. Put one ABN disc into the broth culture and incubate again at $35 \pm 2^\circ\text{C}$ for 4 hours.
6. Add 2 drops of Listeria Mono Reagent to the broth.
7. Examine for the appearance of a yellow color.

Note: Positive and negative controls should be run simultaneously with the organism to be tested (see QUALITY CONTROL).

INTERPRETING RESULTS

A positive D-alanyl-aminopeptidase reaction is indicated by the development of a yellow colour, which means that the organism belongs to *Listeria* species other than *Listeria monocytogenes*.

No colour development indicates that the organism is a presumptive *Listeria monocytogenes* (negative D-alanyl-aminopeptidase reaction).

STORAGE

Store at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

Control strain		Incubation	Characteristic reactions
<i>Listeria monocytogenes</i>	ATCC® 35152	18-24 h/ 35 ± 2°C	Negative reaction: No yellow colour development
<i>Listeria innocua</i>	ATCC® 33090		Positive reaction: Yellow colour development
<i>Listeria seeligeri</i>	ATCC® 35967		
<i>Listeria gray</i>	ATCC® 25401		

PERFORMANCE CHARACTERISTICS

Performance testing of Listeria Mono Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

This test provides presumptive identification of *Listeria monocytogenes* but does not replace full biochemical testing.

Due to the small size of *Listeria* colonies on primary isolation media, it is not recommended to carry out the test directly from these plates as use of multiple colonies from the primary isolation plate may lead to a mixed culture and an incorrect result.

Occasionally certain *Listeria* when grown on non-selective media form colonies that are difficult to remove with an inoculating loop. Ensure that there is sufficient colonial material on the inoculating loop before inoculating the test tube. Failure to pick up sufficient material on the loop may lead to false-negative results.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Listeria Mono Test	20 tests	88010

Revision History

Revision	Release Date	Change Summary
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

This IFU document and the SDS are available from the online Support Center:

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