Urea Rapid Test

Rapid test for detection of urease production

INTENDED PURPOSE

Rapid test to evidence urease production on bacterial isolates. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Urea Rapid Test is a rapid test to detect the production of urease by *Proteus*, *Klebsiella* and certain yeasts, such as an identification test for the differentiation of *Salmonella* spp (urease negative).

KIT CONTENT

- 30 tubes of Urea Rapid Test

METHOD PRINCIPLE

Microorganisms that produce enzyme urease hydrolyze urea in ammonia, water and carbon dioxide. The presence of alkaline products determines colour change of phenol red indicator from yellow to red-fuchsia.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: Incubator, inoculating loop, pipettes, culture media, physiological solution (0.85% saline), test tubes, 0.5 McFarland turbidity standard, quality control organisms.

REAGENTS

- A tube of Urea Rapid Test contains desiccated proteic substrata, the pH indicator phenol red and urea.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Take the number of Urea Rapid Test tubes needed from the fridge and allow them to reach room temperature.
2. Add 0.3 mL of physiological solution to each tube.
3. Inoculate heavily with the test organism (*) from a fresh overnight pure culture.
4. Incubate at 35 ± 2 °C for up to 24 hours (**)
5. Examine tubes at intervals of 1,4 and 24 hours for broth colour change.

*Microbial suspension equivalent to 0.5 McFarland turbidity standard (which corresponds to $1-2 \times 10^8$ CFU/ml for *E. coli*)

**Prolonged incubation may result in a false positive test due to hydrolysis of proteins in the medium.

Note: positive and negative controls should be run simultaneously with the organism to be tested (SEE QUALITY CONTROL).

INTERPRETING RESULTS

A positive urease reaction is indicated by a colour change to red Fuchsia.

Further tests are required for confirmation.

STORAGE

2-8 °C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

Appearance of Urea Rapid Test: yellowish powder, yellow solution once reconstituted.

Control strain		Inoculum	Incubation	Characteristic reactions
<i>Escherichia coli</i>	WDCM 00013 ATCC® 25922 NCTC 12241	1-2 x 10 ⁸ CFU/ml	Up to 24 h/ 35 ± 2°C	Negative reaction: No colour change
<i>Proteus mirabilis</i>	WDCM 00023 ATCC® 29906 NCTC 11938			Positive reaction: Colour change to red fuchsia
<i>Salmonella</i> Typhimurium	WDCM 00031 ATCC®14028 NCTC 12023			Negative reaction: No colour change

PERFORMANCE CHARACTERISTICS

Performance testing of Urea Rapid Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Due to nutritional variation, some strains may result in poor growth or fail to grow on this medium providing false negative results.

Allow at least 30 seconds for the color to develop before considering the test negative. Pigmented organisms may yield ambiguous results.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Urea Rapid Test	30 tests	88011

Revision History

Revision	Release Date	Change Summary
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

This IFU document and the SDS are available from the online Support Center:

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References / Riferimenti

1. EN ISO 11133:2014+Amd1:2018+Amd2:2020. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. Public Health England. (2019). Urease test. UK Standards for Microbiology Investigations. TP 36 Issue 4. <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
3. Public Health England. (2018). Indole test. UK Standards for Microbiology Investigations. TP 19 Issue 4
4. Iseberg H.D. and Sundheim L.H.(1980) Indole reactions in bacteria. J. Bacteriol. 75:682-690. Baltimore, Williams & Wilkins.
5. MacFaddin J.F. (1980) Biochemical tests for identification for medical bacteria. 2nd ed. 173-183. Baltimore, Williams & Wilkins.
6. Blazevic D.J. and Ederer, G.M.(1975) Principles of biochemical tests in diagnostic microbiology. 63-67. New York, John Wiley & Sons.

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