H₂S Rapid Test

Rapid test for hydrogen sulfide production.

INTENDED PURPOSE

Rapid test for evidence hydrogen sulfide production on bacterial isolates. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

H₂S Rapid Test is for the identification of microorganisms based on their ability to produce hydrogen sulfide (H₂S).

KIT CONTENT

- 30 Tubes of H₂S Rapid Test

METHOD PRINCIPLE

This test consists of a desiccated medium which includes sodium thiosulfate and ferrous ammonium sulfate as indicators of hydrogen sulfide production. The ferrous ammonium sulfate reacts with H₂S gas to produce ferrous sulfide, a black precipitate

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loop, pipettes, culture media, physiological solution (0.85% NaCl), vaseline (mineral) oil, quality control organisms.

REAGENTS

- A tube of H₂S Rapid Test contains desiccated biochemical substrates and nutritive components, including sodium thiosulfate and ferrous ammonium sulfate

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Take the number of H₂S Rapid Test tubes needed from the fridge and allow them to reach room temperature.
2. Add 0.3 ml of sterile physiological solution to each tube.
3. Inoculate heavily with the test organism from a fresh overnight pure culture.
4. Overlay with 2-3 drops of sterile vaseline oil.
5. Incubate at 35 ± 2°C for up to 24 hours.
6. Examine tubes after 4-hour incubation and at the end of the incubation period for H₂S production.

Note: Positive and negative controls should be run simultaneously with the organism to be tested (see QUALITY CONTROL).

INTERPRETING RESULTS

H₂S formation (positive test) is indicated by a blackening of the medium. A negative H₂S test is denoted by the absence of blackening.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

Appearance of H₂S Rapid Test: Yellowish powder, light amber solution once reconstituted.

Control strain		Incubation	Characteristic reactions
<i>Klebsiella pneumoniae</i>	ATCC® 13883	Up to 24 h/ 35 ± 2°C	Negative reaction: No blackening
<i>Escherichia coli</i>	ATCC® 25922		
<i>Proteus mirabilis</i>	ATCC® 25933		Positive reaction: Blackening of the medium
<i>Salmonella Typhimurium</i>	ATCC® 14028		

PERFORMANCE CHARACTERISTICS

Performance testing of H₂S Rapid Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Due to nutritional variation, some strains may result in poor growth or fail to grow on this medium providing false negative results.

For identification, organisms must be in pure culture. Morphological, biochemical, and/or serological tests should be performed for final identification.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
H ₂ S Rapid Test	30 tests	88013

Revision History

Revision	Release Date	Change Summary
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

This IFU document and the SDS are available from the online Support Center:

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