Arginine Decarboxylase Test

Rapid test for detecting arginine decarboxylase activity.

INTENDED PURPOSE

Rapid test for evidence arginine decarboxylase production on bacterial isolates. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Arginine Decarboxylase Test is used for the identification and differentiation of Enterobacteriaceae and other gram-negative organisms by means of decarboxylase reactions.

KIT CONTENT

- 30 Tubes of Arginine Decarboxylase Test

METHOD PRINCIPLE

The enzyme arginine decarboxylase attacks the carboxylic group of the amino acid arginine with formation (combination with dihydrolase) of amine putrescine. The sugar contained in the medium is fermented by all the enterobacteria with consequent initial colour change of indicator system from purple to yellow. The acid medium supports the amino acid decarboxylation reactions with formation of putrescine. The amines production alkalizes the medium and induces a new indicator colour change from yellow to purple. The negative reaction is shown by the appearance of a yellow colour in the tube, while the positive reaction is shown by a purple colour.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loop, pipettes, culture media, physiological solution, mineral (vaseline) oil, quality control organisms.

REAGENTS

- A tube of Lysine Decarboxylase Test contains desiccated biochemical substrates and nutritive components including arginine, glucose, bromocresol purple and pyridoxal.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Take the number of Arginine Decarboxylase Test tubes needed from the fridge and allow them to reach room temperature.
2. Add 0.3 mL of physiological solution to each tube.
3. Inoculate with the test organism from a fresh overnight pure culture.
4. Overlay each inoculated tube with (2-3 drops) sterile mineral oil to separate the medium from external oxygen and promote fermentation.
5. Incubate tubes with caps tightened at $35 \pm 2^\circ\text{C}$ for up to 96 hours.
6. Examine tubes after 4 hours incubation and then daily for colour reactions.

INTERPRETING RESULTS

A positive result is indicated by the development of a purple colour.

A negative test is given by the development of a yellow colour for glucose-fermenting microorganisms, while non-glucose-fermenters will result in little or no colour change as compared to an uninoculated tube.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

Appearance: Purplish powder, yellowish-red solution once reconstituted

Control strain		Incubation	Characteristic reactions
<i>Enterobacter cloacae</i>	ATCC® 13047	18-24 h/ 35 ± 2°C	Positive reaction: Development of a purple color
<i>Pseudomonas aeruginosa</i>	ATCC® 27853		
<i>Proteus mirabilis</i>	ATCC® 25923		Negative reaction: Development of a yellow color
<i>Klebsiella pneumoniae</i>	ATCC® 13383		

PERFORMANCE CHARACTERISTICS

Performance testing of Arginine Decarboxylase Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Use of a heavy inoculum may be necessary to obtain a positive result in 4 h.

Interpretation earlier than 12-18 h incubation may lead to false-negative results for organisms that ferment glucose.

Exposure to air may cause alkalization of the surface of the medium giving false-positive results.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Arginine Decarboxylase Test	30 tests	88016

Revision History

Revision	Release Date	Change Summary
1	2024-09-16	Added notice to report any malfunction, defect or incident.
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

In case of malfunctions or defects, contact immediately Liofilchem (*) or the local representative.

In case of incident associated with the device, notify immediately Liofilchem (*) or its local representative and the National Competent Authority.

*Please login to <https://www.liofilchemstore.it/login.php> (user ID and password required) and click on Complaint.

This IFU document and the SDS are available from the online Support Center:

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References / Riferimenti

1. Lal A, Cheeptham N (2015) Decarboxylase Broth Protocol. American Society for Microbiology.
2. Lenette EH (1995) Manual of Clinical Microbiology.
3. Blazevic DJ, Ederer GM (1975) Principles of biochemical tests in diagnostic microbiology. 63-67. New York, John Wiley & Sons.
4. Brooks K, Sodeman T (1974) A rapid method for determining decarboxylase and dihydrolase activity. J Clin Pathol. 27(2):148-52. DOI: 10.1136/jcp.27.2.148.
5. Goldschmidt MC, Lockhart BM (1971) Rapid methods for determining decarboxylase activity: arginine decarboxylase. Appl Microbiol. 22(3):350-7. DOI: 10.1128/am.22.3.350-357.1971.

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