Indole Test

Rapid test for detection of indole production.

INTENDED PURPOSE

Rapid test for evidence indole production in bacterial isolates. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Indole test is used for the rapid qualitative determination of indole production as a first step in screening non-lactose fermenting colonies from agar plates.

KIT CONTENT

- 30 tubes of Indole Test
- 1 Kovac's reagent (3.0 ml)

METHOD PRINCIPLE

Bacteria that produce the enzyme tryptophanase can degrade the amino acid tryptophan into pyruvic acid, ammonia and indole. The indole produced reacts with p-dimethylaminobenzaldehyde of the Kovac's (indole) reagent to form a red-violet compound.

The production of indole from tryptophan is a characteristic absent in *Salmonella* but present in *E.coli*, *Morganella* and some species of *Klebsiella*, *Aerobacter* and *Citrobacter*.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: incubator, inoculating loop, pipettes, culture media, physiological solution (0,85% saline), test tubes, 0.5 McFarland turbidity standard, quality control organisms.

REAGENTS

- A tube of Indole Test contains desiccated biochemical substrates and nutritive components such as phenol red and tryptone, which includes tryptophane.
- Kovac's Reagent consist of 5% p-dimethylaminobenzaldehyde dissolved in a solution of 25% hydrochloric acid and 75% isobutyl alcohol.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Take the number of Indole Test tubes needed from the fridge and allow them to reach room temperature.
2. Add 0.3 ml of physiological solution to each tube.
3. Inoculate heavily with the test organism (*) from a fresh overnight pure culture.
4. Incubate at 35 ± 2 °C for up to 24 hours (**).
5. Examine tubes at intervals of 1,4 and 24 hours for broth colour change.
6. Add 2-3 drops of Kovac's Reagent to perform the Indole test.

*Microbial suspension equivalent to 0.5 McFarland turbidity standard (which corresponds to $1-2 \times 10^8$ CFU/ml for *E.coli*).

**Prolonged incubation may result in a false-positive test due to hydrolysis of proteins in the medium.

Note: Positive and negative controls should be run simultaneously with the organism to be tested (see QUALITY CONTROL)

INTERPRETING RESULTS

The appearance of a red ring at the surface of the broth within 30 seconds (after addition of Kovac's Reagent) indicates a positive indole reaction.

Further tests are required for confirmation.

STORAGE

2-8 °C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

Appearance of Indole Test: Yellowish powder, yellow solution once reconstituted.

Appearance of Kovac's Reagent: Yellow-green liquid, clear.

Control strain		Inoculum	Incubation	Characteristic reactions
<i>Escherichia coli</i>	WDCM 00013 ATCC® 25922 NCTC 12241	1-2 x 10 ⁸ CFU/ml	Up to 24 h/ 35 ± 2°C	Positive reaction: Formation of a red ring
<i>Proteus mirabilis</i>	WDCM 00023 ATCC® 29906 NCTC 11938			Negative reaction: Yellow/ brown ring
<i>Salmonella Tiphymurium</i>	WDCM 00031 ATCC® 14028 NCTC 12023			Negative reaction: Yellow ring

PERFORMANCE CHARACTERISTICS

Performance testing of Indole Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Due to nutritional variation, some strains may result in poor growth or fail to grow on this medium providing false negative results.

Allow at least 30 seconds for the color to develop before considering the indole test negative. Pigmented organisms may yield ambiguous results.

Indole production by certain anaerobes, such as *Clostridium* species, can be rapidly degraded resulting in false negative reactions.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.

- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Indole Test	30 tests	88017

Revision History

Revision	Release Date	Change Summary
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

This IFU document and the SDS are available from the online Support Center:

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


References / Riferimenti

1. EN ISO 11133:2014+Amd1:2018+Amd2:2020. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
2. Public Health England. (2019). Urease test. UK Standards for Microbiology Investigations. TP 36 Issue 4. <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
3. Public Health England. (2018). Indole test. UK Standards for Microbiology Investigations. TP 19 Issue 4
4. Iseberg H.D. and Sundheim L.H.(1980) Indole reactions in bacteria. J. Bacteriol. 75:682-690. Baltimore, Williams & Wilkins.
5. MacFaddin J.F. (1980) Biochemical tests for identification for medical bacteria. 2nd ed. 173-183. Baltimore, Williams & Wilkins.
6. Blazevic D.J. and Ederer, G.M.(1975) Principles of biochemical tests in diagnostic microbiology. 63-67. New York, John Wiley & Sons.

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