Sugar Fermentation

For the differentiation of microorganisms.

INTENDED PURPOSE

Tests for sugar fermentation as aid for identification of both Gram-positive and Gram-negative bacteria isolated from clinical and non-clinical specimens. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Sugar Fermentation are for differentiation of microorganisms based on their ability to ferment carbohydrates. The tubes contain a desiccated protein mixture and the sugar substrate of the fermentation.

KIT CONTENT

- 30 Tubes of the specific "Sugar Fermentation" Test

For example, the product ref. number 88208, Adonitol Test, contains 30 tubes each containing adonitol as a single carbohydrate.

The only exception is the product ref. number 88220, Xylose Test, consisting of 30 tubes and 1 cartridge with 30 Xylose Discs.

METHOD PRINCIPLE

Carbohydrate fermentation tests detect the ability of a microorganism to metabolize a specific carbohydrate allowing to differentiate among bacterial groups or species. Fermentation reactions occurring in a culture media are detected by the color change of a pH indicator when acid products are formed.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loop, pipettes, culture media, physiological solution (0.85% saline), quality control organisms.

REAGENTS

- A tube of "Sugar Fermentation" Test contains desiccated biochemical substrates and nutritive components including the specific carbohydrate and bromocresol purple as a pH indicator.

For Xylose Test, the carbohydrate is impregnated in the discs provided with the kit.

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Allow product to reach room temperature before use (take the tubes needed only).
2. Add 0.3 mL of physiological solution to each tube. In case of using Xylose Test, place a xylose disc in the tube containing the just rehydrated substrates.
3. Using a sterile loop, pick up one or more than one well isolated colony from a fresh culture plate (18-24 hours old) to make a suspension of the test organism in the test tube.
4. Cap the tube and incubate at $35 \pm 2^\circ\text{C}$ for 18-24 hours. Incubate Glucose Test for 16-18 hours.
5. Repeat with another tube but skip step 3 as it will serve as a negative control (uninoculated test).

INTERPRETING RESULTS

A positive reaction is indicated by a colour change from purple to yellow.

If there is no color change (the medium remains purple) the test is negative.

STORAGE

Store at 2-8°C in its original package. In such conditions the product is valid until the expiry date shown on the label. Do not use them beyond that date. Dispose of them if they show signs of deterioration.

SHELF LIFE

1 year.

QUALITY CONTROL

Description	Positive control	Negative control
Adonitol Test	<i>Klebsiella pneumoniae</i> ATCC® 13883	<i>Escherichia coli</i> ATCC® 25922
Arabinose Test	<i>Escherichia coli</i> ATCC® 25922	<i>Enterococcus faecalis</i> ATCC® 19433
Arabitol Test	<i>Listeria monocytogenes</i> ATCC® 35152	<i>Proteus mirabilis</i> ATCC® 25933
Dulcitol Test	<i>Klebsiella pneumoniae</i> ATCC® 13883	<i>Escherichia coli</i> ATCC® 25922
Galactose Test	<i>Escherichia coli</i> ATCC® 25922	<i>Pseudomonas aeruginosa</i> ATCC® 27853
Glucose Test	<i>Escherichia coli</i> ATCC® 25922	<i>Plesiomonas shigelloides</i> ATCC® 14029
Inositol Test	<i>Klebsiella pneumoniae</i> ATCC® 13883	<i>Escherichia coli</i> ATCC® 25922
Inulin Test	<i>Streptococcus pneumoniae</i> ATCC® 6303	<i>Streptococcus pyogenes</i> ATCC® 19615
Lactose Test	<i>Escherichia coli</i> ATCC® 25922	<i>Salmonella typhimurium</i> ATCC® 14028
Levulose Test	<i>Streptococcus pneumoniae</i> ATCC® 6303	<i>Neisseria meningitidis</i> ATCC® 13090
Maltose Test	<i>Escherichia coli</i> ATCC® 25922	<i>Pseudomonas aeruginosa</i> ATCC® 27853
Mannitol Test	<i>Staphylococcus aureus</i> ATCC® 25923	<i>Pseudomonas aeruginosa</i> ATCC® 27853
Mannose Test	<i>Salmonella typhimurium</i> ATCC® 14028	<i>Proteus mirabilis</i> ATCC® 25933
Raffinose Test	<i>Klebsiella pneumoniae</i> ATCC® 13883	<i>Proteus mirabilis</i> ATCC® 25933
Rhamnose Test	<i>Salmonella typhimurium</i> ATCC® 14028	<i>Proteus mirabilis</i> ATCC® 25933
Salicin Test	<i>Klebsiella pneumoniae</i> ATCC® 13883	<i>Escherichia coli</i> ATCC® 25922
Sorbitol Test	<i>Escherichia coli</i> ATCC® 25922	<i>Streptococcus pyogenes</i> ATCC® 19615
Sucrose Test	<i>Klebsiella pneumoniae</i> ATCC® 13883	<i>Salmonella typhimurium</i> ATCC® 14028
Threalose Test	<i>Escherichia coli</i> ATCC® 25922	<i>Pseudomonas aeruginosa</i> ATCC® 27853
Xylose Test	<i>Salmonella typhimurium</i> ATCC® 14028	<i>Staphylococcus aureus</i> ATCC® 25923

PERFORMANCE CHARACTERISTICS

Performance testing of Sugars Fermentation was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

It is recommended that biochemical, immunological, molecular testing be performed on colonies from pure culture for complete identification.

Reversal of pH may occur due usually to further breakdown of acid formed from carbohydrates, or to accumulation of ammonia from breakdown of nitrogen compounds.

Frequent observation of results is necessary since organisms do not necessarily ferment all carbohydrates at the same rate.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.

- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Adonitol Test	30 tests	88208
Arabinose Test	30 tests	88209
Arabitol Test	30 tests	88207
Dulcitol Test	30 tests	88210
Galactose Test	30 tests	88201
Glucose Test	30 tests	88202
Inositol Test	30 tests	88211
Inulin Test	30 tests	88212
Lactose Test	30 tests	88203
Levulose Test	30 tests	88213
Maltose Test	30 tests	88204
Mannitol Test	30 tests	88214
Mannose Test	30 tests	88215
Raffinose Test	30 tests	88205
Rhamnose Test	30 tests	88216
Salicin Test	30 tests	88217
Sorbitol	30 tests	88218
Sucrose	30 tests	88206
Trehalose	30 tests	88219
Xylose	30 tests	88220

Revision History

Revision	Release Date	Change Summary
10	2023-05-15	Updated layout and content in compliance with IVDR 2017/746.

This IFU document and the SDS are available from the online Support Center: liofilchem.com/ifu-sds

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References / Riferimenti

1. Reiner K (2012) Carbohydrate Fermentation Protocol. American Society for Microbiology.
2. Bailey WR, Baron EJ (1986) Bayley and Scott's Diagnostic Microbiology. 7th Ed. Mosby.
3. Lennette EH (1985) Manual of Clinical Microbiology. 4th Ed. American Society for Microbiology.
4. Sanders AC, Faber JE, Cook TM (1957) A rapid method for the characterization of enteric pathogens using paper discs. Appl. Microbiol. 5:36–40.

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