Bacitracin Test

Diagnostic discs for identification of *Streptococcus pyogenes*.

INTENDED PURPOSE

Test used for differentiating group A streptococci from other beta-haemolytic streptococci belonging to different Groups. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

Bacitracin Test are paper discs impregnated with bacitracin, a polypeptide antibiotic, and used in a diagnostic workflow for identification of Lancefield group A streptococci based on susceptibility to that antimicrobial agent.

KIT CONTENT

- 2 cartridges of 50 discs.

Cartridges are individually packed in a foil-sealed blister pack with a desiccant.

METHOD PRINCIPLE

A paper disc impregnated with bacitracin is applied directly to an inoculated agar plate before incubation. *Streptococcus pyogenes* or group A Streptococcus (GAS), which are sensitive to bacitracin, produce a zone of inhibition around the disc, whilst other beta-haemolytic streptococci either grow right up to the edge of the disc or show a very small zone of inhibition.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loop, pipettes, culture media, physiological solution, forceps, quality control organisms.

REAGENTS

- Each disc is impregnated with 0.04 IU of Bacitracin

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Allow product to reach room temperature before use for minimizing condensation on the disc.
2. Using a sterile loop or swab inoculate evenly the surface of a blood agar plate with a pure culture of the test organism exhibiting the appropriate haemolysis from the primary isolation plate.
3. Using sterile forceps, apply one disc of Bacitracin Test to the surface of the inoculated agar.
4. Invert and incubate plates at $35 \pm 2^\circ\text{C}$ for 18-24 hours in ambient air (or in atmosphere with 5-10% CO_2).
5. Examine for an inhibition zone around the disc.

INTERPRETING RESULTS

The test organism is considered sensitive (S) to 0.04 IU bacitracin and presumptively *S. pyogenes* if a zone of inhibition is formed around the disc.

No zone of inhibition (growth up to and around disc) indicates that the test organism is resistant (R) to 0.04 IU bacitracin and is not *S. pyogenes*.

STORAGE

Store Bacitracin Test at $-20^\circ\text{C}/+8^\circ\text{C}$ in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. Once opened from the desiccant-containing packaging, discs should be

stored at 2-8°C and used within 7 days. Do not use the product beyond the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

2 years.

QUALITY CONTROL

Control strain		Interpretation
<i>Streptococcus pyogenes</i>	ATCC® 19615	S, zone of inhibition (≥ 10 mm)
<i>Streptococcus agalactiae</i>	ATCC® 13813	R, no inhibition zone

PERFORMANCE CHARACTERISTICS

Performance testing of Bacitracin Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Bacitracin susceptibility is a presumptive test only. Further biochemical and/or serological tests are needed to confirm identification.

Different concentrations of bacitracin produce different zone sizes. Bacitracin test discs containing low level of bacitracin (0.04 UI) are differential tests only. Discs with a higher concentration of bacitracin (10 units) are used for antimicrobial susceptibility testing (AST).

An excessive inoculum concentration may result in the absence of a zone of inhibition (false-negative) with some group A streptococci, while a light inoculum may produce a false zone of inhibition (false-positive). Also, some streptococci may not grow in the absence of CO₂.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.
- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
Bacitracin Test	2 x 50 discs	9502

Revision History

Revision	Release Date	Change Summary
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

This IFU document and the SDS are available from the online Support Center:

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
References / Riferimenti

1. Public Health England. (2021). Identification of *Streptococcus* species, *Enterococcus* species and morphologically similar organisms. UK Standards for Microbiology Investigations. ID 4. Issue 4. <https://www.gov.uk/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
2. Ruoff KL, Whiley RA, Beighton D (1999) Streptococcus. In P.R. Murray, E.J. Baron, M.A. Tenover, F.C. Tenover, R.H. Tenover (ed) Manual of Clinical Microbiology 7th ed. American Society for Microbiology, Washington, D.C.
3. Maxted WR (1953) The use of bacitracin for identifying Group A haemolytic Streptococci. J Clin Path. 6(3):224-226.

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