X Factor Test V Factor Test V+X Factor Test

INTENDED PURPOSE

Test used for *Haemophilus* spp. differentiation. This test is intended as an aid in the diagnosis, requiring further tests to complete the diagnostic results.

DESCRIPTION

X and V Factor Test are paper discs impregnated with hemin and nicotinamide adenine dinucleotide (NAD), respectively. V+X Factor Test discs contain both hemin and NAD. The three discs are used for differentiation of *Haemophilus* species based upon their requirements for X and V growth factors.

KIT CONTENT

- 2 cartridges of 50 discs

METHOD PRINCIPLE

Species of the genus *Haemophilus* have varying requirements for X and V growth factors, which are both present in blood. The significant differences in growth factor requirements of *Haemophilus* spp. allows for their differentiation. The factors are incorporated in filter paper discs which are placed on a blood free medium previously inoculated with the organism under test. After incubation, the presence or absence of growth around the discs indicates a requirement for that particular factor.

MATERIALS REQUIRED BUT NOT PROVIDED

Standard microbiological supplies and equipment such as: inoculating loop, pipettes, swabs, culture media, quality control organisms.

REAGENTS

- Each disc of X Factor Test contains Hemin
- Each disc of V Factor Test contains Nicotinamide Adenine Dinucleotide (NAD)
- Each disc of V+X Factor Test contains Hemin and NAD

SPECIMEN

Collect specimens in sterile containers or with sterile swabs and transport to the laboratory. Process each specimen using procedures appropriate for that sample. This product is recommended for use only with pure cultures.

Refer to specific guidelines for more detailed information.

TEST PROCEDURE

1. Allow product to reach room temperature before use, for minimizing condensation on the disc.
2. Using a sterile swab, evenly inoculate the entire surface of an agar plate (Tryptic Soy Agar, Brain Heart Infusion Agar or Mueller Hinton Agar) with a pure suspension of the test organism.
3. Using aseptic technique, place one disc of X Factor Test, V Factor Test and V+X Factor Test on the inoculated agar plate, in the form of an equilateral triangle with at least 30-35 mm between discs and 1-2 cm apart from the edge of the plate.
4. Incubate at $35 \pm 2^\circ\text{C}$ for 24-48 in 5-10% carbon dioxide atmosphere.

INTERPRETING RESULTS

If the organism requires X factor only, it will grow only around the X and the X+V factors discs.

If the organism requires V factor only, it will grow only around the V and the X+V factors discs.

If both X and V factors are required, the organism will grow only around the X+V factor disc.

See below for example of test results.

Haemophilus species	Growth around disc		
	X	V	V+X
<i>H. aegyptius</i>	-	-	+
<i>H. aphrophilus</i>	-	-	+
<i>H. ducreyi</i>	+	-	+
<i>H. haemolyticus</i>	-	-	+
<i>H. influenzae</i>	-	-	+
<i>H. parahaemolyticus</i>	-	+	+
<i>H. parainfluenzae</i>	-	+	+
<i>Aggregatibacter aphrophilus</i>	-	+	+

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

SHELF LIFE

2 years.

QUALITY CONTROL

Control strain		Growth around disc		
		X	V	V+X
<i>Haemophilus influenzae</i>	ATCC® 19418	-	-	+
<i>Aggregatibacter aphrophilus</i> (formerly <i>H. parainfluenzae</i>)	ATCC® 7901	-	+	+

PERFORMANCE CHARACTERISTICS

Performance testing of X Factor Test, V Factor Test and V+X Factor Test was carried out using the QC strains listed above. The results obtained met the established criteria.

LIMITATIONS

Because similarities exist in growth factor requirements of *Haemophilus* spp., it is not recommended that this procedure be the sole criterion for species identification.

Care should be taken during inoculation of specimens onto culture media to prevent nutrient carryover.

WARNING AND PRECAUTIONS

- 1) **For *in vitro* diagnostic use (IVD).**
- 2) **For laboratory professional use only.**
- 3) Operators must be trained and have certain experience. Please read the instructions carefully before using the product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.
- 4) Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.
- 5) Do not use if the product or packaging appears to be damaged.
- 6) Follow standard precautions. All patient specimens should be considered potentially infectious and handled accordingly.
- 7) Handle all specimens as if infectious using safe laboratory procedures. Dispose of hazardous or biologically contaminated materials according to the practices of your institution.
- 8) Avoid cross-contamination of samples by using disposable tips and changing them after each sample.
- 9) Do not mix reagents of different batches. Please use the product within the validity period.
- 10) Do not eat, drink, smoke, apply cosmetics or handle contact lenses in areas where reagents and human specimens are handled.
- 11) Results should be interpreted by a trained professional in conjunction with the patient's history and clinical signs and symptoms, and epidemiological risk factors.
- 12) Ensure laboratory equipment is calibrated and maintained in accordance with the laboratory's procedure.

- 13) When test results are transmitted from the laboratory to an informatics centre, attention has to be done to avoid erroneous data transfer.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

See the references at the end of this document.

TABLE OF SYMBOLS

See the table of symbols at the end of this document.

ORDER INFORMATION

Product	Packaging	Ref.
X Factor Test	2 cartridges x 50 discs	9503
V Factor Test	2 cartridges x 50 discs	9504
V + X Factor Test	2 cartridges x 50 discs	9505

Revision History

Revision	Release Date	Change Summary
0	2023-05-15	Updated layout and content in compliance with IVDR 2017/746, version reset to revision 0

This IFU document and the SDS are available from the online Support Center:

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
References / Riferimenti

1. Public Health England. (2019). X and V Factor Test. UK Standards for Microbiology Investigations. TP 38 Issue 4. <https://www.gov.uk/guidance/uk-standards-for-microbiology-investigations-smi-quality-and-consistency-in-clinical-laboratories>
2. Kilian M. (1980) Haemophilus. In Manual of Clinical Microbiology. Eds. Lennette et al. Amer. Soc. For Microbiol. 3rd edn. Washington.
3. Hammond GW et al. (1978). Determination of the hemin requirement of Haemophilus ducreyi: evaluation of the porphyrin test and media used in the satellite growth test. J. Clin. Micro. 7(3), 243-246.

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