# CHROMATIC STAPH AUREUS/ MRSA

## FORMULA TIPICA CHROMATIC STAPH AUREUS (g/l)

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


Ref. 18007	Contenuto / Content: 20 piastre a 2 settori / plates with 2 sectors
------------	---

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## CHROMATIC STAPH AUREUS TYPICAL FORMULA (g/l)

Peptone	40.0
Sodium Chloride	25.0
Chromogenic Selective Mix	2.6
Agar	15.0
Final pH 7.0 ± 0.2	

## CHROMATIC MRSA TYPICAL FORMULA (g/l)

Peptone and Yeast Extract	40.0
Sodium Chloride	25.0
Chromogenic Selective Mix	2.5
Agar	15.0
Final pH 6.9 ± 0.2	

### DESCRIPTION

CHROMATIC STAPH AUREUS medium has a clear appearance and light amber colour, and is used for *Staphylococcus aureus* isolation. CHROMATIC MRSA medium has an opaque and cloudy appearance and light beige colour, and is used for methicillin- resistant *Staphylococcus aureus* isolation.

### METHOD PRINCIPLE

Specially selected peptones supply the nutrients. The chromogen mixtures consist of artificial substrates (chromogens), which release an insoluble coloured compound when hydrolyzed by specific enzymes. Selective agents have also been added for the suppression of gram- negative microorganisms and partial suppression of yeasts.

### TEST PROCEDURE

Inoculate the plates by streaking the specimen onto the surface of the media using a sterile loop.  
Incubate at 35±2°C for 18-24 hours.

### INTERPRETING RESULTS

Growth of mauve to orange- mauve colored colonies on both media indicates the presence of methicillin- resistant *Staphylococcus aureus*. Growth of mauve to orange- mauve colored colonies only on CHROMATIC STAPH AUREUS indicates the presence of *Staphylococcus aureus*. Most gram- positive microorganisms, if not inhibited, will produce white, blue, green or blue-green colonies. Gram- negative microorganisms and yeasts are partially to completely suppressed. Final identification must be performed through serological and biochemical tests.

### QUALITY CONTROL

- Control of general characteristics, label and print
- Sterility control  
5 days at 20 ± 2°C  
5 days at 35 ± 2°C
- Microbiological control  
Plates are inoculated with the following microbial strains: *Escherichia coli* ATCC 25922, *Staphylococcus aureus* MRSA ATCC 43300, *Proteus mirabilis* ATCC 25933, *Staphylococcus aureus* ATCC 25923.  
Inoculum for productivity: 10-100 CFU/ml  
Inoculum for specificity: ≤ 10<sup>4</sup> CFU/ml  
Inoculum for selectivity: 10<sup>4</sup>-10<sup>5</sup> CFU/ml  
Incubation conditions: 18-24 h at 35±/-2°C.

### STORAGE

6-12°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

### WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.