SWAB SOLUTION SKELETON

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Rev. 1/ 10.07.2008

SWAB SOLUTION SKELETON

Solution to hydrate samples collection swabs

DESCRIPTION

SWAB SOLUTION SKELETON is a watery sterile solution to hydrating samples swabs to use before performing microbiological surfaces checking.

CONTENT OF THE PACKAGES

Each package contains:

- 20 tubes containing 5 mL of **SWAB SOLUTION SKELETON**
- 1 instruction sheet

PRINCIPLE OF THE METHOD

SWAB SOLUTION SKELETON contains 0,1% of tryptone and 0,85% of NaCl in watery sterile solution, as described by official European Community Gazzetta (Decision 2001/471/CE 08.06.01). Beside it contains Tween 80 and lecithine which are able to neutralize surface disinfectans, particularly quaternary ammonium compounds, phenol disinfectans, hexachlorophene, formalin and ethylic alcohol.

COMPOSITION	(g/ L)
Sodio Cloruro	8.5
Triptone	1.0
Agar	1.0
Tween 80	30.0 mL
Lecitina	3.0
pH= 7.0 ± 0.2	

USE

- Take one tube from the package and aseptically introduce the swab for 10 seconds.
- Eliminate the liquid in excess from the swab lightly pushing it against the walls of the tube.
- Perform the sampling and close the swabs into its container.
- Repeat the procedure for other samples. Do not use the same tube for more than one sample.

QUALITY CONTROL

Each batch of **SWAB SOLUTION SKELETON** is submitted to the quality control using the following microorganisms:

Microorganism		Growth
<i>Escherichia coli</i>	ATCC 25922	Good
<i>Staphylococcus aureus</i>	ATCC 25923	Good
<i>Enterococcus faecalis</i>	ATCC 19433	Good
<i>Streptococcus pyogenes</i>	ATCC 19615	Good

PRECAUTIONS

SWAB SOLUTION SKELETON cannot be classified as being hazardous according to the current legislation, nor does it contain harmful substances in concentrations $\geq 1\%$. It therefore does not require a Safety Data Sheet to be available.

SWAB SOLUTION SKELETON is a disposable device to be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store **SWAB SOLUTION SKELETON** at 2-8°C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, **SWAB SOLUTION SKELETON** will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **SWAB SOLUTION SKELETON** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

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
PRESENTATION

Product	REF	Σ
SWAB SOLUTION SKELETON	20157	20

TABLE OF SYMBOLS

LOT Batch code	 Do not reuse
REF Catalogue number	 Fragile, handle with care
 Manufacturer	 Contains sufficient for <n> tests
 Use by	 Caution, consult accompanying documents
 Temperature limitation	