Mycoplasma Transport Broth

Liquid transport medium for the enrichment and detection of mycoplasma and other pathogenic microorganisms from clinical specimens.

DESCRIPTION

Mycoplasma Transport Broth is a liquid medium used for the collection and transport of clinical specimens (e.g. vaginal and urethral swabs, seminal fluid, etc.) containing *Trichomonas vaginalis*, *Gardnerella vaginalis*, *Candida* spp, mycoplasmas and ureaplasmas from the collection site to the testing laboratory. The broth may be used as enrichment medium as well.

KIT CONTENTS

Each kit contains:

- 20 tubes of Mycoplasma Transport Broth;
- 1 instruction sheet.

PRINCIPLE OF THE METHOD

Mycoplasma Transport Broth can sustain viability of fragile organisms during the period of time elapsed from specimen collection to laboratory processing. The highly nutritive formula allows microbial growth when the broth culture is incubated under appropriate conditions.

COMPOSITION (g/l)

Tryptone	17.0	Sodium Chloride	5.0
Soy Peptone	3.0	Dipotassium Phosphate	2.5
Glucose	2.5	Bovine Albumin	5.0
Final pH 7.3 ± 0.2			

PROCEDURE FOR USE

1. Allow tubes to come to room temperature before using.
2. Collect the specimen preferably from the site where the suspected organism is most likely to be found, with as little external contamination as possible, and at the clinical stage of disease most likely to yield a positive culture.
3. a) Place the specimen into the tube and send to the laboratory for prompt analysis. Mycoplasma Transport Broth will maintain viability of many microorganisms for 72 hours. It is a common practice to refrigerate specimens to minimize growth of contaminant bacteria and fungi, nonetheless certain organisms, such as *T. vaginalis*, seem to have a better recovery when kept at room temperature.
b) Inoculate the medium directly with the clinical specimen and incubate at 36 ± 1°C for 36-72 hours for cultivating mycoplasma.
4. Subculture onto appropriate agar plates.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical sheet of the medium being used.

QUALITY CONTROL

Recovery tests are performed with the following organisms: *Mycoplasma hominis* ATCC® 23114, *Ureaplasma urealyticum* ATCC® 27618, *Gardnerella vaginalis* ATCC® 14018, *Candida albicans* ATCC® 10231. Inoculated tubes are stored at room temperature for 24 hours prior to subculturing onto appropriate media.

FACTORS THAT MAY INVALIDATE THE RESULTS

Condition, timing and volume of specimen collected for culture are significant variable in obtaining reliable culture results.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *in vitro* Diagnostic Use only and must be used by properly trained operators.

STORAGE AND TRANSPORT CONDITIONS

2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the storage or transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, do not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

REFERENCES

- Lennette E.H, A. Balows, W.J. Hausler, H.J. Shadomy (1985) Manual of Clinical Microbiology 4th Ed.

PRESENTATION

Product	Ref.	Content
Mycoplasma Transport Broth	20158	20 x 6.5 ml tubes

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	 Do not reuse

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LIOFILCHEM® S.r.l.