CHRISTENSEN UREA AGAR

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Final pH = 6.8 ± 0.2 at 25°C.

DESCRIPTION

CHRISTENSEN UREA AGAR is a medium for determination of urease enzyme on culture of *Proteus* and other microorganisms, as enterococci and yeasts (*Cryptococcus*). Microorganisms which possess urease activity hydrolyze urea, with production of ammonium ions and consequent colour change of the indicator to purple-red (alkaline pH). Positive reaction of *Proteus* is observed after 6 hours of incubation on slant and on all medium within 24 hours. Urease positive strains of *Enterobacter*, *Citrobacter*, *Klebsiella* hydrolyze urea as rapidly.

TECHNIQUE

Dissolve the content of the bottle in 100°C boiling steam cooker. Check the complete dissolution by inverting the bottle gently. Cool to 45-50°C and aseptically add 10 ml of Urea 40% supplement (code 80059). Mix well avoiding bubbles formation, aseptically distribute into final tubes and allow to solidify in a slant position.

QUALITY CONTROL

Solid medium

Appearance: homogeneous.

Colour: light orange-red.

Prepared medium

Appearance: slightly opalescent.

Colour: reddish-orange.

Incubation conditions: 36±1°C / 6-24 hours.

Microorganism	ATCC	Urease production
<i>Proteus vulgaris</i>	13315	positive, red or cerise medium
<i>Escherichia coli</i>	25922	negative, no colour change

PERFORMANCE AND LIMITATIONS

The alkaline reaction produced in this medium after prolonged incubation may not be caused by urease activity. False positive reactions may occur due to the utilization of peptones or other proteins which raise the pH.

STORAGE

Store at 10-30°C, in a dry environment, in its original container tightly closed and use it before the expiry date on the label or until signs of deterioration or contamination are evident.

Store prepared tubes at 2-8°C.

REFERENCES

- Christensen, W.B. (1946). *J. Bact.* **52**: 461-466.
- Maslen, L.G.C. (1952). *Brit. Med. J.* **2**: 545-546.
- ISO 6785: 2001. IDF 93: 2001.

PACKAGING

Code 412100	
CHRISTENSEN UREA AGAR	5 x 200 ml
UREA 40% supplement.....	1 x 50 ml