## CHROMATIC STREP B

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### DESCRIPTION

**CHROMATIC STREP B** is a selective and differential medium for isolation of Group B streptococci.

### METHOD PRINCIPLE

Peptone is a product of controlled enzymatic hydrolysis of casein and together with heart extract provide peptides and free amino acids. Yeast extract is a source of amino acids and vitamins of group B. The chromogenic selective mix allows the inhibition of Gram-negative and Gram-positive microorganisms except *S. agalactiae*, and its differentiation respect other streptococci. Agar is the solidifying agent.

### PREPARATION

Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes. Allow the medium to solidify.

### TEST PROCEDURE

Inoculate the plates by streaking the specimen onto the medium surface, using a sterile loop.

Incubate at: 36±1°C for 24-48 hours.

### INTERPRETING RESULTS

*Streptococcus agalactiae* shows bright blue-green colonies, with white centre. *Enterococcus faecalis* shows mauve colonies.

### QUALITY CONTROL

1. Control of general characteristics, label and print
2. Microbiological control

Plates are inoculated with the following microbial strains *Streptococcus agalactiae* ATCC 13813, *Enterococcus faecalis* ATCC 19433, *Escherichia coli* ATCC 25922.

Inoculum for productivity: 10-100 UFC/ml

Inoculum for selectivity: 100-1000 UFC/ml

Inoculum for specificity: ≤ 10<sup>4</sup> UFC/ml

Incubation conditions: 24-48 h at 36±1°C.

### STORAGE

2-8°C away from light, until the expiry date on the label or until signs of deterioration or contamination are evident.

### WARNING AND PRECAUTIONS

The product cannot be classified as hazardous under current legislation nor does it contain harmful substances in concentrations ≥1%. It therefore does not require a Safety Data Sheet to be available. The product is designed for *In vitro* diagnostic use and must be used only by properly trained operators.

### DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

Liofilchem s.r.l.

Via Scozia-Zona industriale - 64026 Roseto degli Abruzzi Tel. +39.085.8930745 - Fax +39.085.8930330  
Web site: <http://www.liofilchem.net> E-mail: [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)

Rev. 0 /19.05.2009