HEMO-ANAEROBIC CULTURING

ITALIANO

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HEMO-ANAEROBIC CULTURING

ENGLISH

Diagnostic device for isolation of fastidious microorganisms from venous blood

DESCRIPTION

The bottle **HEMO-ANAEROBIC CULTURING** contains a very nutritive medium, Schaedler + Vit.K1 Broth, for culturing fastidious, anaerobic and facultative anaerobic, microorganisms from venous blood.

The bottle is endowed with a plastic cap red in colour, a perforable rubber subcap, and a transparent protective cap, for 80 and 40 mL bottles.

The bottle **HEMO-ANAEROBIC CULTURING Neonatal** is endowed with a perforable rubber stopper grey in colour and with a metallic screw cap with a lift plastic overcap red in colour.

For an optimum use the product is available in three presentations:

- **HEMO-ANAEROBIC CULTURING** (for adults)
- **HEMO-ANAEROBIC CULTURING Pediatric** (for children)
- **HEMO-ANAEROBIC CULTURING Neonatal** (for new-borns and children under 2 years old)

CONTENT OF THE PACKAGES

On the basis of its use, the content of each package is the following:

HEMO-ANAEROBIC CULTURING

- 6 bottles with 80 mL of Schaedler + Vit.K1 Broth + PABA + SPS
- 1 instructions sheet.

HEMO-ANAEROBIC CULTURING Pediatric

- 6 bottles with 40 mL of Schaedler + Vit.K1 Broth + PABA + SPS
- 1 instructions sheet.

HEMO-ANAEROBIC CULTURING Neonatal

- 6 bottles with 9 mL of Schaedler + Vit.K1 Broth + PABA + SPS
- 1 instructions sheet.

PRINCIPLE OF THE METHOD

Blood samples must be drawn from the patient with rigorously aseptic technique and sterile devices. The so-obtained sample is inoculated into the hemo-culturing bottle and mixed with the liquid medium whose formulation increases the growth of fastidious, anaerobic and facultative anaerobic, microorganisms from venous blood.

Tryptone, peptone and yeast extract in the broth provide amino acids, proteins and vitamins for microorganisms growth. Sodium chloride maintains the osmotic balance of the medium.

Sodium phosphate bibasic is the buffer of the medium and dextrose is a carbon source. p- Aminobenzoic Acid is a forerunner of folic acid and neutralizes the effects of sulphonamides in the inoculum through a competitive inhibition.

Sodium-polyanethol-sulphonate inhibits coagulation, neutralizes bactericide effect of human serum, prevents fagocytosis and partially inactivates some antibiotics (streptomycin, kanamycin, gentamicin and polymyxin B).

L-Cystine, hemin and Vitamin K1 constitute further growth factors for fastidious microorganisms.

CULTURE MEDIUM FORMULA (g/L)

Schaedler + Vit.K1 Broth with PABA and SPS	
Tryptone.....8.0	Potassium Phosphate Bibasic.....0.8
Soy Peptone.....1.0	L- Cystine.....0.4
Peptone.....10.0	Hemin.....0.01
Yeast Extract.....5.0	Sodium-polyanethol-sulphonate..0.3
Sodium Chloride.....1.7	p- Aminobenzoic Acid.....0.05
Dextrose5.7	Vitamin K1.....0.01
Final pH 7.6 ± 0.2	

COLLECTION AND STORAGE OF THE SAMPLE

Blood samples must be collected from the patient by puncturing the vein with strictly aseptic technique and sterile devices.

The zone of the venous puncture must be prepared disinfecting in the usual way, for example with alcohol or other disinfectants.

Do not touch the skin after the disinfection. Perform blood collection caring not to contaminate blood with the disinfectant.

Collect 8-10 mL of venous blood (4-5 mL in case of pediatric collection, 0.5- 1 mL in case of neonatal collection).

TEST PROCEDURE

a) Inoculum technique

1. Observe the broth in the bottle and check there are not signs of contamination.
2. Lift the protective plastic cap and disinfect the external side of the rubber stopper.
3. Penetrate the rubber stopper with the syringe needle and aseptically inject the blood: negative pressure in the bottle will allow the intake of the drawn blood.
4. Close the protective plastic cap and carefully mix blood with culture medium in the bottle.
5. Write patient's name and other informative data on the label and incubate at 36 ± 1 °C.

b) Subculturing technique

1. Gently mix the bottle so to suspend the culture medium.
2. Draw with a sterile syringe an amount of culture medium so to inoculate 1-2 drops on each selective plate prepared for subculture. Use suitable media for the isolation of Gram-positive and Gram-negative fastidious microorganisms. Incubate according to the different needs of the microorganisms to isolate.
3. Even if the culture medium is clear perform the subculture on selective media. Beside perform Gram staining and examine at the microscope (100X) to observe the presence and the morphology of microorganisms.
4. Perform observation at microscope and subcultures after 24, 48 and 72 hours. If there is no growth of microorganisms on selective media, prolong the incubation time of the bottle at 36 ± 1 °C until 7 days, before performing final subculture and considering negative the result.

QUALITY CONTROL

Each batch of **HEMO- ANAEROBIC CULTURING** is subjected to the quality control using the following reference microorganisms:

<i>Streptococcus mutans</i>	ATCC 35668
<i>Neisseria meningitidis</i>	ATCC 13090
<i>Clostridium sporogenes</i>	ATCC 19404
<i>Bacteroides fragilis</i>	ATCC 25285
<i>Fusobacterium nucleatum</i>	ATCC 25586
<i>Haemophilus influenzae</i>	ATCC 10211
<i>Haemophilus parainfluenzae</i>	ATCC 7901
<i>Peptostreptococcus anaerobius</i>	ATCC 27337
<i>Streptococcus pyogenes</i>	ATCC 19615
<i>Streptococcus pneumoniae</i>	ATCC 6305

LIMITS AND WARNINGS

HEMO-AEROBIC CULTURING can promote the growth of a wide variety of microorganisms.

Nevertheless in clinical use some microorganisms can fail to grow for the following reasons:

- the strain is particularly fastidious.
- the microorganism was not present in the blood sample.
- the blood sample contained an inhibiting concentration of a chemoantibiotic.
- some microorganisms do not grow in the recommended incubation time. In some cases early elimination of a culture bottle can lead to a falsely negative result.

PRECAUTIONS

The product **HEMO-ANAEROBIC CULTURING** cannot be classified as hazardous under current legislation, nor does it contain harmful substances in concentrations $\geq 1\%$; it therefore does not require a Safety Data Sheet to be available.

HEMO-ANAEROBIC CULTURING is a disposable device to be used only for diagnostic use *in vitro*. It must be used in the laboratory by properly trained personnel using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store at 10-25 °C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature.

In such conditions, the product will remain valid until the expiry date indicated on the label. Do not use beyond that date.

Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **HEMO-ANAEROBIC CULTURING** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

REFERENCES

1. Murray P.R. Manual of Clinical Microbiology 7th ed. 2005 p. 64- 67. ASM Press, Washington D.C.
2. Finegold S.M. Diagnostic Microbiology 7th ed. 1986. p.205- 224. Published C.V. Mosby Co. St. Louis

PRESENTATION

Product	Code	Package
HEMO-ANAEROBIC CULTURING	490020	6 bottles
HEMO-ANAEROBIC CULTURING <i>Pediatric</i>	490040	6 bottles
HEMO-ANAEROBIC CULTURING <i>Neonatal</i>	490060	6 bottles

TABLE OF SYMBOLS

LOT	Batch Code		Do not reuse
REF	Catalogue Number		Fragile, handle with care
	Manufacturer		Contains sufficient for <n> tests
	Use by		Caution, consult accompanying documents
	Temperature limitation	IVD	<i>In Vitro Diagnostic Medical Device</i>