LISTERIA DEMI FRASER BROTH - BAG

ITALIANO

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LISTERIA DEMI FRASER BROTH - BAG

ENGLISH

Ready to use liquid medium in bag for enrichment of *Listeria monocytogenes* in foods (ISO 11290)

TYPICAL FORMULA (g/l)

| | | | |
|--------------------|------|-------------------------------|---------|
| Tryptone | 5.0 | Sodium Phosphate Bibasic | 12.0 |
| Peptone | 5.0 | Potassium Phosphate Monobasic | 1.35 |
| Lithium Chloride | 3.0 | Aesculin | 1.0 |
| Beef Extract | 5.0 | Ammonium Ferric Citrate | 0.5 |
| Yeast Extract | 5.0 | Nalidixic Acid | 0.01 |
| Sodium Chloride | 20.0 | Acriflavine | 12.5 mg |
| Final pH 7.2 ± 0.2 | | | |

DESCRIPTION

LISTERIA DEMI FRASER BROTH – BAG is a ready to use liquid medium in bag for enrichment of *Listeria monocytogenes* in foods according to ISO 11290. The bags can be placed on special supports, connected to sterile tubes with the aid of specific joints and used with diluters, peristaltic pumps or, simply, you can bring the liquid to fall.

METHOD PRINCIPLE

LISTERIA DEMI FRASER BROTH prepared according to the formulation described by Fraser, is a modification of UVM enrichment broth. Tryptone, peptone and beef and yeast extracts provide nitrogen, carbon, sulphur and vitamins. Lithium chloride and sodium chloride maintain the osmotic balance of the medium. Potassium phosphate monobasic and sodium phosphate bibasic are buffers. High sodium chloride concentration is a selective agent. The presence of *Listeria* spp. is indicated by hydrolysis of aesculin to produce esculetin which reacts with ferric ions causing the blackening of the broth.

TECHNIQUE

Inoculate 25 g or 25 ml of the sample in 225 ml of LISTERIA DEMI FRASER BROTH. Homogenize and incubate at 30°C for 18-24 hours. Take 0.1 ml of the culture, inoculate 10 ml of LISTERIA FRASER BROTH (ref. 24131) and incubate at 36±1°C for 18-24 hours.

RESULTS INTERPRETATION

Compare each inoculated tube with an uninoculated control against a white background. Tubes that darken or turn black must be subcultured onto O.A. LISTERIA AGAR (ref. 10223), LISTERIA PALCAM AGAR (ref. 10041) or LISTERIA AESCULIN AGAR PLATE (ref. 11230). Tubes that retain the original yellow colour must be inoculated on plating media to obtain negative control free from *Listeria* spp.

QUALITY CONTROL

- Control of general characteristics, label and print.
- Sterility control:
 - 2 days at 22 ± 2°C, in aerobiosis.
 - 1 day at 36 ± 2°C, in aerobiosis.
- Microbiological control.
 - Inoculum for productivity: 10-100 CFU/ml.

Inoculum for selectivity: 10^4 - 10^5 CFU/ml.
 Inoculum for specificity: $\leq 10^4$ CFU/ml.
 Incubation conditions: 18-24 h at $35 \pm 2^\circ\text{C}$.

| Microorganism | | Growth | Aesculin reaction |
|-------------------------------|------------|---------------------|--------------------------|
| <i>Enterococcus faecalis</i> | ATCC 19433 | Partially inhibited | Negative = yellow liquid |
| <i>Escherichia coli</i> | ATCC 25922 | Partially inhibited | Negative = yellow liquid |
| <i>Listeria monocytogenes</i> | ATCC 19111 | Good | Positive = black liquid |
| <i>Listeria monocytogenes</i> | ATCC 35152 | Good | Positive = black liquid |

STORAGE AND TRANSPORT CONDITIONS

10-25°C away from light, until the expiry date on the label. Storage or transport at 2-10°C do not alter in any way the performance of the product. Eliminate if signs of deterioration or contamination are evident.

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product must be used only by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

REFERENCES

- Fraser, J.A. And W.H. Sperber (1988). J Food Prot 51: 762-765
- AFNOR (1993). Food Microbiology. Detection of *Listeria monocytogenes* – Routine method V 08-055.
- ISO 11290 (2004) Microbiology of food and animal feeding stuffs – Horizontal method for detection and enumeration of *Listeria monocytogenes*.

PACKAGING

| Product | REF |  |
|----------------------------------|--------|---------------------------------------------------------------------------------------|
| LISTERIA DEMI FRASER BROTH – BAG | 499020 | 3 bags x 3 litres |
| LISTERIA DEMI FRASER BROTH – BAG | 499025 | 3 bags x 5 litres |

TABLES OF SYMBOLS

| | | | | |
|-----------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------|
|  Batch number |  Do not reuse |  Manufacturer |  Use by |  Fragile, handle with care |
|  Catalogue number |  Temperature limits |  Contains sufficient for <n> tests |  Attention, see instructions for use | |

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
 Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net