# VAGITEST

ITALIANO

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# VAGITEST

ENGLISH

## Dip-slide for the screening of urogenital infections

### DESCRIPTION

VAGITEST is a ready-to-use system for haemophils, lactobacilli, *Gardnerella vaginalis* and pathogenic neisseria detection, with culture media solidified onto a plastic support. Such systems are universally known as 'dip-slide'. The results obtained with VAGITEST can be used as an aid in diagnosis of vaginal infections and must be used in conjunction with data on the patient's clinical history. VAGITEST is available in the 10 and 120 test versions.

### CONTENTS OF THE PACKAGES

Each VAGITEST pack contains:

- 10 ready-to-use culture systems (dip-slides) and 10 tablets for CO<sub>2</sub> generation in the 10 tests version, and 120 ready-to-use culture systems (dip-slides) and 10 tablets for CO<sub>2</sub> generation in the 120-test version.
- 1 instruction sheet.

### CONFIGURATIONS

VAGITEST is available in the configuration shown in table n°1.

Table n°1

#### VAGITEST (ref. 50020 and ref. 51020)

1. Chocolate Agar (brown)

2. Rogosa Agar (amber)

3. Thayer Martin Agar (brown)

### PRINCIPLE OF THE METHOD

VAGITEST inoculated with the clinical material (vaginal secretions, urethral secretions etc), permits a search for haemophils, lactobacilli, *Gardnerella vaginalis* and pathogenic neisseria. Chocolate Agar is a medium enriched with defibrinated cooked horse blood and VITALEX Growth SUPPLEMENT, a mixture containing vitamins, cofactors, amino acids and carbohydrates necessary for the isolation of haemophil pathogens, of *Gardnerella vaginalis*, of pathogenic neisseria and other demanding micro-organisms. Thayer Martin Agar maintains the rich medium characteristics of Chocolate Agar, but addition of the VCAT supplement makes it a selective medium for isolation of pathogenic neisseria from clinical samples. The VCAT supplement is a mixture of antibiotics: Vancomycin and Colistine which inhibit the growth of Gram-Positive and Gram-Negative contaminants, the antimycotic Amphotericine prevents fungus growth, while Trimetropim blocks the growth and swarming of proteins.

Rogosa Agar is a selective medium for the isolation of lactobacilli of oral, vaginal and faecal origin. Its selectivity is provided by a low pH and a high concentration of acetates that impede the growth of other bacteria, thus favouring that of lactobacilli.

### TYPICAL FORMULATIONS OF THE MEDIA USED IN VAGITEST (grams/litre)

CHOCOLATE AGAR		ROGOSA AGAR		THAYER MARTIN AGAR	
Peptone	15.0	Tryptose	10.0	Peptone	15.0
Starch	1.0	Yeast Extract	5.0	Starch	1.0
Dibasic Potassium Phosphate	4.0	Glucose	10.0	Dibasic Potassium Phosphate	4.0
Monobasic Potassium Phosphate	1.0	Arabinose	5.0	Monobasic Potassium Phosphate	1.0
Sodium Chloride	5.0	Sucrose	5.0	Sodium Chloride	5.0
Defibrinated Horse Blood	50.0 ml	Sodium Acetate	15.0	Defibrinated Horse Blood	50.0 ml
Vitalex Growth Supplement	10.0 ml	Ammonium Citrate	2.0	Vitalex Growth Supplement	10.0 ml
Agar	17.0	Monobasic Potassium Phosphate	6.0	Vancomycin	0.002
pH 7.2 ± 0.2		Magnesium Sulphate	0.57	Colistine	0.0075
		Manganese Sulphate	0.12	Amphotericin	0.001
		Iron Sulphate	0.03	Trimetoprim	0.003
		Agar	15.0	Agar	15.0
		Tween 80	1.0 ml	pH 7.2 ± 0.2	
		pH 5.4 ± 0.2			

### COLLECTION AND CONSERVATION OF SAMPLES

Collect the vaginal secretion with two swabs of synthetic fibre, one for the Gram colouration, the other for microbial isolation, for collection and storage of samples for gonococcus search, the AMIES transport medium is recommended, Samples for meningococcus detection must not be stored but sent directly to the laboratory.

### TEST PROCEDURE

1. Take a dip-slide from the refrigerator and leave it at ambient temperature for about 5 minutes
2. Identify the sample with an adhesive label or by writing the data directly on the dip-slide container
3. Unscrew and extract the slide from its cylindrical container without touching the surface of the culture media
4. Holding the slide by the stopper, seed the surfaces of the culture with the clinical material (vaginal secretion, urethral secretion etc) taken with a suitable collection swab
5. Introduce a CO<sub>2</sub> tablet moistened with a drop of sterile water into the slide cup
6. Put the slide back into the container and screw it shut
7. Incubate at 36 ± 1°C for 24-48 ore.

### INTERPRETATION OF THE RESULTS

Following incubation, the presence of bacteria is revealed by the presence of colonies on the surfaces of the culture media.

The culture media are used to estimate the quality and quantity of the various types of microbes that support the infection or constitute the normal vaginal flora, such as lactobacilli. Gonococcus grows on Chocolate Agar and on Thayer Martin Agar in colonies of 0.5 to 1 mm, white or grey in colour, opaque and raised. Presumed colonies of gonococcus are identified by Gram colouration, oxidase and sugar fermentation tests. Meningococcal colonies are larger than those of gonococci, and round and convex in form. *Gardnerella vaginalis* grows on Chocolate Agar in greyish semi-opaque colonies. The haemophils grow on Chocolate Agar in greyish, semi-opaque colonies. Lactobacilli grow on Rogosa Agar in white-grey colonies with regular margins which may be flattened or raised, smooth or rough. Not all lactobacilli grow on this medium.

**QUALITY CONTROL**

For control of the quality of VAGITEST, the dip-slides are inoculated with the strains of microbes listed in table 2, diluted in sterile physiological solution so as to obtain:

- Dimension of the inoculum for productivity: 10-100 CFU/ml
- Dimension of the inoculum for selectivity:  $10^4$ - $10^5$  CFU/ml
- Dimension of the inoculum for specificity:  $\leq 10^4$  CFU/ml
- Conditions for incubation: 24-48 hours at  $36\pm 1^\circ\text{C}$  in a  $\text{CO}_2$  enriched atmosphere (tablet supplied with the pack).

The results are interpreted as shown in Table n°2. The appearance of the media is also checked; for the colour of the media see Table n°1 of the Configuration section.

Table n°2: Strains used and theoretical results in the quality control.

Microorganism	CHOCOLATE AGAR	ROGOSA AGAR	THAYER MARTIN AGAR
<i>Escherichia coli</i> ATCC 25922	growth	growth inhibited	growth inhibited
<i>Enterococcus faecalis</i> ATCC 19433	growth	growth inhibited	growth inhibited
<i>Staphylococcus aureus</i> ATCC 25923	growth	growth inhibited	growth inhibited
<i>Gardnerella vaginalis</i> ATCC 14018	growth	growth inhibited	growth inhibited
<i>Neisseria gonorrhoeae</i> ATCC 43069	growth	growth inhibited	growth
<i>Lactobacillus casei</i> ATCC 393	growth inhibited	growth	growth inhibited

**LIMITS AND NOTES**

For collection of the clinical material for examination, use swabs of synthetic fibre, not cotton ones. All suspect colonies cultivated on VAGITEST media must be subjected to identification.

*N. lactamica* also grows on the Thayer Martin Agar medium. It has been reported that 3-10% of gonococci (normally AHU strains) are sensitive to Vancomycin and do not grow on Thayer Martin Agar (Mirrett *et al.*, 1981). These strains can be recovered on Chocolate Agar. Prolonged incubation (over 24 hours) may give rise to excessive growth of concomitant bacterial flora on the Thayer Martin Agar medium, and this may mask the growth of the gonococcus. Do not use the  $\text{CO}_2$  generation tablets if they are broken or damp, or if the pack is not intact.

**PRECAUTIONS**

The media described here are not classified as hazardous by current legislation, nor do they contain harmful substances in concentrations  $\geq 1\%$ . No safety data sheets are therefore required for them.

VAGITEST are disposable devices to be used only for diagnostic purposes *in vitro*. VAGITEST is intended for use in a professional environment and must be used in a laboratory by adequately trained personnel using approved asepsis and safety methods for dealing with pathogenic agents.

**DIP-SLIDE CONSERVATION**

Store at 2-8°C in the original packaging. Keep away from sources of heat and avoid excessive changes of temperature. In such conditions, the products described here will be valid until the expiry date shown on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration (conspicuous changes in the colour of the media, drying or detachment of the medium from its support, presence of fungal or bacterial colonies). DO NOT FREEZE.

**DISPOSAL OF USED MATERIAL**

After use, VAGITEST and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

**BIBLIOGRAPHY**

- Kellogg D.S., Holmes, K.K., Hill G.A. (1976). Laboratory Diagnosis of Gonorrhoea. Cumitech 4, American Society for Microbiology. Washington, D.C.
- Martin, J.E. and Lester, A. (1971) HSMHA Health Reports. **86**, 30-33.
- Martin, J.E., Armstrong, J.H., Smith, P.B. (1974). App. Microbiology **27**, 802-805.
- Mickelsen P. *et al.*, (1977) J. Clin. Microbiol, **5**, 488.
- Mirrett, S.L., Reller, L.B., Knapp, J.S. (1981) J. Clin. Microbiol. **14**, 94-99.
- Thayer, J.D. and Martin, J.E. (1966). Public Health Reports. **81**, 559-562.

**PRESENTATION**

Product	REF	
VAGITEST	50020	10 Dip Slide
VAGITEST	51020	120 Dip Slide

**TABLE OF SYMBOLS**

 In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	 Batch code

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**LIOFILCHEM® S.r.l.**