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
# INTEGRAL SYSTEM ENTEROBATTERI

*Colorimetric system for biochemical identification and susceptibility testing of enterobacteria*

## DESCRIPTION

INTEGRAL SYSTEM ENTEROBATTERI is a 24-well system containing desiccated biochemical and antibiotic substrates for the biochemical identification and susceptibility testing of enterobacteria.

The system is inoculated with the bacterial suspension of the micro-organism under examination and incubated at  $36\text{ }^{\circ}\text{C} \pm 1^{\circ}\text{C}$  for 18-24 hours. The tests for identification and for the susceptibility testing are interpreted by assessing the change in colour of the various wells.

## KIT CONTENTS

The kit contains:

- 20 systems of INTEGRAL SYSTEM ENTEROBATTERI
- 20 Vials of Suspension Medium (7.0 ml/vial)
- 20 Oxidase Disc
- 1 Instruction sheet

## ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

- Vaseline oil for microbiology use (Vaselin oil 2 bottles of 50 ml) .....cod. 80278
- Physiological Solution (20 vials of 7.0 ml/vial).....cod. 20095
- Kovac's Reagent (2 bottles of 25 ml).....cod. 80270
- Sundry microbiology laboratory materials

## CONFIGURATION

The configuration of the system is shown in Table 1.

Table 1

Well	BIOCHEMICAL IDENTIFICATION
<b>1-LDC</b>	Decarboxylation of lysine
<b>2-ODC</b>	Decarboxylation of ornithine
<b>3-H<sub>2</sub>S</b>	Production of hydrogen sulphide
<b>4-IND</b>	Indole test
<b>5-LAC</b>	Fermentation of lactose
<b>6-DUL</b>	Fermentation of dulcitol
<b>7-PA</b>	Deamination of phenylalanine
<b>8-UR</b>	Hydrolysis of urea
<b>9-CIT</b>	Utilisation of citrate
<b>10-OX</b>	Oxidase test
Well	SUSCEPTIBILITY TESTING (*)
<b>11-AK</b>	Amikacin - 32 µg/ml
<b>12-CN</b>	Gentamycin - 8 µg/ml
<b>13-TOB</b>	Tobramycin - 8 µg/ml
<b>14-TZP</b>	Piperacillin+Tazobactam - 128/4 µg/ml
<b>15-FOS</b>	Fosfomicin - 200 µg/ml
<b>16-CFP</b>	Cefoperazone - 64 µg/ml
<b>17-CTX</b>	Cefotaxime - 64 µg/ml
<b>18-CAZ</b>	Ceftazidime - 32 µg/ml
<b>19-AMS</b>	Ampicillin + Sulbactam – 32/16 µg/ml
<b>20-NA</b>	Nalidixic acid - 32 µg/ml
<b>21-CIP</b>	Ciprofloxacin - 4 µg/ml
<b>22-LEV</b>	Levofloxacin - 8 µg/ml
<b>23-SXT</b>	Co-Trimoxazole - 8 µg/ml
<b>24-C</b>	Growth control for the susceptibility testing

(\*): The concentration of each antibiotic is according to the updated NCCLS-January 2004, Vol 24 N°1.<sup>(1)</sup>

## PRINCIPLE OF THE METHOD

INTEGRAL SYSTEM ENTEROBACTERI makes possible simultaneous biochemical identification and assessment of the susceptibility testing of enterobacteria isolated from clinical samples.

- The identification is based on biochemical tests performed on culture media containing specific substrates in the wells from **1-LDC** to **10-OX**.
- The susceptibility testing is assessed on the basis of the growth or inhibition of the micro-organisms in the media containing the antibiotic and a growth indicator in the wells from **11-AK** to **23-SXT**. Well **24-C** does not contain antibiotics but only the culture medium and the indicator; it serves as a microbe growth control for assessment of the susceptibility testing.

## COMPOSITION

Table n°2

<i>Well</i>	<i>Content</i>
<b>1-LDC</b>	Culture medium for showing lysine decarboxylase
<b>2-ODC</b>	Culture medium for showing ornithine decarboxylase
<b>3-H<sub>2</sub>S</b>	Culture medium for showing production of hydrogen sulphide
<b>4-IND</b>	Culture medium for showing production of indole
<b>5-LAC</b>	Culture medium for showing fermentation of lactose
<b>6-DUL</b>	Culture medium for showing fermentation of dulcitol
<b>7-PA</b>	Culture medium for showing phenylalanine deamination
<b>8-UR</b>	Culture medium for showing hydrolysis of urea
<b>9-CIT</b>	Culture medium for showing utilisation of citrate
<b>10-OX</b>	Culture medium for showing enzyme cytochrome-oxidase
<b>11-AK</b>	Culture medium containing Amikacin - 32 µg/ml
<b>12-CN</b>	Culture medium containing Gentamycin - 8 µg/ml
<b>13-TOB</b>	Culture medium containing Tobramycin - 8 µg/ml
<b>14-TZP</b>	Culture medium containing Piperacillin+Tazobactam - 128/4 µg/ml
<b>15-FOS</b>	Culture medium containing Fosfomycin - 200 µg/ml
<b>16-CFP</b>	Culture medium containing Cefoperazone - 64 µg/ml
<b>17-CTX</b>	Culture medium containing Cefotaxime - 64 µg/ml
<b>18-CAZ</b>	Culture medium containing Ceftazidime - 32 µg/ml
<b>19-AMS</b>	Culture medium containing Ampicillin + Sulbactam - 32/16 µg/ml
<b>20-NA</b>	Culture medium containing Nalidixic acid - 32 µg/ml
<b>21-CIP</b>	Culture medium containing Ciprofloxacin - 4 µg/ml
<b>22-LEV</b>	Culture medium containing Levofloxacin - 8 µg/ml
<b>23-SXT</b>	Culture medium containing Co-Trimoxazole - 8 µg/ml
<b>24-C</b>	Medium without antibiotics

**Suspension Medium (g/l):** Yeast extract 5 g; Beef peptone 3 g; Glucose 2 g; Distilled water 1000.0 ml;  
pH 6.8±0.2

## COLLECTION AND STORAGE OF THE SAMPLES

The colonies to be subjected to the biochemical identification test and susceptibility testing assessment with INTEGRAL SYSTEM ENTEROBACTERI must be taken from a culture medium – selective or non selective – used for the isolation of enterobacteria.

## TEST PROCEDURE

1. Take a system from the kit.
2. After checking that the colonies that have developed on the culture media presumably belong to the enterobacteria group, take one or more morphologically similar, well isolated, colonies from a solid culture medium and suspend them in 5 ml of Physiological Solution for microbiological use in such a way as to obtain a turbidity equivalent to 0.5 McFarland (**Bacterial suspension**).
3. Trasfer
  - 0.2 ml of **Bacterial suspension** into each of the first 10 wells and to wells **1-LDC**, **2-ODC**, **3-H<sub>2</sub>S** and **8-UR** add 2 drops of Vaseline Oil for microbiological use (BIOCHEMICAL IDENTIFICATION).
  - 0.01 ml of **Bacterial suspension** to the vial of **Suspension Medium** contained into the kit and distribute 0.2 ml into the wells from **11-AK** to **24-C** (SUSCEPTIBILITY TESTING).
4. Cover the system with the lid provided and incubate at 36 °C ± 1°C for 18-24 hours.
5. After incubation. add:
  - 2 drops of Kovac's Reagent to well **4-IND** (INDOLE TEST).
  - 1 OXIDASE DISC to well **10-OX** (OXIDASE TEST).

## INTERPRETATION OF THE RESULTS

### BIOCHEMICAL IDENTIFICATION

Interpret the results of the first 10 wells using Table n° 3 and form the 4-digit code following the instructions given to the paragraph **FORMATION OF NUMERICAL CODE**. Then use the code and Table n°5 to identify the bacteria.

### SUSCEPTIBILITY TESTING

Observe the changes of colour of the wells from **11-AK** to **23-SXT** and interpret the results using Table n°3. The control well (**24-C**) must be positive (yellow). If it is negative (blue or grey), it is necessary to check the vitality of the inoculum and whether it has been correctly prepared; then to repeat the tests using a new system.

Table n°3

Well	BIOCHEMICAL IDENTIFICATION	Well colour	
		Positive reaction	Negative reaction
<b>1-LDC</b>	Decarboxylation of lysine	Violet	Yellow-Brown
<b>2-ODC</b>	Decarboxylation of ornithine	Violet	Yellow-Brown
<b>3 -H<sub>2</sub>S</b>	Production of hydrogen sulphide	Black	Yellow
<b>4-IND</b>	Indole test	Pink-Red ring	Yellow
<b>5-LAC</b>	Fermentation of lactose	Yellow	Blue-Green
<b>6-DUL</b>	Fermentation of dulcitol	Yellow	Blue-Green
<b>7-PA</b>	Deamination of phenylalanine	Black-Brown	Yellow
<b>8-UR</b>	Hydrolysis of urea	Red-Fuchsia	Yellow-Orange
<b>9-CIT</b>	Utilisation of citrate	Blue-Dark Green	Green
<b>10-OX</b>	Oxidase test	Blue-Purple (immediate reaction)	Colourless

SUSCEPTIBILITY TESTING		
WELL COLOUR	BACTERIAL GROWTH	INTERPRETATION
Blue	Inhibited	S = Sensitive
Grey	Intermediate	I = Intermediate sensitivity
Yellow	Good	R = Resistant

## FORMATION OF THE NUMERICAL CODE

1) The 10 biochemical tests are divided into 3 groups each containing 3 tests and one group containing 1 test, each one is indicated with a positiveness value of 1,2,4.

- Value 1 : first test positive in each group  
(LDC, IND, PA,OX)
- Value 2 : second test positive in each group  
(ODC,LAC,UR)
- Value 4 : third test positive in each group  
(H<sub>2</sub>S, DUL,CIT)
- Value zero: negative reaction in every group

2) Adding the number of positive reactions in each group, one obtains a 4-digit code which is used to identify the micro-organism from the table of numerical codes as in the example.

	Group I			Group II			Group III			Group IV
Well	LDC	ODC	H <sub>2</sub> S	IND	LAC	DUL	PA	UR	CIT	OX
Positiveness code	1	2	4	1	2	4	1	2	4	1
Results	+	+	-	+	+	-	-	-	-	-
Codes addition	3			3			0			0
NUMERICAL CODE: 3300				IDENTIFICATION: <i>E.coli</i>						

## TABLE OF BIOCHEMICAL REACTIONS

Table n°4

Micro-organisms	Group I			Group II			Group III			Group IV
	LDC	ODC	H <sub>2</sub> S	IND	LAC	DUL	PA	UR	CIT	OX
<i>E.coli</i>	V	V	-	+	+	V	-	-	-	-
<i>E.coli inactive</i>	V	V	-	+	-	V	-	-	-	-
<i>Shigella spp.</i>	-	V	-	V	-	-	-	-	-	-
<i>Edwardsiella spp.</i>	+	+	+	+	-	-	-	-	-	-
<i>Citrobacter spp.</i>	-	V	+	V	V	V	-	V	V	-
<i>Salmonella spp.</i>	+	+	V	-	-	+	-	-	V	-
<i>Arizona spp.</i>	+	+	+	V	V	-	-	-	+	-
<i>K. pneumoniae</i>	+	-	-	-	+	V	-	V	V	-
<i>Klebsiella oxytoca</i>	+	-	-	+	+	V	-	+	V	-
<i>Enterobacter cloacae</i>	V	+	-	-	V	-	-	V	+	-
<i>E. aerogenes</i>	+	+	-	-	+	-	-	-	+	-
<i>Enterobacter hafnia</i>	+	+	-	-	V	-	-	-	+	-
<i>Serratia spp.</i>	+	+	-	-	-	-	-	V	+	-
<i>Proteus mirabilis</i>	-	+	+	-	-	-	+	+	V	-
<i>Proteus vulgaris</i>	-	-	V	+	-	-	+	+	V	-
<i>Morganella morganii</i>	-	+	-	+	-	-	+	+	-	-
<i>Providencia stuartii</i>	-	-	-	+	-	-	+	-	+	-
<i>Providencia rettgeri</i>	-	-	-	+	-	-	+	+	+	-
<i>Pseudomonas spp.</i>	+	+	-	-	-	-	-	-	+	+

+ : Positive Reaction

- : Negative Reaction

V: Variable Reaction

\*Dulcitol: Some strains of *Salmonella* ferment dulcitol slowly; they may therefore be identified initially as *Arizona spp.* A confirmation test is thus required using a polyvalent antiserum for *Salmonella*.

## TABLE OF NUMERICAL CODES

Table n°5

Code	Micro-organism	Code	Micro-organism	Code	Micro-organism	Code	Micro-organism
0000	<i>Shigella</i> spp.	2240	<i>Enterobacter cloacae</i>	4220	<i>Citrobacter</i> spp.	6160	<i>Citrobacter</i> spp.
0100	<i>E.coli</i> inactive 45%	2260	<i>Enterobacter cloacae</i>	4240	<i>Citrobacter</i> spp.	6200	<i>Citrobacter</i> spp.
	<i>Shigella</i> spp. 55%	2300	<i>E.coli</i>	4260	<i>Citrobacter</i> spp.	6220	<i>Citrobacter</i> spp.
0130	<i>Proteus vulgaris</i>	2500	<i>E.coli</i> inactive	4300	<i>Citrobacter</i> spp.	6240	<i>Citrobacter</i> spp.
0150	<i>Providencia stuartii</i>	2700	<i>E.coli</i>	4320	<i>Citrobacter</i> spp.	6260	<i>Citrobacter</i> spp.
0170	<i>Proteus vulgaris</i> 0.5%	3040	<i>Enterobacter hafniae</i> 60%	4340	<i>Citrobacter</i> spp.	6300	<i>Citrobacter</i> spp.
	<i>Providencia rettgeri</i> 99.5%		<i>Serratia</i> spp. 37%	4360	<i>Citrobacter freundii</i>	6320	<i>Citrobacter</i> spp.
0300	<i>E.coli</i>		<i>Enterobacter cloacae</i> 3%	4400	<i>Citrobacter</i> spp.	6340	<i>Citrobacter</i> spp.
0500	<i>E.coli</i> inactive	3041	<i>Pseudomonas</i> spp.	4420	<i>Citrobacter</i> spp.	6360	<i>Citrobacter</i> spp.
0700	<i>E.coli</i>	3060	<i>Serratia</i> spp. 98.9 %	4440	<i>Citrobacter</i> spp.	6400	<i>Citrobacter</i> spp.
1100	<i>E.coli</i> inactive		<i>Enterobacter cloacae</i> 1.1%	4460	<i>Citrobacter</i> spp.	6420	<i>Citrobacter</i> spp.
1200	<i>Klebsiella pneumoniae</i>	3100	<i>E.coli</i> inactive	4500	<i>Citrobacter</i> spp.	6440	<i>Citrobacter</i> spp.
1220	<i>Klebsiella pneumoniae</i>	3240	<i>Enterobacter cloacae</i> 7.7%	4520	<i>Citrobacter</i> spp.	6460	<i>Citrobacter</i> spp.
1240	<i>Klebsiella pneumoniae</i>		<i>Enterobacter aerogenes</i> 75%	4540	<i>Citrobacter</i> spp.	6500	<i>Citrobacter</i> spp.
1260	<i>Klebsiella pneumoniae</i>		<i>Enterobacter hafniae</i> 17.3%	4560	<i>Citrobacter</i> spp.	6520	<i>Citrobacter</i> spp.
1300	<i>E.coli</i>	3260	<i>Enterobacter cloacae</i>	4600	<i>Citrobacter</i> spp.	6540	<i>Citrobacter</i> spp.
1320	<i>Klebsiella oxytoca</i>	3300	<i>E.coli</i>	4620	<i>Citrobacter</i> spp.	6560	<i>Citrobacter</i> spp.
1360	<i>Klebsiella oxytoca</i>	3400	<i>Salmonella</i> spp.	4640	<i>Citrobacter</i> spp.	6600	<i>Citrobacter</i> spp.
1500	<i>E.coli</i> inactive	3440	<i>Salmonella</i> spp.	4660	<i>Citrobacter</i> spp.	6620	<i>Citrobacter</i> spp.
1600	<i>Klebsiella pneumoniae</i>	3500	<i>E.coli</i> inactive	4700	<i>Citrobacter</i> spp.	6640	<i>Citrobacter</i> spp.
1620	<i>Klebsiella pneumoniae</i>	3700	<i>E.coli</i>	4720	<i>Citrobacter</i> spp.	6660	<i>Citrobacter</i> spp.
1640	<i>Klebsiella pneumoniae</i>	4000	<i>Citrobacter</i> spp.	4740	<i>Citrobacter</i> spp.	6700	<i>Citrobacter</i> spp.
1660	<i>Klebsiella pneumoniae</i>	4020	<i>Citrobacter</i> spp.	4760	<i>Citrobacter</i> spp.	6720	<i>Citrobacter</i> spp.
1700	<i>E.coli</i>	4040	<i>Citrobacter</i> spp.	6000	<i>Citrobacter</i> spp.	6740	<i>Citrobacter</i> spp.
1720	<i>Klebsiella oxytoca</i>	4060	<i>Citrobacter</i> spp.	6020	<i>Citrobacter</i> spp.	6760	<i>Citrobacter</i> spp.
1760	<i>Klebsiella oxytoca</i>	4100	<i>Citrobacter</i> spp.	6030	<i>Proteus mirabilis</i>	7040	<i>Arizona</i> spp.
2000	<i>Shigella</i> spp.	4120	<i>Citrobacter</i> spp.	6040	<i>Citrobacter</i> spp.	7100	<i>Edwardsiella</i> spp.
2040	<i>Enterobacter cloacae</i>	4130	<i>Proteus vulgaris</i>	6060	<i>Citrobacter</i> spp.	7140	<i>Arizona</i> spp.
2060	<i>Enterobacter cloacae</i>	4140	<i>Citrobacter</i> spp.	6070	<i>Proteus mirabilis</i>	7240	<i>Arizona</i> spp.
2100	<i>E.coli</i> inactive 44.6%	4160	<i>Citrobacter</i> spp.	6100	<i>Citrobacter</i> spp.	7340	<i>Arizona</i> spp.
	<i>Shigella</i> spp. 55.4%	4170	<i>Proteus vulgaris</i>	6120	<i>Citrobacter</i> spp.	7400	<i>Salmonella</i> spp.
2130	<i>Morganella morganii</i>	4200	<i>Citrobacter</i> spp.	6140	<i>Citrobacter</i> spp.	7440	<i>Salmonella</i> spp.

## QUALITY CONTROL

Each batch of INTEGRAL SYSTEM ENTEROBATTERI is subjected to quality control using the following ATCC reference strains.

<i>Enterobacter cloacae</i> .....	ATCC 13047	<i>Proteus mirabilis</i> .....	ATCC 25933
<i>Escherichia coli</i> .....	ATCC 25922	<i>Pseudomonas aeruginosa</i> .....	ATCC 27853
<i>Klebsiella pneumoniae</i> .....	ATCC 13883	<i>Salmonella typhimurium</i> .....	ATCC 14028

## FACTORS THAT MAY INVALIDATE THE RESULTS

Poor standardisation of the inoculum; application of the method to micro-organisms not in the Enterobacteria group; mixed or contaminated cultures; use of expired systems or expired supplementary reagents; incorrect application of the technique.

## LIMITS AND WARNINGS

The identification of *Salmonella spp.* and of *Shigella spp.* must be confirmed using appropriate antisera. In some cases, supplementary confirmation tests may be required for definitive identification of the micro-organisms.

## PERFORMANCE

The bacterial identification results obtained with the INTEGRAL SYSTEM ENTEROBATTERI agree with those obtained using traditional methods as indicated by Piccolomini et al.<sup>(4)</sup>.

The antibiogram results obtained with the INTEGRAL SYSTEM ENTEROBATTERI agree with those obtained using the radial diffusion methods according to Bauer et al.<sup>(5)</sup> recommended by the Food and Drug Administration (FDA) and by the National Committee for Clinical Laboratory Standards (NCCLS) U.S.A <sup>(1,8)</sup>.

## PRECAUTIONS

The product, INTEGRAL SYSTEM ENTEROBATTERI, cannot be classified as hazardous under current legislation nor does it contain harmful substances in concentrations  $\geq 1\%$ . It therefore does not require a Safety Data Sheet to be available. INTEGRAL SYSTEM ENTEROBATTERI is a disposable device to be used only for diagnostic use *in vitro*; it is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

## STORAGE

Store at 2-8°C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.



## DISPOSAL OF USED MATERIAL

After use, INTEGRAL SYSTEM ENTEROBATTERI and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

## PRESENTATION

Product	Code	Kit
INTEGRAL SYSTEM ENTEROBATTERI	71714	20 tests

## TABLE OF SYMBOLS

SYMBOL	MEANING	SYMBOL	MEANING
	<i>In Vitro</i> Diagnostic Medical Device		Do not reuse
	Manufacturer		Contains sufficient for <n> tests
	Catalogue number		Fragile, handle with care
	Use by		Caution, consult accompanying documents
	Temperature limitation		Batch code

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