Mycoplasma System VET

System for detection, count and susceptibility testing
of mycoplasmas in the veterinary field.

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DESCRIPTION

Mycoplasma System VET is a 24-well system containing desiccated biochemical and antibiotic substrates for detection, semi-quantitative determination and susceptibility testing of *Mycoplasma* spp. and *Ureaplasma* spp.

NOTE 1: Only *Mycoplasma* species that hydrolyze arginine can be detected using this system.

NOTE 2: Some of the main species of *Mycoplasma* that affect farm animals and pets and that can be detected with the system are shown in Table 1 at the end of this document.

The system is inoculated with the sample suspension and incubated at $36 \pm 1^\circ\text{C}$ for 24-48 hours.

The tests results are interpreted by assessing the change in color of the various wells and performing microscope examination.

CONTENT OF THE KIT

- 20 Mycoplasma System VET
- 20 Vials of Physiological Solution (7.0 mL/vial)
- Instructions Sheet (includes Test Results Form), also available from liofilchem.com/ifu-sds

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

- Vaseline oil (ref. 87006)
- Sundry microbiology laboratory materials
- Cover slides and slides
- Microscope

CONFIGURATION

Well	Count and identification of mycoplasmas / ureaplasmas
1-Uu 10 ³	<i>Ureaplasma</i> spp. (titre = 10 ³ CFU/mL)
2-Uu 10 ⁴	<i>Ureaplasma</i> spp. (titre = 10 ⁴ CFU/mL)
3-Uu ≥ 10 ⁵	<i>Ureaplasma</i> spp. (titre ≥ 10 ⁵ CFU/mL)
4-Mh 10 ⁴	<i>Mycoplasma</i> spp. (titre = 10 ⁴ CFU/mL)
5-Mh ≥ 10 ⁵	<i>Mycoplasma</i> spp. (titre ≥ 10 ⁵ CFU/mL)
Well	Detection of <i>T. vaginalis</i> and <i>Candida</i> spp.
6-TR/YE	<i>Trichomonas vaginalis</i> / <i>Candida</i> spp.

Susceptibility testing				
Wells		Antibiotic	Concentrations (µg/mL)	
7-TE	8-TE	Tetracycline	4	8
9-ENR	10-ENR	Enrofloxacin	5	10
11-MAR	12-MAR	Marbofloxacin	10	20
13-DO	14-DO	Doxycycline	4	8
15-E	16-E	Erythromycin	8	16
17-CLA	18-CLA	Clarithromycin	8	16
19-TY	20-TY	Tylosin	30	60
21-CD	22-CD	Clindamycin	4	8
23-AZM	24-AZM	Azithromycin	4	8

PRINCIPLE OF THE METHOD

Mycoplasma System VET allows the detection, semi-quantitative count, presumptive identification and susceptibility testing of Mycoplasmataceae.

The system also provides the detection of *Trichomonas vaginalis* and *Candida* spp.

- The presence and semi-quantitative count of *Ureaplasma* spp. is indicated by a color change from yellow to red in wells **1-Uu 10³**, **2-Uu 10⁴** and **3-Uu $\geq 10^5$**
- The presence and semi-quantitative count of *Mycoplasma* spp. is evidenced by a color change from yellow to red in wells **4-Mh 10⁴** and **5-Mh $\geq 10^5$**
- The presence of *Trichomonas vaginalis* and *Candida* spp. is determined by microscope (40x) examination of a drop of culture liquid taken from well **6-TR/YE** to assess the presence of mobile ciliated trophozoites for identification of *Trichomonas vaginalis* and/or the presence of chlamyospores and hyphae for identification of *Candida* spp.
- The susceptibility testing of urogenital mycoplasmas is evaluated by observing a color change from yellow to red in wells **7-TE** to **24-AZM** and following the instructions under INTERPRETATION OF THE RESULTS.

COLLECTION AND CONSERVATION OF THE SAMPLES

The samples must be sent for inoculation in the Mycoplasma System VET immediately after they are taken. Do not store in the refrigerator for even limited periods of time, since low temperatures may damage the vitality of some particularly sensitive microorganisms such as *Trichomonas vaginalis* prejudicing the final result.

TEST PROCEDURE

PREPARATION OF THE SAMPLE

Carefully open a vial of Physiological Solution* using the ampoule cracker provided with the kit.

1. SWAB

Immerse the swab (after obtaining the clinical material) in the vial of physiological solution and wait 5 minutes. Carefully squeeze the swab against the vial wall so that the specimen is dispersed homogeneously in the Physiological Solution.

NOTE: Keep the swab, dipped in a nutritive broth, until the conclusion of the test.

2. BIOLOGICAL FLUIDS

Add 0.2 mL of sample to the vial of Physiological Solution, shake gently and wait 5 minutes before inoculating.

3. SAMPLE MATERIAL FROM MYCOPLASMA TRANSPORT BROTH (ref. 20158)

Add 1 mL of Mycoplasma Transport Broth, containing the sample under examination, to the vial of Physiological Solution, shake gently and wait 5 minutes before inoculating.

NOTE: Conserve the specimen in the vial of Mycoplasma Transport Broth until the conclusion of the test.

*Physiological Solution (g/L): Sodium chloride 9 g; Distilled water 1000 mL; pH 6.8 \pm 0.2

INOCULATION OF THE SYSTEM

1. Take a system from its wrap and bring it to room temperature.
2. Write down sample info and the date of the start of the examination.
3. Transfer 0.2 mL of suspension of the sample into each well of the system.
4. Overlay with 1 drop vaseline oil for microbiological use all the wells except for well 6-TR/YE.
5. Cover the system with the lid provided and incubate at 36 \pm 1°C for 24 hours, and up to 48 hours.
6. Observe and confirm the presence of mycoplasmas within 48 hours incubation.

INTERPRETATION OF THE RESULTS**Count and identification of *Mycoplasmas* / *Ureaplasmas***

Watch for the color change in the wells **1-Uu 10³** to **5-Mh ≥ 10⁵** and interpret the results as shown below.

Detection of *Trichomonas vaginalis* and *Candida* spp.

Take a drop of liquid from the well, deposit it on a glass slide, place a cover slip on top and examine at the microscope (40x) for the presence of *T. vaginalis* and *Candida* spp. Interpret as shown below.

Antimicrobial Susceptibility Testing of *Mycoplasmas* / *Ureaplasmas*

Watch for the color change in the wells **7-TE** to **24-AZM** and interpret the results as shown below. For each antibiotic, the results of two wells are evaluated together. Refer to wells **1-Uu 10³** and **4-Mh 10⁴** as growth control for the antibiogram of *Ureaplasma* spp and *Mycoplasma* spp, respectively.

Note the results on the TEST RESULTS FORM (copy as many forms as necessary).

Well	COUNT OF MYCOPLASMAS AND UREAPLASMAS	Well color	
		Positive reaction*	Negative reaction
1-Uu 10³	<i>Ureaplasma</i> spp. (titre = 10 ³ CFU/mL) ⁽¹⁾	red	yellow
2-Uu 10⁴	<i>Ureaplasma</i> spp. (titre = 10 ⁴ CFU/mL) ⁽²⁾	red	yellow
3-Uu ≥ 10⁵	<i>Ureaplasma</i> spp. (titre ≥ 10 ⁵ CFU/mL) ⁽³⁾	red	yellow
4-Mh 10⁴	<i>Mycoplasma</i> spp. (titre = 10 ⁴ CFU/mL)	red	yellow
5-Mh ≥ 10⁵	<i>Mycoplasma</i> spp. (titre ≥ 10 ⁵ CFU/mL)	red	yellow

CFU: Colony Forming Units

(1) = Equivalent to 5-20 colonies grown on Mycoplasma Agar A7

(2) = Equivalent to 20-50 colonies grown on Mycoplasma Agar A7

(3) = Over 50 colonies grown on Mycoplasma Agar A7

*An orangey colouring should be considered as test positive.

Well	T.VAGINALIS AND CANDIDA spp.	Microscope examination (40x)
6-TR/YE	<i>Trichomonas vaginalis</i> / <i>Candida</i> spp.	<i>T.vaginalis</i> : mobile ciliated protozoites <i>Candida</i> spp.: chlamydo spores and hyphae

SUSCEPTIBILITY TESTING OF MYCOPLASMAS / UREAPLASMAS

Wells	Antibiotic	Concentrations (µg/mL)
7-TE and 8-TE	Tetracycline	4 and 8
9-ENR and 10-ENR	Enrofloxacin	5 and 10
11-MAR and 12-MAR	Marbofloxacin	10 and 20
13-DO and 14-DO	Doxycycline	4 and 8
15-E and 16-E	Erythromycin	8 and 16
17-CLA and 18-CLA	Clarithromycin	8 and 16
19-TY and 20-TY	Tylosin	30 and 60
21-CD and 22-CD	Clindamycin	4 and 8
23-AZM and 24-AZM	Azithromycin	4 and 8

Well Colour A red* colouring indicates growth (positive) at the tested antibiotic concentration.
A yellow colouring indicates no growth (negative) at the tested antibiotic concentration.

Colour changes in two successive wells containing the same antibiotic are used to categorise results into three susceptibility categories:

Sensitive: Growth is inhibited (YELLOW) by the higher and lower critical concentration of the antibiotic.

Intermediate: Growth is inhibited (YELLOW) by the higher critical concentration of the antibiotic, but not (RED) by the lower critical concentration.

Resistant: Growth is not inhibited (RED) by the higher and lower critical concentration of the antibiotic.

*An orangey colouring should be considered as test positive/growth.

QUALITY CONTROL

Each batch of Mycoplasma System VET is subjected to quality control using the following reference strains:

Mycoplasma hominis ATCC® 23114

Ureaplasma urealiticum ATCC® 27618

Candida albicans ATCC® 10231.

FACTORS THAT MAY INVALIDATE THE RESULTS

Poor standardization of the inoculum; unsuitable material to examine; use of expired systems or expired supplementary reagents; non compliance with temperatures and/or times of incubation.

PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. Mycoplasma System VET is a disposable device to be used for professional use only. It is intended for professional use by properly trained personnel in the laboratory, using approved asepsis and safety methods for handling pathogenic.

STORAGE

Store at 2-8°C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, Mycoplasma System VET and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

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BIBLIOGRAPHY / BIBLIOGRAFIA

1. Razin, S. and Tully, J.G. (1995) Methods and Diagnostic Procedures in Mycoplasma. *California*, 29-31.
2. Rycroft, A.N., Tsounakou, E. and Chalker, V. (2007) Serological Evidence of *Mycoplasma cynos* Infection in Canine Infectious Respiratory Disease. *Veterinary Microbiology*, **120**, 358-362. <http://dx.doi.org/10.1016/j.vetmic.2006.11.011>
3. Hong, S. and Kim, O. (2012) Molecular Identification of *Mycoplasma cynos* from Laboratory Beagle Dogs with Respiratory Disease. *Laboratory Animal Research*, **28**, 61-66. <http://dx.doi.org/10.5625/lar.2012.28.1.61>
4. Chalker, V.J. (2005) Canine Mycoplasmas. *Research Veterinary Science*, **79**, 1-8. <http://dx.doi.org/10.1016/j.rvsc.2004.10.002>
5. Spersger, J., Langer, S., Muck, S., Macher, K., Szostak, M., Rosengarten, R. and Busse, H.J. (2011) *Mycoplasma mucosicanis* sp. nov., Isolated from the Mucosa of Dogs. *International Journal of Systematic Evolution Microbiology*, **61**, 716-721. <http://dx.doi.org/10.1099/ijs.0.015750-0>
6. Rosendal, S. (1974) Canine Mycoplasmas II: Biochemical Characterization and Serological Identification. *Acta Pathologica Microbiologica Scandinavica Section B: Microbiology and Immunology*, **82**, 25-32.
7. Rosendal, S. (1975) Canine Mycoplasmas: Cultural and Biochemical Studies of Type and Reference Strains. *Acta Pathologica Microbiologica Scandinavica Section B: Microbiology and Immunology*, **83**, 457-462.
8. Rosendal, S. (1978) Canine Mycoplasmas: Pathogenicity of Mycoplasmas Associated with Distemper Pneumonia. *The Journal of Infectious Diseases*, **138**, 203-210. <http://dx.doi.org/10.1093/infdis/138.2.203>
9. Rosendal, S. (1982) Canine Mycoplasmas: Their Ecologic Niche and Role in Disease. *Journal of the American Veterinary Medical Association*, **15**, 1212-1214.
10. DaMassa, A.J., Wakenell, P.S. and Brooks, D.L. (1992) Mycoplasmas of Goats and Sheep: Review Article. *Journal of Veterinary Diagnostic Investigation*, **4**, 101-113. <http://dx.doi.org/10.1177/104063879200400126>
11. Corrales, J.C., Esnal, A., De la Fe, C., Sánchez, A., Assunção, P., Poveda, J.B. and Contreras, A. (2007) Contagious Agalactia in Small Ruminants. *Small Ruminant Research*, **68**, 154-166. <http://dx.doi.org/10.1016/j.smallrumres.2006.09.010>
12. Santos, S.B., Nascimento, E.R., Faccini, J.L.H., Barreto, M.L. and Mota, R.A. (2012) Mycoplasmas Species Associated with the Ear Mites in Cattle. *Pesquisa Veterinária Brasileira*, **32**, 293-296. <http://dx.doi.org/10.1590/S0100-736X2012000400004>
13. Edward, D.G. (1955) A Suggested Classification and Nomenclature for Organisms of the Pleuropneumonia Group. *International Journal of Systematic and Evolutionary Microbiology*, **5**, 85-93. <http://dx.doi.org/10.1099/0096266x-5-2-85>
14. Armstrong, D., Tully, J.G., Yu, B., Morton, V., Friedman, M.H. and Steger, L. (1970) Previously Uncharacterized Mycoplasma Isolates from an Investigation of Canine Pneumonia. *Infection and Immunity*, **1**, 1-7.
15. Taylor, R.R., Mohan, K., Miles, R.J. (1996) Diversity of energy-yielding substrates and metabolism in avian mycoplasmas. *Veterinary Microbiology*, **51**, 291-304.

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Mycoplasma System VET

Sample data
Notes

Mark the obtained results in the relevant square

Well	Count and Identification of mycoplasmas/ureaplasmas	Well colour			
		Positive reaction*		Negative reaction	
1-Uu 10³	<i>Ureaplasma</i> spp. (titre = 10 ³ CFU/mL)	red		yellow	
2-Uu 10⁴	<i>Ureaplasma</i> spp. (titre = 10 ⁴ CFU/mL)	red		yellow	
3-Uu ≥ 10⁵	<i>Ureaplasma</i> spp. (titre ≥ 10 ⁵ CFU/mL)	red		yellow	
4-Mh 10⁴	<i>Mycoplasma</i> spp. (titre = 10 ⁴ CFU/mL)	red		yellow	
5-Mh ≥ 10⁵	<i>Mycoplasma</i> spp. (titre ≥ 10 ⁵ CFU/mL)	red		yellow	

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TEST RESULTS FORM

Mycoplasma System VET

Indicate microbial growth with **+**
 Indicate the inhibition of growth with **-**

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Sensitive Growth is inhibited (YELLOW) by the higher and lower critical concentration of the antibiotic.

Intermediate Growth is inhibited (YELLOW) by the higher critical concentration of the antibiotic, but not (RED) by the lower critical concentration.

Resistant Growth is not inhibited (RED) by the higher and lower critical concentration of the antibiotic.

Date of Test	Operator
---------------------	-----------------