Sensi Test ***gram-negative*** **REF. 76010**

- ***Sistema per l'antibiogramma dei batteri gram-negativi e degli enterococchi***
- ***System for the susceptibility testing of gram-negative bacteria and enterococci***
- ***Système pour l'antibiogramme des bactéries gram négatives et des entérocoques***
- ***Sistema para el antibiograma de las bacterias gramnegativas y de los enterococos***
- ***Sistema para o antibiograma das bactérias gram-negativos e dos enterococos***
- ***Σύστημα για αντιβιογράμμα των βακτηριδίων gram-αρνητικό και των εντεροκόκκων***

Code F10613
Rev. 0 / 24.08.2004

Sensi Test gram-negative

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Sensi Test gram-negative

System for the susceptibility testing of gram-negative bacteria and enterococci

DESCRIPTION

Sensi Test gram-negative is a 24-well system containing desiccated antibiotics for the susceptibility testing of gram-negative bacteria and enterococci. The system is inoculated with the bacterial suspension of the micro-organism under examination and incubated at $36^{\circ}\text{C} \pm 1^{\circ}\text{C}$ for 18-24 hours. The susceptibility testing is interpreted by assessing the change in colour of the various wells in the system, and the strain of bacteria is defined as Sensitive (S), Intermediate (I) or Resistant (R).

KIT CONTENTS

The kit contains:

- 20 **Sensi Test gram-negative** systems
- 1 Instruction sheet
- 20 vials of Suspension Medium (7 ml / vial)
- 20 ANTIBIOGRAMMA/SUSCEPTIBILITY TESTING Forms

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

- SUSPENSION BROTHcode 20089
- Sundry microbiology laboratory materials
- PHYSIOLOGICAL SOLUTION.....code 20095

CONFIGURATION

The configuration of the system is shown in Table n°1.

Table n°1

Well	ANTIBIOTIC (*)
1-AK	Amikacin – 32 µg / ml
2-CN	Gentamicin – 8 µg / ml
3-TOB	Tobramycin - 8 µg / ml
4-AMS	Ampicillin / Sulbactam - 32/16 µg / ml
5-AUG	Amoxicillin / Clavulanic acid - 32/16 µg / ml
6-MEZ	Mezlocillin – 128 µg / ml
7-TZP	Piperacillin / Tazobactam - 128/4 µg / ml
8-TCC	Ticarcillin / Clavulanic acid - 128/2 µg / ml
9-FOS	Fosfomicin – 200 µg / ml
10-F	Nitrofurantoin – 128 µg / ml
11-CFP	Cefoperazone – 64 µg / ml
12-CTX	Cefotaxime – 64 µg / ml
13-CAZ	Ceftazidime – 32 µg / ml
14-CRO	Ceftriaxone – 64 µg / ml
15-CXM	Cefuroxime – 32 µg / ml
16-NA	Nalidixic acid – 32 µg / ml
17-LOM	Lomefloxacin – 8 µg / ml
18-CIP	Ciprofloxacin – 4 µg / ml
19-LEV	Levofloxacin – 8 µg / ml
20-PEF	Pefloxacin – 8 µg / ml
21-TEC	Teicoplanin – 32 µg / ml
22-ATM	Aztreonam – 32 µg / ml
23-SXT	Co-trimoxazole – 8 µg / ml
24-C	Growth control

(*):The concentration of each antibiotic is according to the updated NCCLS–January2004; Vol 24 N°1⁽¹⁾

PRINCIPLE OF THE METHOD

Sensi Test gram-negative makes it possible to execute the susceptibility testing of :

- Urinary enterococci and enterobacteria;
- Enterobacteria isolated from other clinical samples;
- *Pseudomonas spp.* and non-fermenting bacteria;
- Gram-negative bacteria of veterinary origin.

The susceptibility testing is assessed on the basis of the growth or inhibition of the micro-organisms in media containing antibiotic and a growth indicator in wells from **1-AK** to **23-SXT**. Well **24-C** does not contain antibiotics, only the culture medium, and the indicator therefore serves as a control for bacterial growth.

COMPOSITION

Table n°2

Well	Content
1-AK	Culture medium with Amikacin - 32 µg/ml
2-CN	Culture medium with Gentamicin - 8 µg/ml
3-TOB	Culture medium with Tobramycin - 8 µg/ml
4-AMS	Culture medium with Ampicillin / Sulbactam - 32/16 µg/ml
5-AUG	Culture medium with Amoxicillin / Clavulanic acid - 32/16 µg/ml
6-MEZ	Culture medium with Mezlocillin – 128 µg/ml
7-TZP	Culture medium with Piperacillin / Tazobactam - 128/4 µg/ml
8-TCC	Culture medium with Ticarcillin / Clavulanic acid - 128/2 µg/ml
9-FOS	Culture medium with Fosfomycin - 200 µg/ml
10-F	Culture medium with Nitrofurantoin - 128 µg/ml
11-CFP	Culture medium with Cefoperazone - 64 µg/ml
12-CTX	Culture medium with Cefotaxime - 64 µg/ml
13-CAZ	Culture medium with Ceftazidime - 32 µg/ml
14-CRO	Culture medium with Ceftriaxone - 64 µg/ml
15-CXM	Culture medium with Cefuroxime - 32 µg/ml
16-NA	Culture medium with Nalidixic acid - 32 µg/ml
17-LOM	Culture medium with Lomefloxacin - 8 µg/ml
18-CIP	Culture medium with Ciprofloxacin - 4 µg/ml
19-LEV	Culture medium with Levofloxacin - 8 µg/ml
20-PEF	Culture medium with Pefloxacin - 8 µg/ml
21-TEC	Culture medium with Teicoplanin - 32 µg/ml
22-ATM	Culture medium with Aztreonam – 32 µg/ml
23-SXT	Culture medium with Co-trimoxazole - 8 µg/ml
24-C	Culture medium without antibiotics

Suspension Medium (g/l): Yeast extract 5 g; Beef peptone 3 g; Glucose 2 g; Distilled water 1000 ml;
pH 6.8 ± 0.2

COLLECTION AND STORAGE OF THE SAMPLES

The colonies to be subjected to the susceptibility testing with **Sensi Test gram-negative** must be taken from a selective or non-selective culture medium used for isolation of the micro-organism under examination.

TEST PROCEDURE

- Take a system from its wrapper and bring it to ambient temperature.
- Write down the name, date and source of the micro-organism.
- Prepare a bacterial suspension of opacity equivalent to that of the standard 0.5 McFarland, taking one or more morphologically similar, well-isolated, colonies from a solid culture medium and suspending them in 5 ml of SUSPENSION BROTH (code 20089) or in 5 ml of PHYSIOLOGICAL SOLUTION (code 20095).
- Take a vial of the Suspension Medium contained in the kit.
- Transfer 0.01 ml of bacterial suspension in SUSPENSION BROTH (code 20089) into the vial of Suspension Medium contained in the kit and distribute 0.2 ml (4 drops) into each well of the system
- Cover the system with the lid provided and incubate at $36^{\circ}\text{C} \pm 1^{\circ}\text{C}$ for 18-24 hours.

INTERPRETATION OF THE RESULTS

At the end of the incubation, observe the change in colour of the wells and interpret the results using Table 3. Transcribe the data onto the ANTIBIOGRAMMA/SUSCEPTIBILITY TESTING form provided.

The control well (**24-C**) must be positive (yellow). In the event that it is negative (blue or grey), the vitality of the sample colonies must be checked and the test repeated using a new system and a recently grown bacterial culture.

Table n°3

WELL COLOUR	BACTERIAL GROWTH	INTERPRETATION
Blue	Inhibited	S = Sensitive
Grey	Partially inhibited	I = Intermediate
Yellow	Good	R = Resistant

QUALITY CONTROL

Each batch of **Sensi Test gram-negative** is subjected to quality control using the following reference strains of bacteria:

<i>Escherichia coli</i>	ATCC 25922
<i>Pseudomonas aeruginosa</i>	ATCC 27853
<i>Enterobacter cloacae</i>	ATCC 13048
<i>Enterococcus faecalis</i>	ATCC 19433
<i>Klebsiella pneumonite</i>	ATCC 13883

FACTORS THAT MAY INVALIDATE THE RESULTS

- Use of mixed cultures
- Inaccurate standardisation of the inoculum
- Use of inoculum solutions other than those suggested
- Application of the method to slow growing or anaerobic bacteria
- Use of expired systems
- Temperature and times for incubation other than those suggested
- Lack of growth in well **24-C**

LIMITS AND WARNINGS

The results obtained with the system **Sensi Test gram-negative** must be interpreted by a clinical microbiologist before proceeding to therapeutic treatment of the patient.

PERFORMANCE

The system **Sensi Test gram-negative** makes it possible to execute the susceptibility testing according to the method of micro-dilution in broth. The results obtained with the system **Sensi Test gram-negative** agree with those obtained following the standards indicated by the National Committee for Clinical Laboratory Standards (NCCLS) U.S.A. ^(2,4) and with those from the radial diffusion method according to Bauer et al. ^(3,4) recommended by the Food and Drug Administration (FDA).

PRECAUTIONS

The product, **Sensi Test gram-negative**, is not classified as hazardous under current legislation nor does it contain harmful substances in concentrations $\geq 1\%$. It therefore does not require a Safety Data Sheet to be available. **Sensi Test gram-negative** is a disposable device to be used only for diagnostic use *in vitro*. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8°C in the original packaging. In such conditions the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.


DISPOSAL OF USED MATERIAL

After use, **Sensi Test gram-negative** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Code	Kit
Sensi Test gram-negative	76010	20 tests

TABLE OF SYMBOLS

<u>SYMBOL</u>	<u>MEANING</u>	<u>SYMBOL</u>	<u>MEANING</u>
	<i>In Vitro</i> Diagnostic Medical Device		Do not reuse
	Manufacturer		Contains sufficient for <n> tests
	Catalogue number		Fragile, handle with care
	Use by		Caution, consult accompanying documents
	Temperature limitation		Batch code

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IVD

1. NCCLS - Performance Standards for antimicrobial Susceptibility Testing, fourteen Informational Supplement- January 2004, Vol 24 N°1.
2. NCCLS -Methods for Dilution Antimicrobial Susceptibility tests for Bacteria that grow aerobically; Approved Standard – Sixth Edition (2003).
3. Bauer, A.W., W.M.M. Kirby, J.C. Sherris and M. Turk. 1966. Antimicrobial susceptibility testing by a standardized single disk method. *AM. J. Clin. Pathol.* **45**:493-496.
4. Dati di archivio Liofilchem s.r.l."SENSI TEST GRAM NEGATIVE" (Aprile 2003).

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