Sensi Test gram-positive

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Sensi Test gram-positive

System for the susceptibility testing of gram-positive bacteria

DESCRIPTION

Sensi Test gram-positive is a 24-well system containing desiccated antibiotics for the susceptibility testing of gram-positive bacteria. The system is inoculated with the bacterial suspension of the micro-organism under examination and incubated at 36°C ± 1°C for 18-24 hours. The susceptibility testing is interpreted by assessing the change in colour of the various wells in the system, and the strain of bacteria is defined as Sensitive (S), Intermediate (I) or Resistant (R).

KIT CONTENTS

The kit contains:

20 Sensi Test gram-positive systems	20 vials of <i>Inoculum Broth</i> (7 ml / vial)
Instruction sheet	20 ANTIBIOGRAMMA forms

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

SUSPENSION BROTHcode 20089	PHYSIOLOGICAL SOLUTION.....code 20095
Sundry microbiology laboratory materials	

CONFIGURATION

The configuration of the system is shown in Table n°1.

Table n°1

Well	ANTIBIOTIC (*)
1-CN	Gentamicin – 8 µg / mL
2-NET	Netilmicin – 32 µg / mL
3-TOB	Tobramycin - 8 µg / mL
4- AUG	Amoxicillin / Clavulanic acid - 8/4 µg / mL
5- AMS	Ampicillin / Sulbactam - 32/16 µg / mL
6- PIP	Piperacillin – 16 µg / mL
7-CEC	Cefaclor - 32 µg / mL
8-CFM	Cefixime - 32 µg / mL
9-CID	Cefonicid – 32 µg / mL
10-CAZ	Ceftazidime – 32 µg / mL
11-CRO	Ceftriaxone – 64 µg / mL
12-CXM	Cefuroxime – 32 µg / mL
13-CIP	Ciprofloxacin – 4µg / mL I
14-LEV	Levofloxacin – 8 µg / mL
15-PEF	Pefloxacin – 8 µg / mL
16-AZM	Azithromycin – 8 µg / mL
17-CLR	Clarithromycin – 8 µg / mL
18-ERY	Erythromycin – 8 µg / mL
19-MK	Miokamycin – 8 µg / mL
20-ROX	Roxitromycin – 8 µg / mL
21-FOS	Fosfomicin – 200 µg / mL
22-RD	Rifampicin – 4 µg / mL
23-SXT	Co-trimoxazole – 8 µg / mL
24-C	Growth control

(*):The concentration of each antibiotic is according to the updated NCCLS–January2004; Vol 24 N°1⁽¹⁾

PRINCIPLE OF THE METHOD

Sensi Test gram-positive makes it possible to execute the susceptibility testing of :

- Staphylococci, streptococci
- Gram-positive bacteria of veterinary origin.

The susceptibility testing is assessed on the basis of the growth or inhibition of the micro-organisms in media containing antibiotic and a growth indicator in wells from 1-CN to 23-SXT. Well 24-C does not contain antibiotics, only the culture medium, and the indicator therefore serves as a control for bacterial growth.

COMPOSITION

Each well contains an essiccated culture medium and a growth indicator. After the inoculum of the system with the microorganism in *Inoculum broth* the obtained antibiotic concentrations of the wells are those indicated on the Table n°1.

Inoculum broth (g/L): *Glucose anhydrous 30.0 g; Tryptone 5.0g; Chloride Sodium 1000 ml; Extract yeast 3.0g; Agar 1.0g; Distilled water 1000 mL; pH 7.5 ± 0.1*

COLLECTION AND STORAGE OF THE SAMPLES

The colonies to be subjected to the susceptibility testing with **Sensi Test gram-positive** must be taken from a selective or non-selective culture medium used for isolation of the micro-organism under examination.

TEST PROCEDURE

- Take a system from its wrapper and bring it to ambient temperature.
- Write down the name, date and source of the micro-organism.
- Prepare a bacterial suspension of opacity equivalent to that of the standard 0.5 McFarland, taking one or more morphologically similar, well-isolated, colonies from a solid culture medium and suspending them in 5 ml of SUSPENSION BROTH (code 20089) or in 5 ml of PHYSIOLOGICAL SOLUTION (code 20095).
- Take a vial of the *Inoculum Broth* contained in the kit.
- Transfer 0.2 ml of bacterial suspension in SUSPENSION BROTH (code 20089) or in PHYSIOLOGICAL SOLUTION (code 20095) into the vial of *Inoculum Broth* contained in the kit and distribute 0.2 ml (4 drops) into each well of the system.
- Cover the system with the lid provided and incubate at 36°C ± 1°C for 18-24 hours.

INTERPRETATION OF THE RESULTS

At the end of the incubation, observe the change in colour of the wells and interpret the results using Table 2. Transcribe the data onto the ANTIBIOGRAMMA forms provided.

The control well (**24-C**) must be positive (yellow). In the event that it is negative (red or orange), the vitality of the sample colonies must be checked and the test repeated using a new system and a recently grown bacterial culture.

Table n°2

WELL COLOUR	BACTERIAL GROWTH	INTERPRETATION
Red	Inhibited	S = Sensitive
Orange	Partially inhibited	I = Intermediate
Yellow	Good	R = Resistant

QUALITY CONTROL

Each batch of **Sensi Test gram-positive** is subjected to quality control using the following reference strains of bacteria:

Staphylococcus aureus..... ATCC 25923
Streptococcus pneumoniae..... ATCC 49619

FACTORS THAT MAY INVALIDATE THE RESULTS

- Use of mixed cultures
- Inaccurate standardisation of the inoculum
- Use of inoculum solutions other than those suggested
- Application of the method to slow growing or anaerobic bacteria
- Use of expired systems
- Temperature and times for incubation other than those suggested
- Lack of growth in well **24-C**

LIMITS AND WARNINGS

The results obtained with the system **Sensi Test gram-positive** must be interpreted by a clinical microbiologist before proceeding to therapeutic treatment of the patient.

PERFORMANCE

The system **Sensi Test gram-positive** makes it possible to execute the susceptibility testing according to the method of micro-dilution in broth. The results obtained with the system **Sensi Test gram-positive** agree with those obtained following the standards indicated by the National Committee for Clinical Laboratory Standards (NCCLS) U.S.A. ^(2,4) and with those from the radial diffusion method according to Bauer et al. ^(3,4) recommended by the Food and Drug Administration (FDA).

PRECAUTIONS

The product, **Sensi Test gram-positive**, is not classified as hazardous under current legislation nor does it contain harmful substances in concentrations $\geq 1\%$. It therefore does not require a Safety Data Sheet to be available. **Sensi Test gram-positive** is a disposable device to be used only for diagnostic use *in vitro*. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8°C in the original packaging. In such conditions the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **Sensi Test gram-positive** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Code	Kit
Sensi Test gram-positive	76020	20 tests

TABLE OF SYMBOLS

IVD <i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

BIBLIOGRAFIA / BIBLIOGRAPHY

1. NCCLS - Performance Standards for antimicrobial Susceptibility Testing, fourteenth Informational Supplement- January 2004, Vol 24 N°1.
2. NCCLS -Methods for Dilution Antimicrobial Susceptibility tests for Bacteria that grow aerobically; Approved Standard – Sixth Edition (2003).
3. Bauer, A.W., W.M.M. Kirby, J.C. Sherris and M. Turk. 1966. Antimicrobial susceptibility testing by a standardized single disk method. *AM. J. Clin. Pathol.* **45**:493-496.
4. Dati di archivio Liofilchem s.r.l."SENSI TEST GRAM POSITIVE" (Giugno 2005).

LIOFILCHEM
Bacteriology Products

LIOFILCHEM S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto D.A. (TE) - ITALY

Tel +390858930745 Fax +390858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net