SensiQuattro Gram-negative

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SensiQuattro Gram-negative

ENGLISH

System for determination of the Minimum Inhibitory Concentration (MIC)
of gram-negative bacteria with the broth microdilution method.

DESCRIPTION

SensiQuattro Gram-negative is a 40-well plate containing 9 dried antibiotics in four two-fold dilutions and 1 antibiotic in three two-fold dilutions.

The system allows to perform the antibiogram with evaluation of the minimum inhibitory concentration (MIC), for gram-negative bacteria towards antibiotic agents establishing the sensitivity of the microorganism according to the EUCAST (*European Committee on Antimicrobial Susceptibility testing*) or CLSI (*Clinical and Laboratory Standards Institute*) interpretative criteria.

KIT CONTENT

ref. 76031	<ul style="list-style-type: none"> • 20 Systems SensiQuattro Gram-negative • 20 Vials of Physiological Solution • 20 Vials of Suspension Medium 	ref. 79031	<ul style="list-style-type: none"> • 4 Systems SensiQuattro Gram-negative • 4 Vials of Physiological Solution • 4 Vials of Suspension Medium
Instruction Sheet		TEST RESULTS FORM	

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

Mc Farland 0.5 Barium Sulphate Standard (ref.80400)	Tips for multi.channel pipette (ref.96758)
Multi.channel pipette 50-300 µl (ref.96759)	Solution reservoir for multichannel pipette (ref. 96761)

CONFIGURATION

WELL	ANTIBIOTIC	µg/ml	WELL	ANTIBIOTIC	µg/ml
CN	Gentamicin	2	AMS	Ampicillin / Sulbactam	8 / 4
CN	Gentamicin	4	AMS	Ampicillin / Sulbactam	16 / 8
CN	Gentamicin	8	AMS	Ampicillin / Sulbactam	32 / 16
CN	Gentamicin	16	AMS	Ampicillin / Sulbactam	64 / 32
AK	Amikacin	8	CIP	Ciprofloxacin	0.5
AK	Amikacin	16	CIP	Ciprofloxacin	1
AK	Amikacin	32	CIP	Ciprofloxacin	2
AK	Amikacin	64	CIP	Ciprofloxacin	4
TOB	Tobramycin	4	LEV	Levofloxacin	1
TOB	Tobramycin	8	LEV	Levofloxacin	2
TOB	Tobramycin	16	LEV	Levofloxacin	4
TOB	Tobramycin	32	LEV	Levofloxacin	8
TZP	Piperacillin / Tazobactam	16 / 4	CXM	Cefuroxime	4
TZP	Piperacillin / Tazobactam	32 / 4	CXM	Cefuroxime	8
TZP	Piperacillin / Tazobactam	64 / 4	CXM	Cefuroxime	16
TZP	Piperacillin / Tazobactam	128 / 4	CXM	Cefuroxime	32
FOS	Fosfomicin	32	SXT	Co-trimoxazole	2 / 38
FOS	Fosfomicin	64	SXT	Co-trimoxazole	4 / 76
FOS	Fosfomicin	128	SXT	Co-trimoxazole	8 / 152
FOS	Fosfomicin	256	Growth	Growth control	

PRINCIPLE OF THE METHOD

The system is inoculated and rehydrated with a standardized microbial suspension and incubated at $36\pm 1^{\circ}\text{C}$ for 18-24 hours. After incubation the susceptibility test (MIC) is determined on the basis of growth or inhibition of microorganism in the well containing the antibiotic, nutritional factors and a growth indicator. The well **Growth** does not contain antibiotics and is used as a control of microbial growth.

COLLECTION AND STORAGE OF THE SAMPLE

The colonies to be subjected to the susceptibility test (MIC) must be properly isolated and taken from a culture medium, selective or non-selective, used for isolation of the microorganism.

TEST PROCEDURE

1. Take a system **SensiQuattro Gram-negative** from its wrapper and let it reach room temperature.
2. Write down the patient's name, date and origin of the organism.
3. Take several morphologically similar, well-isolated, colonies and suspend them in a vial of *Physiological Solution* from the kit and obtain a suspension of opacity equivalent to the standard 0.5 McFarland (ref. 80400).
4. Transfer 0.2 ml of the bacterial suspension in a vial of Suspension Medium* contained in the kit and drain the solution into the reservoir for multi-channel pipetting (ref. 96761).
5. Dispense 0.15 ml of bacterial suspension into each well of the system using an eight tips multichannel pipette or a single tip pipette.
6. Cover the system with the lid provided and incubate at $36\pm 1^{\circ}\text{C}$ for 18-24 hours.

* Suspension Medium (g/l): Yeast Extract 5.0 g; Beef Peptone 3.0 g; Glucose 2.0 g; Distilled Water 1000 ml;
pH 6.8 ± 0.2

INTERPRETATION OF THE RESULTS

After incubation, observe the color change in the wells and establish the values of minimum inhibitory concentration using table No. 2.

- The control well Growth must be yellow in color. If it turns blue, check the viability of the colonies picked and repeat the test using a new system and a microbial culture of recent growth.

Table No. 2.

Well colour	Microbial growth
Blue	Inhibited
Yellow/Grey	Good

- The obtained values of minimum inhibitory concentration (MIC) must be interpreted according to the current EUCAST ⁽¹⁾ or CLSI ⁽²⁾ interpretative criteria.
- Note the results on the TEST RESULTS FORM included in the kit (copy as many form as necessary).

QUALITY CONTROL

Each batch of **SensiQuattro Gram-negative** is subjected to quality control using the following microorganisms:

<i>Escherichia coli</i>	ATCC® 25922	<i>Klebsiella pneumoniae</i>	ATCC® 700603
<i>Enterobacter cloacae</i>	ATCC® 13047	<i>Pseudomonas aeruginosa</i>	ATCC® 27853
<i>Klebsiella pneumoniae</i>	ATCC® 13883	<i>Salmonella typhimurium</i>	ATCC® 14028

FACTORS THAT MAY INVALIDATE THE RESULTS

Inaccurate standardization of the inoculum; clinical material not suitable; use of additional systems and reagents; temperature and incubation time not respected.

PRECAUTIONS

The product, **SensiQuattro Gram-negative**, is not classified as hazardous under current legislation, however see the safety data sheet for a correct use. **SensiQuattro Gram-negative** is a disposable device to be used only for diagnostic use *in vitro* and must be used in the laboratory by appropriately trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE AND TRANSPORT CONDITIONS

Store **SensiQuattro Gram-negative** at 2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the storage or transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, do not alter in any way the performance of the product. Discard the panel if there are visible signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **SensiQuattro Gram-negative** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with guidelines used in the laboratory for decontamination and disposal of potentially infected material.

REFERENCES

1. EUCAST. Breakpoint tables for interpretation of MICs and zone diameters.
2. CLSI M100-S23. Performance Standards for Antimicrobial Susceptibility Testing.

PRESENTATION

Product	Ref.	Packaging
SensiQuattro Gram-negative	76031	20 tests
SensiQuattro Gram-negative	79031	4 tests

TABLE OF SYMBOLS

IVD <i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code