SensiQuattro Gram-positive

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SensiQuattro Gram-positive

ENGLISH

System for determination of the Minimum Inhibitory Concentration (MIC)
of gram-positive bacteria with the broth microdilution method.

DESCRIPTION

SensiQuattro Gram-positive is a 40-well plate containing 9 dried antibiotics in four two-fold dilutions and 1 antibiotic in three two-fold dilutions.

The system allows to perform the antibiogram with evaluation of the minimum inhibitory concentration (MIC), for gram-positive bacteria towards antibiotic agents establishing the sensitivity of the microorganism according to the EUCAST (*European Committee on Antimicrobial Susceptibility testing*) or CLSI (*Clinical and Laboratory Standards Institute*) interpretative criteria.

KIT CONTENT

ref. 76032	<ul style="list-style-type: none"> • 20 Systems SensiQuattro Gram-positive • 20 Vials of Physiological Solution • 20 Vials of Inoculum Broth 	ref. 79032	<ul style="list-style-type: none"> • 4 Systems SensiQuattro Gram-positive • 4 Vials of Physiological Solution • 4 Vials of Inoculum Broth
Instruction Sheet		- TEST RESULTS FORM	

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

Mc Farland 0.5 Barium Sulphate Standard (ref.80400)	Tips for multi.channel pipette (ref.96758)
Multi.channel pipette 50-300 µl (ref.96759)	Solution reservoir for multichannel pipette (ref. 96761)

CONFIGURATION

WELL	ANTIBIOTIC	µg/ml	WELL	ANTIBIOTIC	µg/ml
CN	Gentamicin	2	CRO	Ceftriaxone	8
CN	Gentamicin	4	CRO	Ceftriaxone	16
CN	Gentamicin	8	CRO	Ceftriaxone	32
CN	Gentamicin	16	CRO	Ceftriaxone	64
LEV	Levofloxacin	1	AUG	Amoxicillin / Clavulanic acid	4 / 2
LEV	Levofloxacin	2	AUG	Amoxicillin / Clavulanic acid	8 / 4
LEV	Levofloxacin	4	AUG	Amoxicillin / Clavulanic acid	16 / 8
LEV	Levofloxacin	8	AUG	Amoxicillin / Clavulanic acid	32 / 16
AMS	Ampicillin / Sulbactam	8 / 4	CLR	Clarithromycin	1
AMS	Ampicillin / Sulbactam	16 / 8	CLR	Clarithromycin	2
AMS	Ampicillin / Sulbactam	32 / 16	CLR	Clarithromycin	4
AMS	Ampicillin / Sulbactam	64 / 32	CLR	Clarithromycin	8
VA	Vancomycin	2	RD	Rifampicin	0.5
VA	Vancomycin	4	RD	Rifampicin	1
VA	Vancomycin	8	RD	Rifampicin	2
VA	Vancomycin	16	RD	Rifampicin	4
LNZ	Linezolid	2	SXT	Co-trimoxazole	2 / 38
LNZ	Linezolid	4	SXT	Co-trimoxazole	4 / 76
LNZ	Linezolid	8	SXT	Co-trimoxazole	8 / 152
LNZ	Linezolid	16	Growth	Growth control	

PRINCIPLE OF THE METHOD

The system is inoculated and rehydrated with a standardized microbial suspension and incubated at $36\pm 1^{\circ}\text{C}$ for 18-24 hours. After incubation the susceptibility test (MIC) is determined on the basis of growth or inhibition of microorganism in the well containing the antibiotic, nutritional factors and a growth indicator. The well **Growth** does not contain antibiotics and is used as a control of microbial growth.

COLLECTION AND STORAGE OF THE SAMPLE

The colonies to be subjected to the susceptibility test (MIC) must be properly isolated and taken from a culture medium, selective or non-selective, used for isolation of the microorganism.

TEST PROCEDURE

1. Take a system **SensiQuattro Gram-positive** from its wrapper and let it reach room temperature.
2. Write down the patient's name, date and origin of the organism.
3. Take several morphologically similar, well-isolated, colonies and suspend them in a vial of Physiological Solution from the kit and obtain a suspension of opacity equivalent to the standard 0.5 McFarland (ref. 80400).
4. Transfer 0.2 ml of the bacterial suspension in a vial of Inoculum Broth* contained in the kit and drain the solution into the reservoir for multi-channel pipetting (ref. 96761).
5. Dispense 0.15 ml of bacterial suspension into each well of the system using an eight tips multichannel pipette or a single tip pipette.
6. Cover the system with the lid provided and incubate at $36\pm 1^{\circ}\text{C}$ for 18-24 hours.

* Inoculum Broth (g/l): Glucose 30.0 g; Tryptone 5.0 g; NaCl 9.0 g; Yeast Extract 5.0 g; Agar 1.0 g; Distilled Water 1000 ml; pH 7.0 ± 0.2

INTERPRETATION OF THE RESULTS

After incubation, observe the color change in the wells and establish the values of minimum inhibitory concentration using table No. 2.

- The control well Growth must be yellow in color. If it turns blue, check the viability of the colonies picked and repeat the test using a new system and a microbial culture of recent growth.

Table No. 2.

Well colour	Microbial growth
Red	Inhibited
Yellow/Orange	Good

- The obtained values of minimum inhibitory concentration (MIC) must be interpreted according to the current EUCAST ⁽¹⁾ or CLSI ⁽²⁾ interpretative criteria.
- Note results on the TEST RESULTS FORM included in the kit (copy as many form as necessary).

QUALITY CONTROL

Each batch of **SensiQuattro Gram-positive** is subjected to quality control using the following microorganisms:

<i>Staphylococcus aureus</i>	ATCC® 29213	<i>Streptococcus pneumoniae</i>	ATCC® 49619
<i>Staphylococcus xylosus</i>	ATCC® 29971	<i>Streptococcus agalactiae</i>	ATCC® 13813
<i>Streptococcus pyogenes</i>	ATCC® 19615	<i>Enterococcus faecalis</i>	ATCC® 29212

FACTORS THAT MAY INVALIDATE THE RESULTS

Inaccurate standardization of the inoculum; clinical material not suitable; use of additional systems and reagents; temperature and incubation time not respected.

PRECAUTIONS

The product, **SensiQuattro Gram-positive**, is not classified as hazardous under current legislation, however see the safety data sheet for a correct use. **SensiQuattro Gram-positive** is a disposable device to be used only for diagnostic use *in vitro* and must be used in the laboratory by appropriately trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE AND TRANSPORT CONDITIONS

Store **SensiQuattro Gram-positive** at 2-8°C away from light, until the expiry date on the label. However, our stability studies have shown that the storage or transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, do not alter in any way the performance of the product. Discard the panel if there are visible signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **SensiQuattro Gram-positive** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with guidelines used in the laboratory for decontamination and disposal of potentially infected material.

REFERENCES

1. EUCAST. Breakpoint tables for interpretation of MICs and zone diameters.
2. CLSI M100-S23. Performance Standards for Antimicrobial Susceptibility Testing.

PRESENTATION

Product	Ref.	Packaging
SensiQuattro Gram-positive	76032	20 tests
SensiQuattro Gram-positive	79032	4 tests

TABLE OF SYMBOLS

IVD <i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code