URIN SYSTEM Plus

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URIN SYSTEM Plus

ENGLISH

System for microbial count, identification and susceptibility testing directly from urine specimens.

DESCRIPTION

URIN SYSTEM Plus is a 24-well system containing desiccated biochemical and antibiotic substrates for microbial count, identification and susceptibility testing directly from urine specimens. The system is inoculated with a dilution of urine after the presence of microorganisms has been ascertained by microscope examination of the urinary sediment, and incubated at 36 ± 1 °C for 18-24 hours. The tests results are interpreted by assessing the color change in the various wells.

KIT CONTENTS

The kit contains:

Ref. 74160	20 URIN SYSTEM Plus	20 Vials of Suspension Medium (7 mL/vial)	20 Vials of Physiological Solution (4.5 mL/vial)
Ref. 79160	4 URIN SYSTEM Plus	4 Vials of Suspension Medium (7 mL/vial)	4 Vials of Physiological Solution (4.5 mL/vial)
1 Instruction sheet		-	1 TEST RESULTS FORM

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

Single-channel pipette, 10 µL	Solution reservoir for multichannel pipette (ref. 96761)
Multichannel pipette, 50-300 µL (ref. 96759)	Glass slides and covers
Tips for multichannel or single-channel pipettes (ref. 96758)	Microscope

PRINCIPLE OF THE METHOD

URIN SYSTEM Plus allows the microbial count, identification and susceptibility testing directly from urine specimens in 18-24 hours.

- The total bacterial count is evaluated by the microbial growth in the wells **1-GR+**, and **2-GR++** which contain culture media and specific growth indicators.
- The microbial identification is based on the growth in differential media that contain selective agents for various microorganisms, in the wells **3-ESC** to **9-CAN**.
- The susceptibility testing is assessed on the basis of the growth or inhibition of the microorganisms in the media containing an antimicrobial agent and a growth indicator in the wells **10-AK** to **23-SXT**. The well **24-C** is the microbial growth control. It contains the culture medium (the same used in the susceptibility testing wells) and a growth indicator.

CONFIGURATION

The configuration of the system is shown in Table 1.

Table 1.

Well	TOTAL MICROBIAL COUNT	
1-GR+	$10^5 \leq \text{CFU/mL} \leq 10^6$	
2-GR++	CFU/mL > 10^6	
Well	MICROBIAL IDENTIFICATION	
3-ESC	<i>Escherichia coli</i>	
4-PRO	<i>Proteus</i> spp. / <i>Providencia</i> spp.	
5-PSE	<i>Pseudomonas</i> spp.	
6-KES	KES Group (<i>Klebsiella</i> , <i>Enterobacter</i> , <i>Serratia</i>)	
7-STR	<i>Enterococcus</i> spp.	
8-STA	<i>Staphylococcus aureus</i>	
9-CAN	<i>Candida</i> spp.	
Well	SUSCEPTIBILITY TESTING	Concentration (µg/mL)
10-AK	Amikacin	32
11-CN	Gentamicin	8
12-TOB	Tobramycin	8
13-TZP	Piperacillin / Tazobactam	128/4
14 FOS	Fosfomicin	200
15-CFP	Cefoperazone	64
16-CTX	Cefotaxime	64
17-CAZ	Ceftazidime	32
18-AMS	Ampicillin / Sulbactam	32/16
19-TE	Tetracycline	16
20-CIP	Ciprofloxacin	4
21-LEV	Levofloxacin	8
22-AUG	Amoxicillin / Clavulanic acid	32/16
23-SXT	Co-trimoxazole	8
24-C	Growth control for susceptibility testing	

COLLECTION AND CONSERVATION OF THE SAMPLES

Collect the urine aseptically in a sterile container at least 6 hours after the previous urination, it is in any case best to collect the first urine in the morning. The collection most used methods are that of *mitto intermedio* for adults who urinate on command and that of the adhesive plastic bag in early infancy. In special cases it may be necessary to obtain samples by catheter from the bladder or by supra-pubic puncture. The samples of urine must be sent to the laboratory to seed the **URIN SYSTEM Plus** immediately after collection. If that is not possible, store the sample in the refrigerator at 2-8 °C for up to 24 hours. Sterile containers, containing preservatives, are available commercially for the storage of urine that cannot be processed within a short time. If these containers are used, always follow the manufacturer's recommendations. The urine must be collected before the start of any antibacterial treatment or at least 48 hours after the last administration of antimicrobial drug.

TEST PROCEDURE**1) PREPARATION OF THE SAMPLE**


- Examine the urinary sediment under the microscope to determine whether there is a significant presence of microorganisms.
- Transfer 0.5 mL of urine into a vial of Physiological Solution contained in the kit (**Suspension A**).
- Transfer 0.2 mL of Suspension A into a vial of Suspension Medium* contained in the kit (**Suspension B**).

*Suspension Medium (g/L):

Mueller Hinton Broth 21g, Yeast extract 5g;
Beef peptone 3g; Glucose 2g;
Distilled water 1000 mL; pH 6,8 ± 0,2

2) INOCULATION OF THE SYSTEM

- Take a system from its wrap and bring it to room temperature.
- Write down the name of the patient and the date of the start of the examination.
- Transfer the **Suspension A** into a reservoir solution (ref. 96761) and using a multichannel or single channel pipette with sterile tips, dispense 0.2 mL of Suspension A into the wells from **1-GR+** to **9-CAN** (**Microbial count and Identification**).
- Transfer the **Suspension B** into another reservoir solution (ref. 96761) and using a multichannel or single channel pipette with sterile tips, dispense 0.2 mL of Suspension B into the wells from **10-AK** to **24-C** (**Susceptibility Testing**).
- Cover the system with the lid provided and incubate at 36 ± 1 °C for 18-24 hours.
- After incubation, watch for the color change in the wells and interpret the results.

WORKFLOW**INTERPRETATION OF THE RESULTS**

- At the end of the incubation, watch for the color change in the wells and interpret the results using Table 2. The control well (**24-C**) must be positive (grey-yellow). In the event that it is negative (blue) the examination must be repeated using a new system.
- Note the results on the TEST RESULTS FORM (copy as many forms as necessary).

Table 2.

Well	TOTAL MICROBIAL COUNT	Well colour	
		Positive reaction	Negative reaction
1-GR+	$10^5 \leq \text{CFU/mL} \leq 10^6$ (**)	yellow	blue
2-GR++	$\text{CFU/mL} > 10^6$ (**)	yellow	blue

(**) Check that the result conforms with that from the microscope examination of the urinary sediment.

Well	MICROBIAL IDENTIFICATION	Well colour	
		Positive reaction	Negative reaction
3-ESC	<i>Escherichia coli</i>	blue	grey-red
4-PRO	<i>Proteus</i> spp. / <i>Providencia</i> spp.	brown	yellow
5-PSE	<i>Pseudomonas</i> spp.	turbid green	yellow
6-KES	KES Group (<i>Klebsiella</i> spp., <i>Enterobacter</i> spp., <i>Serratia</i> spp.)	yellow	violet
7-STR	<i>Enterococcus</i> spp.	black	yellow
8-STA	<i>Staphylococcus aureus</i>	black bottom	yellow
9-CAN	<i>Candida</i> spp.	yellow	green

Well	SUSCEPTIBILITY TESTING		
	Well colour	Microbial growth	Interpretation
10-AK to 23-SXT	blue	inhibited	sensitive
	grey-yellow	good	resistant

Well	GROWTH CONTROL	Well colour	
		Positive reaction	Negative reaction
24-C	Growth control for the susceptibility testing	grey-yellow	blue

QUALITY CONTROL

Each batch of **URIN SYSTEM Plus** is subjected to quality control using the following reference strains of bacteria at concentrations of 10^4 , 10^5 , 10^6 CFU/ mL:

Escherichia coli ATCC® 25922, *Proteus mirabilis* ATCC® 25933, *Pseudomonas aeruginosa* ATCC® 27853, *Klebsiella pneumoniae* ATCC® 13883, *Staphylococcus aureus* ATCC® 25923, *Enterococcus faecalis* ATCC® 19433, *Candida albicans* ATCC® 10231.

FACTORS THAT MAY INVALIDATE THE RESULTS

Poor standardisation of the inoculum; unsuitable clinical material; use of expired systems or expired supplementary reagents; non compliance with temperatures and times of incubation.

LIMITS AND WARNINGS

For definitive microbial identification of the micro-organisms, confirm the microorganisms using further biochemical confirmatory tests or other identification methods.

PERFORMANCE

The results obtained with the **URIN SYSTEM Plus** agree with those obtained using the URITEST N dip slides (ref. 500232, 51023) and biochemical tests for confirmation of microbial identification. The URITEST N dip slides contain the culture media Cled Agar, MacConkey Agar and Cetrimide Agar.

The results of the susceptibility testing obtained with the system URIN SYSTEM Plus agree with those obtained with the method of radial diffusion according to Bauer et al.⁽⁶⁾ recommended by the Food and Drug Administration (FDA) and from the Clinical and Laboratory Standards Institute (CLSI)⁽¹⁾.

Application of the Wilcoxon non parametric statistical test shows that the two groups of values are not significantly discordant: $P=0.0625$. *Enterococcus* spp. is isolated more frequently with **URIN SYSTEM Plus** than with URITEST N dip slides⁽³⁾.

PRECAUTIONS

The product, **URIN SYSTEM Plus**, is not classified as hazardous under current legislation, nor does it contain harmful substances in concentrations $\geq 1\%$. It therefore does not require a Safety Data Sheet to be available. **URIN SYSTEM Plus** is a disposable device to be used only for diagnostic use *in vitro*. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

CONSERVATION

Store at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.


DISPOSAL OF USED MATERIAL

After use, **URIN SYSTEM Plus** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Ref.	Packaging
URIN SYSTEM Plus	74160	20 tests
URIN SYSTEM Plus	79160	4 tests

TABLE OF SYMBOLS

IVD	<i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF	Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

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This page is intentionally left blank

This page is intentionally left blank

BIBLIOGRAFIA / BIBLIOGRAPHY / BIBLIOGRAPHIE
BIBLIOGRAFIA / BIBLIOGRAFIA / ΒΙΒΛΙΟΓΡΑΦΙΑ

1. CLSI M100-S23 (2013), Performance Standards for antimicrobial Susceptibility Testing.
2. CLSI M7-A9 (2012), Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically.
3. Dati di archivio Liofilchem s.r.l. 'Micro Urin Sistem' (Giugno 2003).
4. Murray, Baron, Pfaller, Tenorev and Tenover: Manual of Clinical Microbiology (1999), 7th Edition.
5. Edwin H.Lenette (1985), Manual of Clinical Microbiology, 4th Edition (ASM).
6. Bauer, A.W., W.M.M. Kirby, J.C. Sherris and M. Turk (1966), Antimicrobic susceptibility testing by a standardized single disk method. AM. J. Clin. Pathol. 45:493-496.

Microbiology Products

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

LIOFILCHEM® s.r.l.

Via Scozia zona ind.le, 64026 Roseto degli Abruzzi (Te) Italy
Tel. +39 0858930745 Fax +39 0858930330 www.liofilchem.net liofilchem@liofilchem.net

F12613
Rev.7 / 24.09.2013