Urin System Chrom

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ENGLISH

Urin System Chrom

System for microbial count, identification and susceptibility testing, directly from urine specimens.

DESCRIPTION

Urin System Chrom is a 32-wells system containing desiccated biochemical, chromogenic and antibiotic substrates for microbial count, identification and susceptibility testing of microorganisms from urine specimens. The system is inoculated with a dilution of urine after the presence of microorganism has been ascertained by microscope examination of urinary sediment and incubated at $36\pm1^{\circ}\text{C}$ for 18-24 hours. The tests results are interpreted by assessing the color change in the various wells.

CONTENT OF THE KIT

Ref. 74161	Ref. 79161
<ul style="list-style-type: none"> • 20 Urin System Chrom • 20 Vials of Suspension Medium (7,0 mL) • 20 Vials of Physiological Solution (4,5 mL) 	<ul style="list-style-type: none"> • 4 Urin System Chrom • 4 Vials of Suspension Medium (7,0 mL) • 4 Vials of Physiological Solution (4,5 mL)

PRINCIPLE OF THE METHOD

Urin System Chrom allows the microbial count, identification and susceptibility testing of pathogenic microorganisms most frequently isolated from urine specimens:

Escherichia coli, *Proteus* spp, *Providencia* spp, *Morganella* spp, *Citrobacter* spp, *Pseudomonas* spp, *Klebsiella* spp, *Enterobacter* spp, *Serratia* spp, *Enterococcus faecalis*, *Staphylococcus aureus*, *Candida albicans*.

Microbial count is evaluated by the microbial growth in the wells **1-GR+** and **2-GR++**.

Identification is obtained with biochemical and chromogenic reactions in the wells **3-ESC** to **13-CAN**.

Susceptibility testing is assessed on the basis of the growth or inhibition of the microorganisms in the wells **14-AK** to **31-SXT** containing antibiotics.

COLLECTION AND CONSERVATION OF THE SAMPLES

Collect aseptically the first urine in the morning in a sterile container.

The collection methods most used are that of *mitto intermedio* for adults who urine on command and that of the adhesive plastic bag in early infancy.

The samples of urine must be sent to the laboratory immediately after collection. If that is not possible, store the sample in the refrigerator at $2\text{-}8^{\circ}\text{C}$ for no more than 24 hours.

The urine must be collected before the start of any antibacterial treatment.

CONFIGURATION

Well	MICROBIAL COUNT		
1-GR+	$10^5 \leq \text{CFU/mL} \leq 10^6$		
2-GR++	$\text{CFU/mL} > 10^6$		
Well	MICROBIAL IDENTIFICATION		
3-ESC	<i>Escherichia coli</i>		
4-IND *	Indole test for confirming <i>Escherichia coli</i>		
5-PD	Phenylalanine deaminase test	<i>Proteus, Providencia, Morganella, Citrobacter</i> spp	
6-UR <input type="checkbox"/>	Urea hydrolysis test		
7-H₂S <input type="checkbox"/>	Hydrogen sulfate production test		
8-PSE	<i>Pseudomonas</i> spp		
9-KES	<i>Klebsiella, Enterobacter, Serratia</i> spp (KES group)		
10-VP *	Voges-Proskauer test for confirming KES group microorganisms		
11-STR	<i>Enterococcus faecalis</i> (D group)		
12-STA	<i>Staphylococcus aureus</i>		
13-CAN	<i>Candida</i> spp		
Well	SUSCEPTIBILITY TESTING ($\mu\text{g/mL}$)		
14-AK	Amikacin	32	
15-CN	Gentamicin	8	
16-TOB	Tobramycin	8	
17-TZP	Piperacillin/Tazobactam	128/4	
18-FOS	Fosfomycin	200	
19-IMI	Imipenem	16	
20-CFM	Cefixime	32	
21-CXM	Cefuroxime	32	
22-CTX	Cefotaxime	64	
23-CAZ	Ceftazidime	32	
24-F	Nitrofurantoin	64	
25-AMS	Ampicillin/Sulbactam	32/16	
26-TE	Tetracycline	16	
27-CIP	Ciprofloxacin	4	
28-LEV	Levofloxacin	8	
29-OFL	Oflloxacin	8	
30-VA	Vancomycin	4	
31-SXT	Co-trimoxazole	8	
32-Growth	Microbial growth control		

: overlay the well with vaseline oil

* : after incubation, add the indicated reagent for performing the test

TEST PROCEDURE

1) PREPARATION OF THE CLINICAL SPECIMEN

- Examine the urinary sediment under the microscope to determine the presence of urinary tract infection.
- Transfer 0.5 mL of urine into a vial of Physiological Solution contained in the kit (**Suspension A**).
- Transfer 0.2 mL of **Suspension A** into a vial of Suspension Medium* contained in the kit (**Suspension B**).

* **Suspension Medium (g/L):** Mueller Hinton Broth 21g, Yeast Extract 5g, Peptone 3g, Glucose 2g, Distilled Water 1000 mL; pH 6.8 ± 0.2

2) INOCULATION OF THE SYSTEM

- Take a system from its wrap and bring it to room temperature.
- Write down the name of the patient and the date of the start of the examination.
- Dispense 0.2 mL of **Suspension A** into the wells from **1-GR+** to **13-CAN (Microbial count and Identification)**.
- Dispense 0.2 mL of **Suspension B** into the wells from **14-AK** to **32-Growth (Susceptibility testing)**.
- Add 1 drop of vaseline oil (**Vaseline Oil Droppers**, ref. 87006) to the wells **6-UR** and **7-H₂S**.
- Cover the system with the lid provided and incubate at 36±1°C for 18-24 hours.
- After incubation watch for the color change in the wells and interpret the results.

Laboratory workflow is illustrated at the end of this document.

INTERPRETATION OF THE RESULTS

IDENTIFICATION

- 1-GR+** **Microbial count** $10^5 \leq \text{CFU/mL} \leq 10^6$: watch for the color change of the well from blue to yellow.
- 2-GR++** **Microbial count** $\text{CFU/mL} > 10^6$: watch for the color change of the well from blue to yellow.
- 3-ESC** ***Escherichia coli***: watch for the color change of the well from colorless to green.
- 4-IND** Confirmation: add 1 drop of Kovac's reagent (**Kovac's Reagent Droppers**, ref. 87001) to the well to carry out the indole test. A red ring indicates a positive reaction.
- 5-PD** ***Proteus* spp**: watch for the color change of the well **5-PD** from yellow to brown, the well **6-UR** from yellow to red-fuchsia and the well **7-H₂S** from yellow to black.
- 6-UR** ***Providencia, Morganella* spp**: watch for the color change of the well **5-PD** from yellow to brown, the well **6-UR** from a yellow to red-fuchsia. The well **7-H₂S** is yellow.
- 7-H₂S** ***Citrobacter* spp**: watch for the color change of the well **7-H₂S** from yellow to black. Wells **5-PD** and **6-UR** are yellow.
- 8-PSE** ***Pseudomonas* spp**: watch for the color change of the well from yellow to turbid green.
Confirmation: take a drop of the well broth culture and perform the oxidase test (**OxidaseTest Stick**, ref. 88029).
- 9-KES** ***Klebsiella, Enterobacter, Serratia* spp**: watch for the color change of the well from colorless to mauve.
- 10-VP** Confirmation: to perform the Voges-Proskauer test, add 2 drops of alpha-naphthol and 1 drop of NaOH or KOH (**VP Reagent Droppers**, ref. 87002 or 87007) to the well. Read within 15-20 minutes, a pink-red color indicates a positive reaction. .
- 11-STR** ***Enterococcus faecalis*** (D group): watch for the color change of the well from yellow to black.
- 12-STA** ***Staphylococcus aureus***, watch for the color change of the well from colorless to mauve.
Confirmation: take a drop of the well broth culture and perform the agglutination test (**Staph Latex Kit**, ref. 96017).
- 13-CAN** ***Candida* spp**, watch for the color change of the well from green to turbid yellow.
Confirmation: take a drop of the well broth culture and examine at the microscope (400x) for the presence of chlamydospores and fungal hyphae.

SUSCEPTIBILITY TESTING

14-AK Watch for the color change of the wells from blue to yellow-gray.

to **31-SXT** Blue = no growth, **Sensitive**.

Yellow-gray = growth, **Resistant**.

32-Growth Well without antibiotics. Any color change from blue and/or the presence of turbidity indicates microbial growth.

Note the results on the TEST RESULTS FORM.

Interpretative table.

Well	MICROBIAL COUNT	Well color	
		Positive reaction	Negative reaction
1-GR+	$10^5 \leq \text{CFU/mL} \leq 10^6$	yellow	blue
2-GR++	$\text{CFU/mL} > 10^6$	yellow	blue

Well	IDENTIFICATION	Well color	
		Positive reaction	Negative reaction
3-ESC	<i>Escherichia coli</i>	green	yellow
4-IND	Indole test (confirmation of <i>Escherichia coli</i>)	red ring	yellow
5-PD	<i>Proteus, Providencia, Morganella</i> spp. Phenylalanine deaminase test	brown-black	yellow
6-UR	<i>Proteus, Providencia, Morganella</i> spp. Urea hydrolysis test	pink-fuchsia	yellow
7-H₂S	<i>Proteus, Citrobacter</i> spp. Hydrogen sulfate production test	brown-black	yellow
8-PSE	<i>Pseudomonas</i> spp	turbid green	yellow
9-KES	KES group (<i>Klebsiella, Enterobacter, Serratia</i> spp.)	mauve	colorless
10-VP	VP test (confirmation of KES group)	pink-red	yellow
11-STR	<i>Enterococcus faecalis</i> (D group)	black	yellow
12-STA	<i>Staphylococcus aureus</i>	mauve	colorless
13-CAN	<i>Candida</i> spp.	yellow	green

Well	SUSCEPTIBILITY TESTING		
	Well color	Microbial growth	Interpretation
14 to 31	blue	inhibited	Sensitive
	yellow-gray	good	Resistant

Well	GROWTH CONTROL	Well color	
		Positive reaction	Negative reaction
32-Growth	Microbial growth control	yellow-gray	blue

QUALITY CONTROL

Each batch of Urin System Chrom is subjected to quality control using reference strains of bacteria at concentrations of 10^4 , 10^5 , 10^6 CFU/ml:

Escherichia coli ATCC® 25922, *Proteus mirabilis* ATCC® 25933, *Klebsiella pneumoniae* ATCC® 13883, *Pseudomonas aeruginosa* ATCC® 27853, *Candida albicans* ATCC® 10231, *Citrobacter freundii* ATCC® 8090, *Enterococcus faecalis* ATCC® 19433, *Staphylococcus aureus* ATCC® 25923, *Enterobacter cloacae* ATCC® 13047.

PERFORMANCE

The results obtained with the Urin System Chrom agree with those obtained using other microbiological and biochemical tests for microbial identification.

The susceptibility testing results obtained with Urin System Chrom agree with those obtained with the broth dilution method recommended by CLSI⁽¹⁾.

FACTORS THAT MAY INVALIDATE THE RESULTS

Poor standardization of the inoculum; unsuitable material to examine; use of expired systems or expired supplementary reagents; non compliance with temperatures and/or times of incubation.

PRECAUTIONS

The product Urin System Chrom does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. Urin System Chrom is a disposable device to be used for *in vitro* diagnostic use only. It is intended for professional use by properly trained personnel in the laboratory, using approved asepsis and safety methods for handling pathogenic.

CONSERVATION

Store Urin System Chrom at 2-8°C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, Urin System Chrom and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Ref.	Packaging
Urin System Chrom	74161	20 tests
Urin System Chrom	79161	4 tests

TABLE OF SYMBOLS

IVD <i>In Vitro Diagnostic medical device</i>		Do not reuse		Manufacturer		Contains sufficient for <n> tests	Temperature limitation
REF Catalogue number		Fragile, handle with care		Use by		Attention, consult accompanying documents	LOT Batch code

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
WORKFLOW

This page is intentionally left blank

BIBLIOGRAPHY

1. CLSI M100-S27, 2017. Performance Standards for Antimicrobial Susceptibility Testing.
2. CLSI M07-A10, 2015. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically.
3. Dati di archivio Liofilchem s.r.l. 'Micro Urin Sistem'. (Giugno 2003).
4. Murray, Baron, Pfaffer, Tenorev and Yolken: Manual of Clinical Microbiology (1999), 7th Edition.
5. Edwin H.Lenette: Manual of Clinical Microbiology (1985), 4th Edition (ASM).
6. Bauer, A.W., W.M.M. Kirby, J.C. Sherris and M. Turk. 1966. Antimicrobic susceptibility testing by a standardized single disk method. AM. J. Clin. Pathol. 45:493-496.

Microbiology Products

Liofilchem® and the Liofilchem company logo are registered trademarks of LIOFILCHEM s.r.l.

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© Copyright LIOFILCHEM 2017