PATHOGENIC SYSTEM

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System for the detection and presumptive identification of pathogenic microorganisms from pharyngeal, cutaneous, auricular and ophthalmic swabs, expectorate and other clinical specimens

DESCRIPTION

PATHOGENIC SYSTEM is a 18 wells system containing desiccated biochemical substrata for the detection and presumptive identification of pathogenic microorganisms and dermatophytes from pharyngeal, cutaneous, auricular and ophthalmic swabs, expectorate, exudate from abscess and injury, and other clinical specimens.

The system is directly inoculated with the suspension of the clinical specimen and incubated at 36 ± 1 °C for 18-24 hours.

The tests for the detection and identification of microorganisms present in the specimen are interpreted by evaluating the color change in the wells and performing immunologic and microscopic confirmation tests.

CONTENT OF THE PACKAGE

The package contains:

20 Systems PATHOGENIC SYSTEM	1 vial containing 20 lyophilized pellet COAGULASE TEST
20 Physiological Solution vials (7.0 mL/vial)	1 TEST RESULTS FORM
1 cartridge containing 20 discs PYR Disc	1 Instruction sheet
1 PYR Reagent (5 mL)	

NECESSARY ITEMS NOT CONTAINED IN THE PACKAGE

STREPTO A CARD (ref.96418)	OXIDASE TEST DISCS (ref.88004)
ENTEROSYSTEM 18R REAGENT (cod.80252) (Vasellin oil, Indole test reagent; VP test reagents)	TRYPTIC SOY BROTH (ref.21104) or BRAIN HEART INFUSION (ref.20104)
	Sundry material for microbiological laboratory

CONFIGURATION

The system has the configuration indicated in table no.1.

Table no.1

Well	PRESUMPTIVE IDENTIFICATION of Group A Streptococcus
1-PYR *	Pyroglutamic-β-naphthylamide test
2-SGA *	Immunochromatographic test for confirmation <i>Streptococcus pyogenes</i> (Group A)
Well	PRESUMPTIVE IDENTIFICATION of Staphylococcus aureus
3-STA*	<i>Staphylococcus aureus</i>
4-COA*	Coagulase test for confirmation of <i>Staphylococcus aureus</i>
Well	PRESUMPTIVE IDENTIFICATION of Haemophilus spp.
5-HAE	<i>Haemophilus spp.</i>
6-X+V	Cultural medium for confirmation of <i>Haemophilus spp.</i>
Well	PRESUMPTIVE IDENTIFICATION of Escherichia coli
7-ESC	<i>Escherichia coli</i>
8-IND *	Indole test for confirmation of <i>Escherichia coli</i>
Well	PRESUMPTIVE IDENTIFICATION of Proteus spp. / Providencia spp.
9-PRO	<i>Proteus spp / Providencia spp</i>
10- <input type="checkbox"/> UR <input type="checkbox"/>	Urea test for confirmation of <i>Proteus spp. / Providencia spp.</i>
Well	PRESUMPTIVE IDENTIFICATION of Pseudomonas spp
11-PSE	<i>Pseudomonas spp</i>
12-OX *	Cytochrome-oxidase test for confirmation of <i>Pseudomonas spp.</i>
Well	PRESUMPTIVE IDENTIFICATION of Group KES microorganisms
13-KES	Group KES (<i>Klebsiella spp, Enterobacter spp, Serratia spp</i>)
14-VP*	Voges Proskauer test for confirmation of Group KES microorganisms
Well	PRESUMPTIVE IDENTIFICATION of Enterococcus faecalis
15-STR	<i>Enterococcus faecalis</i>
16-SB	Cultural medium for confirmation of <i>Enterococcus faecalis</i>
Well	PRESUMPTIVE IDENTIFICATION of Candida spp.
17-CAN	<i>Candida spp.</i>
Well	PRESUMPTIVE IDENTIFICATION of DERMATOPHYTES
18-DTM	Dermatophytes

* : After incubation, add the indicated reagent to perform the test

• : Performe immunochromatographic test

: Add Vaseline Oil

PRINCIPLE OF THE METHOD

PATHOGENIC SYSTEM allows detection and identification of microorganisms frequently isolated from pharyngeal, auricular and ophthalmic swabs, expectorate, exudate from abscess and injury, and other clinical specimens, such as: *Streptococcus pyogenes* (Group A), *Staphylococcus aureus*, *Haemophilus spp.*, *Escherichia coli*, *Proteus spp./Providencia spp.*, *Pseudomonas spp.*, KES group microorganisms, *Enterococcus faecalis*, *Candida spp.* and dermatophytes of the genus *Microsporum spp.*, *Trichophyton spp.*, *Epidermophyton spp.*

- The presence of *Streptococcus pyogenes* is detected by the colorless to pink color change of the well **1-PYR** subsequent to addition of PYR reagent¹ and confirmed through the immunochromatographic test **STREPTO A CARD**² (ref. 96418) carried out from microbial culture of well **2-SGA**.
- The presence of *Staphylococcus aureus* is detected by the red to yellow-orange color change of well **3-STA** and confirmed with the clot presence in the well **4-COA** (Positive Coagulase test).
Strains of *Enterococcus spp.* can cause the red to yellow color change of well **3-STA** without clot presence (Negative Coagulase test).
- The presence of *Haemophilus spp.* is detected by the yellow-orange to red color change of well **5-HAE**⁴ and confirmed through isolation on selective medium (Chocolate Bacitracin Agar, ref. 11023) and biochemical tests from the well **6-X+V**.
- The presence of *Escherichia coli* is detected by the grey-red to blue color change of well **7-ESC** and by the formation of a pink-red ring subsequent to the addition of Kovac's Reagent in the well **8-IND**⁴.
- The presence of *Proteus spp./Providencia spp.* is detected by the yellow to brown-black color change of well **9-PRO**⁴ and confirmed with the yellow to red-fuchsia color change of the well **10-UR**.
- The presence of *Pseudomonas spp.* is detected by the yellow to turbid green color change of well **11-PSE** and confirmed through biochemical test with **OXIDASE TEST DISCS** (ref. 88004) carried out from microbial culture of well **12-OX**⁴.
- The presence of KES group microorganisms (*Klebsiella spp.*, *Enterobacter spp.*, *Serratia spp.*) is detected by the purple to yellow color change of well **13-KES** and by the Voges Proskauer test performed from the well **14-VP**.
The identification of various microorganisms has to be confirmed through isolation on selective media for enterobacteria and biochemical tests⁴.
- The presence of *Enterococcus faecalis* is detected by the yellow to black color change of well **15-STR** and confirmed with the colorless to red color change of the well **16-SB**.
- The presence of *Candida spp.* is detected by the green to yellow color change and by the formation of a sediment in the well **17-CAN**⁵. Watch for the presence of chlamydospores and mycelial hyphae at the microscope.
- The presence of dermatophytes is detected by the yellow to red color change of the well **18-DTM** after the incubation of the system at 36 ± 1 °C for 72 hours. Watch for the presence of spores and mycelial hyphae at the microscope.

COMPOSITION

Table no.2

Well	Content
1- PYR	Culture medium to detect the presence of the enzyme L-pyroglutamyl-aminopeptidase (PYR) of <i>Streptococcus</i> β-hemolytic of Group A
2- SGA	Culture medium to detect <i>Streptococcus</i> β-hemolytic of Group A
3- STA	Culture medium to detect the growth of <i>Staphylococcus aureus</i>
4- COA	Culture medium for coagulase test
5- HAE	Culture medium for selective isolation of <i>Haemophilus spp.</i>
6- H+V	Culture medium to detect the growth of <i>Haemophilus spp.</i>
7- ESC	Culture medium to detect the growth of <i>Escherichia coli</i>
8- IND	Culture medium for the performance of indole test
9- PRO	Culture medium to detect the growth of <i>Proteus spp./Providencia spp.</i>
10- UR	Culture medium for urea hydrolysis
11- PSE	Culture medium to detect the growth of <i>Pseudomonas spp.</i>
12- OX	Culture medium for cytochrome oxidase test
13- KES	Culture medium for the growth of Group KES (<i>Klebsiella spp.</i> , <i>Enterobacter spp.</i> , <i>Serratia spp.</i>) microorganisms
14- VP	Culture medium for Voges Proskauer test
15- STR	Culture medium for selective isolation of <i>Enterococcus faecalis</i>
16- SB	Culture medium for the growth of <i>Enterococcus faecalis</i>
17- CAN	Culture medium for the growth of <i>Candida spp.</i>
18-DTM	Culture medium for the growth of dermatophytes

Physiological Solution (g/L): Sodium chloride 0,9 g; Distilled water 1000.0 mL; pH 7.0 ± 0.2

COLLECTION AND STORAGE OF THE SAMPLE

The clinical materials to be tested from pharyngeal, cutaneous, auricular and ophthalmic swabs, expectorate and other clinical specimens have to be sent to the laboratory for the inoculation of **PATHOGENIC SYSTEM** immediately after the specimen collection. Avoid storage in the fridge even for short time as low temperatures may damage the vitality of some sensitive microorganisms altering the final result.

TEST PROCEDURE

PREPARATION OF THE CLINICAL SPECIMEN

- Subsequent to the collection of the clinical specimen, dip the swab into a tube of nutrient broth (ex. Tryptic Soy Broth ref. 21104, Brain Heart Infusion Broth ref. 20104, Nutrient Broth ref. 20103).
- Leave the swab immersed for 5 minutes.
- Incubate the tube in thermostat at 36 ± 1 °C for 2-4-6 hours.

Note: Inoculating the clinical specimen under examination on a blood agar culture medium is recommended

INOCULATION OF THE SYSTEM

1. Take a system from its pouch and leave it at room temperature.
2. Note name of the patient, date of beginning of the examination and type of clinical material.
3. Transfer 0.5 mL of culture broth of the clinical specimen into a vial of *Physiological Solution* contained in the kit .
4. Put a Pyr Disc, contained in the kit, into the well **1-PYR**.
5. Inoculate 0.2 mL (4 drops) of the Inoculum Solution into each well of **PATHOGENIC SYSTEM**.
6. Cover the well **10-UR** with two drops of VASELINE OIL (ref. 80252)
7. Cover the system with the lid provided and incubate at 36 ± 1 °C for 18-24 hours.
8. Extend the incubation to 72 hours for detection of dermatophytes.

INTERPRETATION OF RESULTS

Subsequently to the incubation:

1. Add 2 drops of PYR Reagent into the well **1-PYR**. Wait 5 minutes and watch for the formation of pink-red colouring of the disc.
2. Perform the rapid immunochromatografic test **STREPTO A CARD**, following the test procedure as in the pack insert (ref.96418), of the microbial culture in the well **2-SGA** to confirm the presence of *Streptococco* β -hemolytic of Group A.
3. Watch for the color change in the well **3-STA**, and interpret the result using the table no.3. In the case of positive reaction (well color yellow -orange) put a lyophilized pellet Coagulase Test, contained into the kit, in the well **4-COA**. Incubate the system for 4-6 hours and watch for the clot formation in the well **4-COA** (coagulase test positive) to confirm the presence of *Staphylococcus aureus*.
4. Watch for the color change in the well **5-HAE**, and interpret the result using the table no.3.
5. Perform the microbial isolation test: take a loop of culture broth from well **6-X+V** and inoculate it on the surface of a selective plate (Chocolate Bacitracin Agar, ref. 11023). Incubate the plate at 36 ± 1 °C in anaerobiosis, and carry out biochemical tests for the cultivated colonies to confirm the presence of *Haemophilus spp.*
6. Watch for the color change in the well **7-ESC**, and interpret the result using the table no.3.
7. Perform indole test adding 2 or 3 drops of Kovac's Reagent (ref. 80270) into the well **8-IND**. The positive reaction is given by the formation within 1 minute of a pink-red ring in the well **8-IND**.
8. Watch for the color change in the well **9-PRO**, **10-UR** and interpret the result using the table no.3.
9. Watch for the color change in the well **11-PSE**, and interpret the result using the table no.3.
10. Perform the oxidase test by adding a disc of **OXIDASE TEST DISCS** (ref.88004) into the well **12-OX** and watch for the formation of a blue-purple coloring within 2 minutes to confirm the presence of *Pseudomonas spp.*
11. Watch for the color change in the well **13-KES**, and interpret the result using the table no.3.
12. Perform the Voges Proskauer Test adding 2 drops of alpha-naphtol and 1 drop of NaOH 40% (ref. 88035) to the well **14-VP** and wait for 15-20 minutes for the appearance of pink-red colour.
13. Watch for the color change of the well **15-STR** and **16-SB**, and interpret the result using the table no.3.
14. Watch for the color change of the well **17-CAN**, and interpret the result using the table no.3. Take a drop from the well **17-CAN**, deposit it on a glass slide, cover it, and watch for chlamydospores and mycelial hyphae at the microscope (40x) to confirm the presence of *Candida spp.*
15. Watch for the color change of the well **18-DTM** and interpret the result using the table no.3 and watch for spores and mycelial hyphae at the microscope (40x) to confirm the presence of dermatophytes.
16. Note the results on the TEST RESULTS FORM (copy as many forms as necessary)

Table no.3

Well	Presumptive identification of <i>Streptococcus</i> Group A	Well color	
		Positive reaction	Negative reaction
1-PYR	Pyroglutamic- β -naphthylamide test	pink	colorless
2-SGA	Immunochromathographic test to confirm <i>Streptococcus pyogenes</i> (Group A)	two red lines	one red lines
Well	Presumptive identification of <i>Staphylococcus aureus</i>	Well color	
		Positive reaction	Negative reaction
3-STA	<i>Staphylococcus aureus</i>	yellow-orange	red
4-COA	Coagulase test to confirm <i>Staphylococcus aureus</i>	clot presence	clot absence
Well	Presumptive identification of <i>Haemophylus spp.</i>	Well color	
		Positive reaction	Negative reaction
5-HAE	<i>Haemophylus spp.</i>	red	yellow-orange
6-X+V	Cultural medium to confirm <i>Haemophylus spp.</i>		
Well	Presumptive identification of <i>Escherichia coli</i>	Well color	
		Positive reaction	Negative reaction
7-ESC	<i>Escherichia coli</i>	blue	grey-red
8-IND	Indole test for confirming <i>Escherichia coli</i>	pink-red	yellow
Well	Presumptive identification of <i>Proteus spp.</i> / <i>Providencia spp.</i>	Well color	
		Positive reaction	Negative reaction
9-PRO	<i>Proteus spp.</i> / <i>Providencia spp.</i>	brown-black	yellow
10-UR	Urea test to confirm <i>Proteus spp.</i> / <i>Providencia spp.</i>	yellow	red-fuchsia
Well	Presumptive identification of <i>Pseudomonas spp.</i>	Well color	
		Positive reaction	Negative reaction
11-PSE	<i>Pseudomonas spp.</i>	turbid green	yellow
12-OX	Cytochrome-oxidase test for confirming <i>Pseudomonas spp.</i>	blue	colorless
Well	Presumptive identification of microorganisms group KES (<i>Klebsiella spp.</i> , <i>Enterobacter spp.</i> , <i>Serratia spp.</i>)	Well color	
		Positive reaction	Negative reaction
13-KES	Group KES (<i>Klebsiella spp.</i> , <i>Enterobacter spp.</i> , <i>Serratia spp.</i>)	yellow	violet
14-VP	Voges Proskauer test to confirm Group KES microorganisms	pink-red	colorless
Well	Presumptive identification of <i>Enterococcus faecalis</i>	Well color	
		Positive reaction	Negative reaction
15-STR	<i>Enterococcus faecalis</i>	brown-black	colorless
16-SB	Cultural medium to confirm <i>Enterococcus faecalis</i>	red	colorless
Well	Presumptive identification of <i>Candida spp.</i>	Well color	
		Positive reaction	Negative reaction
17-CAN	<i>Candida spp.</i>	yellow with sediment	green
	Microscopical examination (40x)	presence of chlamydo spores and mycelial hyphae	absence of chlamydo spores and mycelial hyphae
Well	Presumptive identification of DERMATHOPHYTES	Well color	
		Positive reaction	Negative reaction
18-DTM	Dermatophytes	red	yellow
	Microscopical examination (40x)	presence of spores and mycelial hyphae	absence of spores and mycelial hyphae

Note:

It is possible to carry out the susceptibility test of the various isolated microorganisms, by taking a drop of culture broth from the well and proceeding according to the laboratory methods.

Do not use the wells where supplement reagents are added, like 1-PYR, 8-IND, 10-UR, 12-OX and 14-VP for the susceptibility test.

QUALITY CONTROL

Every batch of **PATHOGENIC SYSTEM** is subjected to the quality control using the following reference microorganisms:

<i>Streptococcus pyogenes</i>	ATCC 19615	<i>Klebsiella pneumoniae</i>	ATCC 13883
<i>Staphylococcus aureus</i>	ATCC 25923	<i>Enterobacter cloacae</i>	ATCC 23355
<i>Haemophilus parainfluenzae</i>	ATCC 7901	<i>Serratia marcescens</i>	ATCC 8100
<i>Escherichia coli</i>	ATCC 25922	<i>Enterococcus faecalis</i>	ATCC 19433
<i>Proteus mirabilis</i>	ATCC 25933	<i>Candida albicans</i>	ATCC 10231
<i>Pseudomonas aeruginosa</i>	ATCC 27853	<i>Trichophyton mentagrophytes</i>	ATCC 9533

FACTORS THAT MAY INVALIDATE THE RESULTS

Poor standardization of the inoculum; clinical material unsuitable; use of expired systems and/or reagents; non compliance with temperatures and times of incubation.

CAUTIONS

The product, **PATHOGENIC SYSTEM** is not classified as hazardous under current legislation, however refer to the safety datasheet for a correct use. **PATHOGENIC SYSTEM** is a disposable device to be used only for diagnostic use *in vitro*. It must be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **PATHOGENIC SYSTEM** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Ref.	Package
PATHOGENIC SYSTEM	71679	20 tests

TABLE OF SYMBOLS

IVD <i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

Rev.2/18.07.2008

BIBLIOGRAFIA / BIBLIOGRAPHY

1. Ellner PD, Williams DA, Hosmer MA, Cohenford MA. Preliminary evaluation of a rapid colorimetric method for the presumptive identification of group A streptococci and enterococci. J. Clin. Microbiol. **22**, 880-881, 1985.
2. Appleton RS, Victoria BE, Tamer D and Ayouf EM. 1985. J.Lab. Clin. Med. **105**:114-119.
3. Philips W, Kloos W. (1981). J.Clin.Microbiol; **14**: 671.
4. Murray, Baron, Pfaller, Tenev. Tenover. Manual of Clinical Microbiology 7th Edition (1999). American Society Microbiology.
5. *Davise H. Larone. Medically important fungi. A Guide to Identification. 2th Edition (1987). Elsevier.*

Bacteriology Products

LIOFILCHEM S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto D.A. (TE) - ITALY

Tel +390858930745 Fax +390858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net