PATHOGENIC SYSTEM AST

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
Prodotto	Ref.	Confezione
PATHOGENIC SYSTEM AST	71681	20 test
PATHOGENIC SYSTEM AST	79681	4 test

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PATHOGENIC SYSTEM AST

ENGLISH

System for the identification and susceptibility testing of the main pathogenic microorganisms directly from several clinical specimens.

DESCRIPTION

PATHOGENIC SYSTEM AST is a 32-well system containing desiccated biochemical and antibiotic substrates for the detection, identification and susceptibility testing of pathogenic microorganisms from pharyngeal, cutaneous, auricular and ophthalmic swabs, expectorate, exudates and other clinical specimens.

The system is directly inoculated with the suspension of the clinical specimen and incubated at $36\pm 1^\circ\text{C}$ for 18-24 hours.

The tests are interpreted by assessing the colour change of the various wells and performing immunoserological and microscope confirmation tests.

CONTENT OF THE PACKAGE

Ref. 71681	Ref. 79681
<ul style="list-style-type: none"> • 20 PATHOGENIC SYSTEM AST • 20 Vials of Physiological Solution (7.0 mL) • 20 Vials of Inoculum Broth 2 (7.0 mL) • 1 Vial of Reagent A (2.0 mL) • 1 Vial of Reagent B (2.0 mL) • 1 Cartridge containing 20 PYR-Disc 	<ul style="list-style-type: none"> • 4 PATHOGENIC SYSTEM AST • 4 Vials of Physiological Solution (7.0 mL) • 4 Vials of Inoculum Broth 2 (7.0 mL) • 1 Vial of Reagent A (2.0 mL) • 1 Vial of Reagent B (2.0 mL) • 1 Cartridge containing 4 PYR-Disc

PRINCIPLE OF THE METHOD

PATHOGENIC SYSTEM AST allows the detection, identification and susceptibility testing of the microorganisms most frequently isolated from clinical specimens, such as:

Streptococcus pyogenes (Group A), *Staphylococcus aureus*, *Haemophylus* spp, *Escherichia coli*, *Streptococcus agalactiae* (Group B), *Proteus* spp/*Providencia* spp, *Pseudomonas* spp, KES Group (*Klebsiella* spp, *Enterobacter* spp, *Serratia* spp), *Enterococcus faecalis* (Group D), *Candida* spp .

The identification is showed based by using biochemical and cromogenic substrates present in the various wells of the system.

The susceptibility testing is assessed on the basis of either the microorganisms growth or inhibition in the wells containing the antibiotic and a growth indicator.

COLLECTION AND STORAGE OF THE SAMPLE

The clinical materials to examined should be sent to the laboratory using swabs with transport media (BIO-CULT AMIES, ref. 85307; BIO-CULT STUART, ref. 85308).

Avoid storage in the fridge for a long time as low temperatures may damage the vitality of some sensitive microorganisms altering the final result.

CONFIGURATION

Well	IDENTIFICATION	
1-PYR *	<i>Streptococcus pyogenes</i> (group A) - PYR test (Pirrolidonil- β -naftilammide hydrolysis)	
2-SGA •	Immunoserological test for confirming <i>Streptococcus pyogenes</i>	
3-STA	<i>Staphylococcus aureus</i>	
4-HAE	<i>Haemophylus</i> spp	
5-NES *	<i>Neisseria</i> spp	
6-STG	<i>Streptococcus agalactiae</i> (group B)	
7-STR	<i>Enterococcus faecalis</i> (group D)	
8-PSE	<i>Pseudomonas</i> spp	
9-PRO	<i>Proteus</i> spp / <i>Providencia</i> spp	
10-ESC	<i>Escherichia coli</i>	
11-IND *	Indole test for confirming <i>Escherichia coli</i>	
12-KES	<i>Klebsiella</i> , <i>Enterobacter</i> , <i>Serratia</i> (KES group)	
13-VP *	Voges-Proskauer test for confirming <i>Klebsiella</i> , <i>Enterobacter</i> , <i>Serratia</i>	
14-CAN	<i>Candida</i> spp	
Well	SUSCEPTIBILITY TESTING	Concentration ($\mu\text{g/mL}$)
15-CN	Gentamicin	16
16-NET	Netilmicin	32
17-AZM	Azithromycin	2
18-CLA	Clarithromycin	8
19-AUG	Amoxicillin / Clavulanic Acid	32/16
20-AMS	Ampicillin / Sulbactam	32/16
21-TZP	Piperacillin / Tazobactam	128/4
22-CRO	Ceftriaxone	64
23-CFP	Cefoperazone	64
24-CAZ	Ceftazidime	32
25-CXM	Cefuroxime	32
26-PEF	Pefloxacin	8
27-LEV	Levofloxacin	8
28-OFX	Ofloxacin	8
29-RD	Rifampicin	4
30-MRP	Meropenem	16
31-SXT	Co-trimoxazole	16
32-Growth	Growth control for the susceptibility tests	

* : After incubation, add the indicated reagents to carry out the test

• : After incubation, perform the immunoserological agglutination test

TEST PROCEDURE**WORKFLOW****PREPARATION OF THE CLINICAL SPECIMEN****1. SWAB: Pharyngeal, Cutaneous, Auricular, Ophthalmic**

Dip the swab into a tube containing 3 mL of nutrient broth*

2. LIQUID SAMPLE: Expectorate, Exudate, Blood culture

Dispense 0.2 mL specimen in a tube containing 3 mL nutrient broth*.

- Incubate the broth culture in thermostat at $36\pm 1^\circ\text{C}$ for 2-3 hours.
- Transfer 0.2 mL of broth culture into a vial of Physiological Solution (**Suspension A**).
- Transfer 0.01 mL of the broth culture into a vial of Inoculum Broth 2** (**Suspension B**).


* Use one of the follow nutrient broths:
Tryptic Soy Broth (ref. 24513),
Brain Heart Infusion Broth (ref. 24104),
Nutrient Broth (ref. 24103),
Todd Hewitt Broth (ref. 24111).

** Inoculum Broth 2 (g/L):
Glucose 15g; Tryptone 2.5g; Sodium Chloride 2.5g;
Yeast Extract 2.5g; Agar 0.5g;
Distilled water 1000 mL; pH 7.0 ± 0.2

INOCULATION OF THE SYSTEM

- Take a system from its wrapper and bring it to room temperature.
- Write down the name of the patient, the date of start of the examination and the type of clinical material.
- Put a **PYR-Disc** into the well **1-PYR**.
- Transfer the **Suspension A** into a reservoir solution (ref. 96761) and using a multichannel or single channel pipette with sterile tips, dispense 0.2 mL of Suspension A into the wells from **1-PYR** to **14-CAN (Identification)**.
- Transfer the **Suspension B** into a reservoir solution and using a multichannel or single channel pipette with sterile tips, dispense 0.2 mL of Suspension B into the wells from **15-CN** to **32-Growth (Antimicrobial Susceptibility Testing)**.
- Cover the system with the lid provided and incubate at $36\pm 1^\circ\text{C}$ for 18-24 hours.
- After incubation, observe the color change in the wells and interpret the results.

Note: Inoculating the clinical specimen under examination on a blood agar culture medium is recommended (Columbia Agar, ref. 11025).

INTERPRETATION OF THE RESULTS

- 1-PYR** *Streptococcus pyogenes* (group **A**): add 1 drop of **Reagent A** into the well and incubate at 36±1°C for 10 minutes. After incubation add 1 drop of **Reagent B** into the well and wait 3 minutes.
Disc turning pink-fuchsia indicates the microorganism presence.
- 2-SGA** *Streptococcus pyogenes* confirmation: take a drop of the well broth culture and perform the **STREPTO A latex kit** agglutination test (ref. 96154).
- 3-STA** *Staphylococcus aureus*: watch the color change of the well from red to yellow-orange.
Confirmation: take a drop of the well broth culture and perform the **STAPH latex kit** agglutination test (ref. 96017).
- 4-HAE** *Haemophilus spp*: watch for the color change of the well from yellow-orange to red.
Confirmation: take a drop of well culture broth and inoculate onto a selective plate of Chocolate Bacitracin Agar (ref. 11023).
- 5-NES** *Neisseria spp*: take a drop of the well broth culture and put it onto a **OXIDASE TEST STICK** (ref. 88029); watch for the development of a blue color in about 2 minutes (positive test).
Confirmation: take a drop of the well broth culture and inoculate onto a selective plate of Thayer Martin Agar (ref.11040).
- 6-STG** *Streptococcus agalactiae* (group **B**): watch for the color change of the well from colorless to green.
Confirmation: take a drop of the well broth culture and perform the **STREPTO B latex kit** agglutination test (ref. 96153).
- 7-STR** *Enterococcus faecalis* (group **D**): watch for the color change of the well from colorless to black.
- 8-PSE** *Pseudomonas spp*: watch for the color change of the well from yellow to turbid green.
- 9-PRO** *Proteus spp/Providencia spp*: watch for the color change of the well from yellow to brown-black.
- 10-ESC** *Escherichia coli*: watch for the color change of the well from colorless to green.
- 11-IND** Confirmation: perform the indole test adding 1 drop KOVAC'S *Reagent* (ref. 80271) into the well and wait for the appearance in about 1 minute of a red ring (positive test).
- 12-KES** **KES group** (*Klebsiella spp, Enterobacter spp, Serratia spp*): watch for the color change of the well from colorless to mauve.
- 13-VP** Confirmation: perform the Voges-Proskauer test adding the reagents (ref. 87002) and wait for a red color development in about 15 minutes.
- 14-CAN** *Candida spp*: watch for the color change of the well from green to turbid yellow.
Confirmation: take a drop of well broth culture and examine at the microscope (400X) for the presence of chlamydospores and fungal hyphae.

SUSCEPTIBILITY TESTING

- 15-31** Susceptibility testing of glucose fermenting microorganisms.
Red = no growth, **Sensitive** ; yellow-orange = growth, **Resistant**.
Susceptibility testing of glucose non-fermenting microorganisms (*Pseudomonas spp* and other species).
Red = no growth, **Sensitive** ; turbid red-fuchsia = growth, **Resistant**.
- 32-Growth** Well broth culture without antibiotics, turning yellow indicates microbial growth.

Note the results on the TEST RESULT FORM included as appendix.

Table 1.

Well	IDENTIFICATION	Well color	
		Positive reaction	Negative reaction
1-PYR	<i>Streptococcus pyogenes</i> (group A) - PYR test	pink-fuchsia disc	light pink disc
2-SGA	Immunoserological test (<i>S. pyogenes</i> confirmation)	agglutination	no agglutination
3-STA	<i>Staphylococcus aureus</i>	yellow-orange	red
4-HAE	<i>Haemophilus</i> spp	red	yellow-orange
5-NES	<i>Neisseria</i> spp (oxidase test)	blue stick	colorless stick
6-STG	<i>Streptococcus agalactiae</i> (group B)	green	colorless
7-STR	<i>Enterococcus faecalis</i> (group D)	brown-black	colorless
8-PSE	<i>Pseudomonas</i> spp	turbid green	yellow-orange
9-PRO	<i>Proteus</i> spp / <i>Providencia</i> spp	brown-black	yellow
10-ESC	<i>Escherichia coli</i>	green	colorless
11-IND	Indole test (<i>Escherichia coli</i> confirmation)	red ring	colorless
12-KES	KES group (<i>Klebsiella</i> spp, <i>Enterobacter</i> spp, <i>Serratia</i> spp)	mauve	colorless
13-VP	VP test (KES group confirmation)	red	colorless
14-CAN	<i>Candida</i> spp	turbid yellow	green

well	SUSCEPTIBILITY TESTING	Well color	
		Positive reaction	Negative reaction
15 to 31	Glucose fermenting microorganisms	yellow-orange	red
	Glucose non-fermenting microorganisms	turbid red-fuchsia	red

32-Growth	Microbial growth control	yellow	red
------------------	--------------------------	--------	-----

QUALITY CONTROL

PATHOGENIC SYSTEM AST is subjected to the quality control using the following reference microorganisms:

Streptococcus pyogenes ATCC® 19615, *Staphylococcus aureus* ATCC® 25923, *Haemophilus parainfluenzae* ATCC® 7901, *Escherichia coli* ATCC® 25922, *Proteus mirabilis* ATCC® 25933, *Pseudomonas aeruginosa* ATCC® 27853, *Klebsiella pneumoniae* ATCC® 13883, *Streptococcus agalactiae* ATCC® 13813, *Neisseria gonorrhoeae* ATCC® 19424, *Enterococcus faecalis* ATCC® 19433, *Candida albicans* ATCC® 10231.

FACTORS THAT MAY INVALIDATE THE RESULTS

Poor standardization of the inoculum; clinical material unsuitable; use of expired systems or expired supplementary reagents; non compliance with temperatures and times of incubation.

LIMITS AND WARNINGS

For definitive identification of the microorganisms, it is necessary to make use of biochemical confirmation tests.

PRECAUTIONS

The product, **PATHOGENIC SYSTEM AST** is not classified as hazardous under current legislation, however refer to the safety data sheet for a correct use. **PATHOGENIC SYSTEM AST** is a disposable device to be used only for diagnostic use *in vitro*. The product must be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE AND TRANSPORT CONDITIONS

Store at 2-8°C in the original packaging. Stability studies have shown that the storage or transport at 18-25°C for 4 days, or at 35-39°C for 48 hours, do not alter in any way the performance of the product. Do not use beyond the expiry date indicated on the label.

DISPOSAL OF USED MATERIAL

After use, **PATHOGENIC SYSTEM AST** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Ref.	Package
PATHOGENIC SYSTEM AST	71681	20 tests
PATHOGENIC SYSTEM AST	79681	4 tests

TABLE OF SYMBOLS

IVD for <i>in vitro</i> diagnostic use	 Do not reuse	 Manufacturer	 Contains sufficient for <n> test	 Temperature limits
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch number

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BIBLIOGRAFIA / BIBLIOGRAPHY

- Ellner PD, Williams DA, Hosmer MA, Cohenford MA. Preliminary evaluation of a rapid colorimetric method for the presumptive identification of group A streptococci and enterococci. *J. Clin. Microbiol.* **22**, 880-881, 1985.
- Appleton RS, Victoria BE, Tamer D and Ayouf EM. 1985. *J.Lab. Clin. Med.* **105**:114-119.
- Philips W, Kloos W. (1981). *J.Clin.Microbiol*; **14**: 671.
- Murray, Baron, Pfaller, Tenorev. Tenover. *Manual of Clinical Microbiology* 7th Edition (1999). American Society Microbioloy.
- *Davise H. Larone. Medically important fungi. A Guide to Identification. 2th Edition (1987). Elsevier.*
- CLSI. Performance Standards for Antimicrobial Susceptibility Testing; Twenty-Fourth Informational Supplement. CLSI document M100-S24. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.
- CLSI. Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria that Grow Aerobically. CLSI document M7-A9. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.
- The European Committee on Antimicrobial Susceptibility Testing. Breakpoint tables for interpretation of MICs and zone diameters. Version 4.0, 2014.

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