Integral System Streptococchi

ITALIANO

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Integral System Streptococci

System for the identification and susceptibility testing of Streptococci.

ENGLISH

DESCRIPTION

Integral System Streptococci is a 24-well system containing desiccated biochemical substrata for the identification and susceptibility testing of streptococci.

The system is inoculated with the bacterial suspension of the microorganism to be examined and incubated at $36\pm 1^\circ\text{C}$ for 18-24 hours.

At the end of the incubation period the color change of the various wells allows to identify the organism by its numerical profile and test susceptibility to 7 different antibiotics.

CONTENT OF THE PACKAGE

Ref. 71720	Ref. 79720
<ul style="list-style-type: none"> • 20 Systems Integral System Streptococci • 20 Vial of Inoculum Broth (7.0 mL) • 1 Cartridge of PYR-Disc • 1 Vial of PYR Reagent (5.0 mL) • 1 Vial of Ninhydrin Reagent (5.0 mL) 	<ul style="list-style-type: none"> • 4 Systems Integral System Streptococci • 4 Vial of Inoculum Broth (7.0 mL) • 1 Cartridge of PYR-Disc • 1 Vial of PYR Reagent (5.0 mL) • 1 Vial of Ninhydrin Reagent (5.0 mL)

NECESSARY ITEMS NOT CONTAINED IN THE PACKAGE

<ul style="list-style-type: none"> • Physiological Solution (ref. 20095).
--

PRINCIPLE OF THE METHOD

Integral System Streptococci allows to simultaneously perform identification and susceptibility testing of the organism to be examined.

Identification is obtained by assessing the combination of positive and negative reactions in the wells **1-PYR** to **9-OPT** for determining the numerical profile. Identification is completed by referring to the Codes table or using the Software on-line.

Susceptibility testing is assessed on the basis of the growth or inhibition of the microorganism in the wells **10-SXT** to **24-PRL**. 7 antibiotics in double concentration are evaluated.

CONFIGURATION

Well	STREPTOCOCCI IDENTIFICATION	
1-PYR	*	Hydrolysis of pirrolidonil- β -naftilammide (PYR test)
2-ESC		Hydrolysis of esculin
3-IPP	*	Hydrolysis of hippurate
4-ONPG		Hydrolysis of ONPG (ortho-nitrophenyl- β -galactoside)
5-ARA		Fermentation of arabinose
6-MAN		Fermentation of mannitol
7-RAF		Fermentation of raffinose
8-BAC		Susceptibility to bacitracin
9-OPT		Susceptibility to optochine
Well	SUSCEPTIBILITY TESTING	
10-SXT	Co-trimoxazole	8 $\mu\text{g/ml}$
11-ERY	Erythromycin	1 $\mu\text{g/ml}$
12-ERY	Erythromycin	4 $\mu\text{g/ml}$
13-PEF	Pefloxacin	1 $\mu\text{g/ml}$
14-PEF	Pefloxacin	4 $\mu\text{g/ml}$
15-KF	Cephalothin	8 $\mu\text{g/ml}$
16-KF	Cephalothin	32 $\mu\text{g/ml}$
17-TEC	Teicoplanin	4 $\mu\text{g/ml}$
18-TEC	Teicoplanin	16 $\mu\text{g/ml}$
19-RD	Rifampicin	4 $\mu\text{g/ml}$
20-RD	Rifampicin	16 $\mu\text{g/ml}$
21-AMP	Ampicillin-sulbactam	8-4 $\mu\text{g/ml}$
22-AMP	Ampicillin-sulbactam	16-8 $\mu\text{g/ml}$
23-PRL	Piperacillin	16 $\mu\text{g/ml}$
24-PRL	Piperacillin	64 $\mu\text{g/ml}$

* : after incubation, add the indicated reagent

TEST PROCEDURE**PREPARATION OF THE BACTERIAL SUSPENSION**

1. The microorganism to be identified must be recently isolated (18-24 h); bacteria cultures older than 48 hours can provide not reliable results.
2. Take one or more morphologically similar, well-isolated colonies from the agar culture medium and suspend in Physiological Solution (**Suspension A**).

INOCULUM OF THE SYSTEM

1. Take a system from its wrapper and bring it to room temperature.
2. Write down patient's name, date and origin of the microorganism.
3. Transfer 0.2 ml of Suspension A into wells **1-PYR** to **7-RAF**.
4. Deposit a **PYR-Disc** into the well **1-PYR**.
5. Transfer 50 µL of Suspension A in Inoculum Broth contained in the packaging (**Suspension B**).
6. Transfer 0.2 mL of Suspension B into wells **8-BAC** to **24-PRL**.
7. Cover the system with the lid provided and incubate at 36±1°C for 18-24 hours.

INTERPRETATION OF RESULTS

At the end of the incubation period:

1. Add 2 drops of PYR Reagent into the well **1-PYR**;
2. Add 2 drops of Ninhydrin Reagent into the well **3-IPP**;
3. Wait for 15-20 minutes and watch for the color change in the wells. Interpret the results using the interpretative table.
4. Note the results and determine the 3-digit code by following instructions provided in the NUMERICAL CODE FORMATION paragraph.
5. Identify the microorganism by referring to the Codes table or using the Software on-line.

Interpretative table.

Well	STREPTOCOCCI IDENTIFICATION	Well color	
		Positive reaction	Negative reaction
1-PYR	PYR test	pink-fuchsia disc	light pink disc
2-ESC	Esculin hydrolysis	black-brown	yellow
3-IPP	Hippurate hydrolysis	blue-violet	colorless
4-ONPG	ONPG hydrolysis	yellow	colorless
5-ARA	Arabinose fermentation	yellow	violet
6-MAN	Mannitol fermentation	yellow	violet
7-RAF	Raffinose fermentation	yellow	violet
8-BAC	Bacitracin test	red	yellow
9-OPT	Optochine test	red	yellow

Well	SUSCEPTIBILITY TESTING		
	Well color with low antibiotic concentration	Well color with high antibiotic concentration	Interpretation
10 to 24	red	red	Sensitive
	yellow	red	Intermediate sensitivity
	yellow	yellow	Resistant

NUMERICAL CODE FORMATION

The biochemical tests are separated into 3 groups of 3 and a value of 1, 2 or 4 is indicated for each:

- Value 1 : first test positive in each group (**PYR, ONPG, RAF**).
- Value 2 : second test positive in each group (**ESC, ARA, BAC**).
- Value 4 : third test positive in each group (**IPP, MAN, OPT**).
- Value 0 : every negative test

A 3-digit code is obtained by adding together the values corresponding to positive reactions within each group. The code allows the identification of the organism under examination by referring to the Codes table. The example below shows how a numerical code can be formed.

Example.

Test	Group 1			Group 2			Group 3		
	PYR	ESC	IPP	ONPG	ARA	MAN	RAF	BAC	OPT
Values	1	2	4	1	2	4	1	2	4
Results	+	-	-	-	-	-	-	+	-
Sum of values	1			0			2		
CODE: 102 IDENTIFICATION: <i>Streptococcus pyogenes</i>									

SCHEME OF REACTIONS

Microorganism	PYR	ESC	IPP	ONPG	ARA	MAN	RAF	BAC	OPT	Hemolysis
<i>S. pyogenes</i>	+	-	-	±	-	-	-	+	-	β
<i>S. agalactiae</i>	-	-	+	-	-	-	-	-	-	α/β
<i>E. faecalis</i>	+	+	±	±	-	+	-	-	-	-/α/β
<i>E. faecium</i>	±	+	±	±	+	+	-	-	-	-/α
<i>S. durans</i>	+	+	-	±	-	-	-	-	-	-/α/β
<i>S. mutans</i>	-	±	-	-	-	+	±	-	-	-/β
<i>S. uberis</i>	-	±	±	±	-	+	±	-	-	-/α
<i>S. sanguis I</i>	-	±	-	-	-	-	±	-	-	α
<i>S. salivarius</i>	-	±	-	±	-	-	±	-	-	-
<i>S. MG intermedius</i>	-	+	-	±	-	-	±	-	-	-/α
<i>S. sanguis II</i>	-	-	-	±	-	-	+	-	-	α
<i>S. mitis</i>	-	-	-	±	-	-	-	-	-	α
<i>S. anginosus constellatus</i>	-	±	-	-	-	-	±	-	-	α
<i>S. morbillorum</i>	-	-	-	-	-	-	-	-	-	-
<i>S. acidominimus</i>	-	-	+	-	-	-	-	-	-	α
<i>S. pneumoniae</i>	-	-	-	+	-	-	+	-	+	α

+ : positive reaction

- : negative reaction

± : variable reaction

Codes table.

Code	Microorganism	Code	Microorganism	Code	Microorganism
000	<i>S.sanguis</i> I <i>S.salivarius</i> <i>S.mitis</i> <i>S.anginosus constellatus</i> <i>S.morbilorum</i>	200	<i>S.sanguis</i> I <i>S.salivarius</i> <i>S. MG intermedius</i> <i>S.anginosus constellatus</i>	400	<i>S.acidominimus</i> CAMP Test S negativo
001	<i>S.sanguis</i> I <i>S.salivarius</i> <i>S.sanguis</i> II <i>S.anginosus constellatus</i>	201	<i>S.sanguis</i> I <i>S.salivarius</i> <i>S. MG intermedius</i> <i>S.anginosus constellatus</i>		<i>S.agalactiae</i> CAMP Test S positivo
010	<i>S.salivarius</i> <i>S.mitis</i>	211	<i>S.salivarius</i> <i>S. MG intermedius</i>	441	<i>S.uberis</i>
011	<i>S.sanguis</i> II <i>S.salivarius</i>	240	<i>S.mutans</i> <i>S.uberis</i>	450	<i>S.uberis</i>
015	<i>S.pneumoniae</i>	241	<i>S.mutans</i> <i>S.uberis</i>	451	<i>S.uberis</i>
040	<i>S.uberis</i> <i>S.mutans</i>	251	<i>S.uberis</i>	640	<i>S.uberis</i>
041	<i>S.uberis</i> <i>S.mutans</i>	260	<i>E.faecium</i>	641	<i>S.uberis</i>
050	<i>S.uberis</i>	270	<i>E.faecium</i>	650	<i>S.uberis</i>
051	<i>S.uberis</i>	300	<i>S.durans</i>	651	<i>S.uberis</i>
102	<i>S.pyogenes</i>	340	<i>E.faecalis</i>	660	<i>E.faecium</i>
112	<i>S.pyogenes</i>	350	<i>E.faecalis</i>	670	<i>E.faecium</i>
		360	<i>E.faecium</i>	740	<i>E.faecalis</i>
		370	<i>E.faecium</i>	750	<i>E.faecalis</i>
				760	<i>E.faecium</i>

QUALITY CONTROL

Every batch of **Integral System Streptococchi** is subjected to the quality control using reference microbial strains: *Streptococcus pyogenes* ATCC® 19615, *Streptococcus pneumoniae* ATCC® 6303, *Enterococcus faecalis* ATCC® 19433, *Streptococcus agalactiae* ATCC® 13813.

FACTORS THAT MAY INVALIDATE THE RESULTS

Poor standardization of the inoculum; unsuitable clinical material; use of expired systems and/or reagents; non compliance with temperatures and times of incubation.

CAUTIONS

The product, **Integral System Streptococchi**, is classified as hazardous under current legislation, refer to the safety data sheet for a correct use. **Integral System Streptococchi** is a disposable device to be used only for diagnostic use *in vitro*. It must be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8°C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **Integral System Streptococchi** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

BIBLIOGRAPHY

- Ronald M. Atlas. Handbook of Microbiological Media, 4th Edition (2010) CRC Press.
- Patrick R. Murray, Ellen Jo Baron, James H. Jorgensen, Michael A. Pfaller, and Robert H. Tenover. Manual of Clinical Microbiology, 8th Edition (2003) ASM Press.
- Edwin H. Lennette. Manual of Clinical Microbiology, 4th Edition (1985) ASM Press.

PRESENTATION

Product	Ref.	Packaging
Integral System Streptococchi	71720	20 Tests
Integral System Streptococchi	79720	4 Tests

TABLE OF SYMBOLS

IVD <i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

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