INTEGRAL SYSTEM

YEASTS *Plus*

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Prodotto	Ref.	Confezione
INTEGRAL SYSTEM YEASTS <i>Plus</i>	71822	20 test
INTEGRAL SYSTEM YEASTS <i>Plus</i>	79822	4 test

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INTEGRAL SYSTEM YEASTS *Plus*

DESCRIPTION

INTEGRAL SYSTEM YEASTS *Plus* is a 24 wells system containing biochemical substrata and dried antimycotics for the identification of the most clinically important yeasts and sensitivity evaluation to antimycotics. The system is inoculated with the cell suspension and incubated at $36 \pm 1^\circ\text{C}$ for 48 hours.

KIT CONTENTS

The kit contains:

Ref. 71822	20 Systems INTEGRAL SYSTEM YEASTS <i>Plus</i>	40 Vials of <i>Physiological Solution</i> (4.5 mL/vial)
Ref. 79822	4 Systems INTEGRAL SYSTEM YEASTS <i>Plus</i>	8 Vials of <i>Physiological Solution</i> (4.5 mL/vial)
1 Cartridge of Xylose Disc - 1 Instructions Sheet - 1 Test Results Form		

ITEMS NECESSARY BUT NOT INCLUDED IN THE KIT

Vaseline oil for microbiological use (VASELINE OIL: 2 bottles x 50 mL, ref.80278)

Various materials for laboratory of microbiology

CONFIGURATION

The configuration of the system is shown in Table no.1:

Table no.1

Well	ASSIMILATION	Well	CHROMATIC REACTION
1-GLU	Glucose	13-CHR	Chromogenic substrate
2-MAL	Maltose	Well	ANTIMYCOTICS
3-SAC	Saccarose	14-NY	Nistatin 1.25 µg/mL
4-LAC	Lactose	15-AMB	Amphotericin 2 µg/mL
5-GAL	Galactose	16-FCY	Flucytosine 16 µg/mL
6-MEL	Melibiose	17-ECN	Econazole 2 µg/mL
7-CEL	Cellobiose	18-KCA	Ketoconazole 0.5 µg/mL
8-INO	Inositol	19-CLO	Clotrimazole 1 µg/mL
9-XYL	Xylose	20-MIC	Miconazole 2 µg/mL
10-RAF	Raffinose	21-ITR	Itraconazole 1 µg/mL
11-TRE	Trehalose	22-VOR	Voriconazole 2 µg/mL
12-DUL	Dulcitol	23-FLU	Fluconazole 64 µg/mL
		24-Growth	Growth control

PRINCIPLE OF THE METHOD

INTEGRAL SYSTEM YEASTS *Plus* allows:

- **Presumptive Identification.** It is based on assimilation reactions of sugars; the tests for the assimilation reactions are interpreted by evaluating the color change of wells **1-GLU** to **12-DUL**. The combination of positive and negative reactions allows the formation of a numerical code which permits to identify the yeasts under examination, through the use of the **TABLE OF CODES**.
The well **13-CHR** contains a chromogenic substrate that permits to differentiate some yeasts by evaluating the color change of the well.
- **Sensitivity to antimycotics.** The tests are evaluated according to growth or inhibition of yeasts in media containing the antimycotic and a growth indicator in the wells **14-NY** to **23-FLU**.
The color change from red to orange in the wells indicates a slow growth of the yeast under examination and an intermediate sensitivity to the concentration of antimycotic in the well.
The color change from red to yellow in the wells indicates a growth of the yeast under examination and resistance to the concentration of antimycotic in the well.
No color change in the well indicates no growth of the yeast under examination and sensitivity to the concentration of antimycotic in the well.
The well **24-Growth** does not contain antimycotics, it contains culture medium and indicator and it works as growth control.

COMPOSITION

Table no.2

Well	Content
1-GLU	Culture medium for the assimilation of glucose
2-MAL	Culture medium for the assimilation of maltose
3-SAC	Culture medium for the assimilation of saccharose
4-LAC	Culture medium for the assimilation of lactose
5-GAL	Culture medium for the assimilation of galactose
6-MEL	Culture medium for the assimilation of melibiose
7-CEL	Culture medium for the assimilation of cellobiose
8-INO	Culture medium for the assimilation of inositol
9-XYL	Culture medium for the assimilation of xylose
10-RAF	Culture medium for the assimilation of raffinose
11-TRE	Culture medium for the assimilation of trehalose
12-DUL	Culture medium for the assimilation of dulcitol
13-CHR	Culture medium containing chromogenic substratum
14-NY	Culture medium containing Nistatin - 1.25 µg/mL
15-AMB	Culture medium containing Amphotericin - 2 µg/mL
16-FCY	Culture medium containing Flucytosine - 16 µg/mL
17-ECN	Culture medium containing Econazole - 2 µg/mL
18-KCA	Culture medium containing Ketaconazole – 0.5 µg/mL
19-CLO	Culture medium containing Clotrimazole - 1 µg/mL
20-MIC	Culture medium containing Miconazole - 2 µg/mL
21-ITR	Culture medium containing Itraconazole - 1 µg/mL
22-VOR	Culture medium containing Voriconazole - 2 µg/mL
23-FLU	Culture medium containing Fluconazole - 64 µg/mL
24-Growth	Culture medium for the growth control
Physiological Solution (g/L): Sodium chloride 9 g; Distilled water 1000 mL; pH 6.8 ± 0.2	

COLLECTION AND STORAGE OF THE SAMPLES

INTEGRAL SYSTEM YEASTS *Plus* is used for the identification of the most clinically important yeasts and sensitivity evaluation to antimycotics isolated on selective culture media for the isolation of yeasts.

TEST PROCEDURE**1) PREPARATION OF THE CELL SUSPENSION**

- Take one or more microbial colonies and suspend in a *Physiological Solution* vial contained in the kit till obtaining a 0.5 Mc Farland turbidity (ref. 80400) (**Suspension A**).
- Dispense 0.02 mL of **Suspension A** into another *Physiological Solution* vial contained in the kit (**Suspension B**).

2) INOCULATION OF THE SYSTEM

- Take a system from its wrapper and bring it to room temperature.
- Take note of name, date and origin of the microorganism.
- Insert one xylose disc in the well **9-XYL**.
- Transfer 0.2 mL (4 drops) of **Suspension A** into wells **1-GLU** to **13-CHR**.
- Transfer 0.2 mL (4 drops) of **Suspension B** into wells **14-NY** to **24-Growth**.
- Cover all wells of the system with one drop of vaseline oil except **13-CHR**.
- Cover the system with the provided lid and incubate at 36 ± 1 °C for 48 hours.

INTERPRETATION OF THE RESULTS**IDENTIFICATION**

- Observe the color change in the wells **1-GLU** to **12-DUL** and interpret the results using the table no.3.
- Write down the obtained results on the **Test Result Form** (copy as many forms as necessary).
- Obtain the 4 digit code following the instructions reported in the paragraph **NUMERICAL CODE FORMATION**.
- Finally, identify the microorganism using the **Table of codes**.
- Observe the color change of the well **13-CHR**.

Table no.3

Well	IDENTICATION OF YEASTS	Well color	
		Positive reaction	Negative reaction
1-GLU	glucose	yellow-grey	purple
2-MAL	maltose	yellow-grey	purple
3-SAC	saccharose	yellow-grey	purple
4-LAC	lactose	yellow-grey	purple
5-GAL	galactose	yellow-grey	purple
6-MEL	melibiose	yellow-grey	purple
7-CEL	cellobiose	yellow-grey	purple
8-INO	inositol	yellow-grey	purple
9-XYL	xylose	yellow-grey	purple
10-RAF	raffinose	yellow-grey	purple
11-TRE	trehalose	yellow-grey	purple
12-DUL	dulcitol	yellow-grey	purple

NUMERICAL CODE FORMATION

- 1) The 12 biochemical tests are divided into 4 groups, each containing 3 tests and every test is indicated with a positivity value of 1, 2, or 4.
- value 1 : first positive test of each group (**1-GLU, 4-LAC, 7-CEL, 10-RAF**)
 - value 2 : second positive test of each group (**2-MAL, 5-GAL, 8-INO, 11-TRE**)
 - value 4 : third positive test of each group (**3-SAC, 6-MEL, 9-XYL, 12-DUL**)
 - value 0 : every negative test
- 2) Through the sum in every group of the positive reaction values, a 4 digit code is obtainable, which allows the identification of the microorganism under examination by using the **Table of codes**, as in the following example.

	Group 1			Group 2			Group 3			Group 4		
Test	1-GLU	2-MAL	3-SAC	4-LAC	5-GAL	6-MEL	7-CEL	8-INO	9-XYL	10-RAF	11-TRE	12-DUL
Code of positivity	1	2	4	1	2	4	1	2	4	1	2	4
Results	+	+	+	-	+	-	-	-	+	-	+	-
Sum of codes	7			2			4			2		
CODE	7242			MICROORGANISM			Candida albicans					

TABLE OF REACTIONS

	glucose	maltose	sucrose	lactose	galactose	melibiose	cellobiose	inositol	xylose	raffinose	trehalose	dulcitol
<i>Candida albicans</i>	+	+	+	-	+	-	-	-	+	-	+	-
<i>Candida catenulata</i>	+	+	-	-	+	-	-	-	+	-	-	-
<i>Candida dubliniensis</i>	+	+	+	-	+	-	-	-	-	-	-	-
<i>Candida famata</i>	+	+	+	V	+	+	+	-	+	+	+	V
<i>Candida glabrata</i>	+	-	-	-	-	-	-	-	-	-	+	-
<i>Candida guilliermondii</i>	+	+	+	-	+	+	+	-	+	+	+	+
<i>Candida kefyr</i>	+	-	+	+	+	-	V	-	V	+	V	-
<i>Candida krusei</i>	+	-	-	-	-	-	-	-	-	-	-	-
<i>Candida lambica</i>	+	-	-	-	-	-	-	-	+	-	-	-
<i>Candida lusitanae</i>	+	+	+	-	+	-	+	-	+	-	+	-
<i>Candida parapsilosis</i>	+	+	+	-	+	-	-	-	+	-	-	-
<i>Candida rugosa</i>	+	-	-	-	+	-	-	-	V	-	-	-
<i>Candida tropicalis</i>	+	+	+	-	+	-	V	-	+	-	+	-
<i>Candida zeylanoides</i>	+	-	-	-	V	-	V	-	-	-	+	-
<i>Candida pseudotropicalis</i>	+	-	+	+	+	-	+	-	V	+	-	-
<i>Candida stellatoidea</i>	+	+	-	-	+	-	-	-	+	-	+	-
<i>Cryptococcus neoformans</i>	+	+	+	-	+	-	+	+	+	V	+	+
<i>Cryptococcus albidus</i>	+	+	+	V	V	+	+	+	+	+	+	V
<i>Cryptococcus laurentii</i>	+	+	+	+	+	V	+	+	+	V	+	+
<i>Cryptococcus luteolus</i>	+	+	+	-	+	+	+	+	+	+	+	+
<i>Cryptococcus terreus</i>	+	V	-	V	V	-	+	+	+	-	+	V
<i>Cryptococcus uniguttulatus</i>	+	+	+	-	V	-	V	+	+	V	V	-
<i>Cryptococcus gastricus</i>	+	+	-	V	+	-	+	+	+	-	+	-
<i>Rhodotorula glutinis</i>	+	+	+	-	V	-	+	-	+	+	+	-
<i>Rhodotorula rubra</i>	+	+	+	-	+	-	V	-	+	+	+	-
<i>Saccharomyces cerevisiae</i>	+	+	+	-	+	-	-	-	-	+	V	-
<i>Hansenula anomala</i>	+	+	+	-	+	-	+	-	+	-	+	-
<i>Geotrichum candidum</i>	+	-	-	-	+	-	-	-	+	-	-	-
<i>Blastoschizomyces capitatus</i>	+	-	-	-	+	-	-	-	-	-	-	-
<i>Prototheca wickerhamii</i>	+	-	-	-	+	-	-	-	-	-	+	-
<i>Trichosporon capitatum</i>	+	-	-	-	+	-	-	-	-	-	-	-
<i>Trichosporon pullulans</i>	+	+	+	+	+	+	+	V	V	+	+	-

TABLE OF CODES

CODE	MICROORGANISM	13-WELL COLOR	DIFFERENTIAL CHARACTERISTICS AND MICROSCOPIC MORPHOLOGY ON CORN MEAL AGAR (ref. 11507) AT 25 °C
1000	<i>Candida krusei</i>	PINK	<i>Candida krusei</i> shows extremely long pseudohyphae, rarely ramified; few blastoconids.
1002	<i>Candida zeylanoides</i>	-	<i>Candida zeylanoides</i> shows pseudohyphae, with blastoconids formation. <i>Candida glabrata</i> shows little spores, without pseudohyphae.
	<i>Candida glabrata</i>		
1012	<i>Candida zeylanoides</i>	-	
1040	<i>Candida lambica</i>	-	
1072	<i>Cryptococcus terreus</i>	-	
1076	<i>Cryptococcus terreus</i>	-	
1172	<i>Cryptococcus terreus</i>	-	
1176	<i>Cryptococcus terreus</i>	-	
1200	<i>Trichosporon capitatum</i>	-	<i>Trichosporum capitatum</i> shows hyphae and pseudohyphae, arthroconids and blastoconids. <i>Blastoschizomyces capitatus</i> shows hyphae and ringconids similar to arthroconids. <i>Candida rugosa</i> shows pseudohyphae and blastoconids.
	<i>Blastoschizomyces capitatus</i>		
	<i>Candida rugosa</i>		
1202	<i>Prototheca wickerhamii</i>	-	<i>Prototheca wickerhamii</i> shows spores, without hyphae formation. <i>Candida zeylanoides</i> shows pseudohyphae, with blastoconids formation.
	<i>Candida zeylanoides</i>		
1212	<i>Candida zeylanoides</i>	-	
1240	<i>Geotrichum candidum</i>	-	<i>Geotrichum candidum</i> shows arthroconids. <i>Candida rugosa</i> shows pseudohyphae and blastoconids.
	<i>Candida rugosa</i>		
1272	<i>Cryptococcus terreus</i>	-	
1276	<i>Cryptococcus terreus</i>	-	
1372	<i>Cryptococcus terreus</i>	-	
1376	<i>Cryptococcus terreus</i>	-	
3072	<i>Cryptococcus terreus</i>	-	
3076	<i>Cryptococcus terreus</i>	-	
3172	<i>Cryptococcus terreus</i>	-	
3176	<i>Cryptococcus terreus</i>	-	
3240	<i>Candida catenulata</i>	-	
3242	<i>Candida stellatoidea</i>	-	
3272	<i>Cryptococcus gastricus</i>	-	<i>Cryptococcus gastricus</i> does not grow at 37 °C on Sabouraud Dextrose Agar.
	<i>Cryptococcus terreus</i>	-	<i>Cryptococcus terreus</i> grows at 37 °C on Sabouraud Dextrose Agar.
3276	<i>Cryptococcus terreus</i>	-	
3372	<i>Cryptococcus gastricus</i>	-	<i>Cryptococcus gastricus</i> does not grow at 37 °C on Sabouraud Dextrose Agar.
	<i>Cryptococcus terreus</i>	-	<i>Cryptococcus terreus</i> grows at 37 °C on Sabouraud Dextrose Agar.
3376	<i>Cryptococcus terreus</i>	-	
5303	<i>Candida kefyri</i>	-	
5311	<i>Candida pseudotropicalis</i>	-	<i>Candida pseudotropicalis</i> shows long blastoconids on the pseudohyphae.
	<i>Candida kefyri</i>	-	<i>Candida kefyri</i> shows many long pseudohyphae with long and oval blastoconids.
5313	<i>Candida kefyri</i>	-	
5341	<i>Candida kefyri</i>	-	
5343	<i>Candida kefyri</i>	-	
5351	<i>Candida pseudotropicalis</i>	-	<i>Candida pseudotropicalis</i> shows long blastoconids on the pseudohyphae.
	<i>Candida kefyri</i>	-	<i>Candida kefyri</i> shows many long pseudohyphae with long and oval blastoconids.
5353	<i>Candida kefyri</i>	-	
7053	<i>Rhodotorula glutinis</i>	-	
7060	<i>Cryptococcus uniguttulatus</i>	-	
7061	<i>Cryptococcus uniguttulatus</i>	-	
7062	<i>Cryptococcus uniguttulatus</i>	-	
7063	<i>Cryptococcus uniguttulatus</i>	-	
7070	<i>Cryptococcus uniguttulatus</i>	-	

CODE	MICROORGANISM	13-WELL COLOR	DIFFERENTIAL CHARACTERISTICS AND MICROSCOPIC MORPHOLOGY ON CORN MEAL AGAR (ref. 11507) AT 25 °C
7071	<i>Cryptococcus uniguttulatus</i>	-	
7072	<i>Cryptococcus uniguttulatus</i>	-	
7073	<i>Cryptococcus uniguttulatus</i>	-	
7200	<i>Candida dubliniensis</i>	-	
7201	<i>Saccharomyces cerevisiae</i>	-	
7203	<i>Saccharomyces cerevisiae</i>	-	
7240	<i>Candida parapsilosis</i>	-	
7242	<i>Candida albicans</i>	GREEN	<i>Candida albicans</i> shows pseudohyphae with terminal chlamidospores.
	<i>Candida tropicalis</i>	BLUE	<i>Candida tropicalis</i> shows blastoconids all along the pseudohyphae.
7243	<i>Rhodotorula rubra</i>	-	
7252	<i>Hansenula anomala</i>	-	<i>Hansenula anomala</i> shows blastoconids with ascospores and without pseudohyphae.
	<i>Candida tropicalis</i>	BLUE	<i>Candida tropicalis</i> shows blastoconids all along the pseudohyphae.
	<i>Candida lusitanae</i>	COLORLESS	<i>Candida lusitanae</i> shows short chains of long blastoconids on the curve pseudohyphae.
7253	<i>Rhodotorula rubra</i>	-	<i>Rhodotorula glutinis</i> and <i>Rhodotorula rubra</i> do not form pseudohyphae. <i>Rhodotorula glutinis</i> uses KNO ₃ <i>Rhodotorula rubra</i> does not uses KNO ₃
	<i>Rhodotorula glutinis</i>		
7260	<i>Cryptococcus uniguttulatus</i>	-	
7261	<i>Cryptococcus uniguttulatus</i>	-	
7262	<i>Cryptococcus uniguttulatus</i>	-	
7263	<i>Cryptococcus uniguttulatus</i>	-	
7270	<i>Cryptococcus uniguttulatus</i>	-	
7271	<i>Cryptococcus uniguttulatus</i>	-	
7272	<i>Cryptococcus uniguttulatus</i>	-	
7273	<i>Cryptococcus uniguttulatus</i>	-	
7276	<i>Cryptococcus neoformans</i>	-	
7277	<i>Cryptococcus neoformans</i>	-	
7376	<i>Cryptococcus laurentii</i>	-	
7377	<i>Cryptococcus laurentii</i>	-	
7473	<i>Cryptococcus albidus</i>	-	
7477	<i>Cryptococcus albidus</i>	-	
7573	<i>Cryptococcus albidus</i>	-	
7577	<i>Cryptococcus albidus</i>	-	
7653	<i>Candida famata</i>	-	
7657	<i>Candida guilliermondii</i>	-	<i>Candida guilliermondii</i> forms short pseudohyphae with groups of blastoconids at the septa.
	<i>Candida famata</i>	-	<i>Candida famata</i> does not form pseudohyphae.
7673	<i>Cryptococcus albidus</i>	-	
7677	<i>Cryptococcus luteolus</i>	-	<i>Cryptococcus luteolus</i> uses KNO ₃ .
	<i>Cryptococcus albidus</i>	-	<i>Cryptococcus albidus</i> does not use KNO ₃ .
7713	<i>Trichosporon pullulans</i>	-	
7733	<i>Trichosporon pullulans</i>	-	
7753	<i>Trichosporon pullulans</i>	-	<i>Trichosporon pullulans</i> shows hyphae and pseudohyphae, arthroconids and blastoconids.
	<i>Candida famata</i>	-	<i>Candida famata</i> shows blastoconids, does not form pseudohyphae.
7757	<i>Candida famata</i>	-	
7773	<i>Trichosporon pullulans</i>	-	<i>Trichosporon pullulans</i> shows hyphae and pseudohyphae, arthroconids and blastoconids.
	<i>Cryptococcus albidus</i>	-	<i>Cryptococcus albidus</i> shows dark round spores; does not form hyphae.
7776	<i>Cryptococcus laurentii</i>	-	
7777	<i>Cryptococcus laurentii</i>	-	<i>Cryptococcus laurentii</i> grows at 37 °C on Sabouraud Dextrose Agar.
	<i>Cryptococcus albidus</i>	-	<i>Cryptococcus albidus</i> does not grow at 37 °C on Sabouraud Dextrose Agar.

SENSITIVITY TO ANTIMYCOTICS

Observe the color change of wells and interpret the results using the table no.4.

Write down data in the Test Result Form.

The control well (**24-Growth**) has to turn yellow (positive result).

Table no.4

WELL COLOR	BACTERIAL GROWTH	INTERPRETATION
Red	Inhibited	S = Sensitive
Orange	Partially inhibited	I = Intermediate
Yellow	Good	R = Resistant

QUALITY CONTROL

Every lot of **INTEGRAL SYSTEM YEASTS *Plus*** is subjected to the quality control using the following reference microorganisms:

<i>Candida albicans</i>	ATCC 90028
<i>Candida krusei</i>	ATCC 6258
<i>Candida tropicalis</i>	ATCC 750
<i>Candida parapsilosis</i>	ATCC 22019

FACTORS THAT MAY INVALIDATE THE RESULTS

Poor standardisation of the inoculum; application of the method to microorganisms not belonging to the yeasts group; improper application of the test procedure; non compliance with temperatures and times of incubation.

PRECAUTIONS

The product **INTEGRAL SYSTEM YEASTS *Plus*** is not classified as hazardous under current legislation, however see the safety data sheet for a correct use.

INTEGRAL SYSTEM YEASTS *Plus* is a disposable device to be used only for diagnostic use *in vitro*. It must be used in the laboratory by properly trained personnel, using approved aseptic and safety methods for handling pathogenic agents.

STORAGE

Store at 2-8 °C in the original packaging. Keep away from sources of heat and avoid excessive changes in temperature. In such conditions, the product will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, **INTEGRAL SYSTEM YEASTS *Plus*** and material that has come into contact with the sample must be decontaminated and disposed of in accordance with the techniques used in the laboratory for decontamination and disposal of potentially infected material.

PRESENTATION

Product	Ref.	Packaging
INTEGRAL SYSTEM YEASTS <i>Plus</i>	71822	20 tests
INTEGRAL SYSTEM YEASTS <i>Plus</i>	79822	4 tests

TABLE OF SYMBOLS

IVD <i>In Vitro</i> Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

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BIBLIOGRAFIA / BIBLIOGRAPHY / BIBLIOGRAPHIE

1. Murray, Baron, Pfaller, Tenorev and Tenover. *Manual of Clinical Microbiology 8th Edition* (2003)
2. Bayley and Scott's. *Diagnostic Microbiology 7th Edition* (1986).
3. Davise H. Larone. *"Medically important fungi: a guide to identification"*. 2th Edition. Elsevier (1987).
4. J.P. Van der Walt, L.Rodrigues De Miranda, S.A.Meyer. *"The Yeasts "*. 3th Edition. Elsevier (1984).
5. Barchiesi, F., A.M. Tortorano, I.F. Di Francesco, M. Cogliati, G.scalise, and M.A.Viviani. *In vitro activity of five antifungal agents against uncommon clinical isolates of Candida spp.* J.Antimicrob.Chemother. **43**.295-299 (1999).
6. NCCLS . *Reference method for broth dilution antifungal susceptibility testing of yeasts. M27-A2.*(2002).

Microbiology Products

LIOFILCHEM S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto D.A. (TE) - ITALY

Tel +39 0858930745 Fax +39 0858930330

www.liofilchem.net E-mail: liofilchem@liofilchem.net