TERGITOL 1.5% Supplement

Selective supplement for coliforms isolation.

DESCRIPTION

TERGITOL 1.5% Supplement is a liquid supplement used with Tergitol Agar Base (ref. 610184 / 620184). The complete medium is used for enumeration and identification of *E. coli* and coliform bacteria in water.

KIT CONTENTS

Each kit contains:

- 4 bottles of TERGITOL 1.5% Supplement

An instruction sheet is available from the online Support Center: liofilchem.com/ifu-sds

PRINCIPLE OF THE METHOD

Tergitol NP-7 inhibits the growth of Gram-positive microorganisms, spore forming Gram-negative bacteria and swarming of *Proteus* spp. favouring the recovery of coliform bacteria.

COMPOSITION (g/l)

Tergitol NP-7 15 g	Deionized Water 1000 ml
--------------------	-------------------------

PROCEDURE FOR USE

Aseptically add 7 ml of TERGITOL 1.5% Supplement to 1 litre of Tergitol Agar Base, autoclaved as directed and cooled to 45-50°C. Mix well and pour into sterile Petri dishes.

NOTE: Another supplement, TTC 1% Supplement (ref. 80300 / 80430), is used for the preparation of Tergitol Agar Base and is added aseptically (2.5 ml TTC 1% to 1 litre base medium) after sterilization of the medium.

TECHNIQUE AND RESULTS INTERPRETATION

Refer to the technical sheet of the medium being prepared.

QUALITY CONTROL

1. Appearance control: Clear to slightly opalescent, colourless viscous liquid
2. Microbiological control.

Prepare the medium per label directions. Inoculate the medium with the strains indicated below and incubate at $36 \pm 2^\circ\text{C}$ for 21 ± 3 h.

Control strain		Growth
<i>Escherichia coli</i>	WDCM 00179	Good
<i>Escherichia coli</i>	WDCM 00013	Good
<i>Enterobacter aerogenes</i>	WDCM 00175	Good
<i>Citrobacter freundii</i>	WDCM 00006	Good
<i>Enterococcus faecalis</i>	WDCM 00009	Inhibited

WARNING AND PRECAUTIONS

For professional use only. Operators must be trained and have certain experience in the laboratory methods. Please read the instructions carefully before using this product. Reliability of assay results cannot be guaranteed if there are any deviations from the instructions in this document.

Consult the Safety Data Sheet (SDS) for information regarding hazards and safe handling practices.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration.

PRESENTATION

Product	REF	Σ
TERGITOL 1.5% Supplement	80042	4 bottles x 40 ml

REFERENCES

- Chapman GH. (1951). A culture medium for detecting and confirming *Escherichia coli* in ten hours. *Am. J. Public Health* 41:1381.
- Chapman GH. (1947). A superior culture medium for the enumeration and differentiation of coliforms. *J. Bacteriol.* 53:504.
- Kulp W, Mascoli C, Tavshanjian O (1953). Use of tergitol-7 triphenyl tetrazolium chloride agar as the coliform confirmatory medium in routine sanitary water analysis. *Am.J.Public Health* 43:1111.

TABLE OF SYMBOLS

LOT Batch code	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Consult instructions for use	

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BIBLIOGRAFIA

- Chapman GH. (1951). A culture medium for detecting and confirming *Escherichia coli* in ten hours. Am. J. Public Health 41:1381.
- Chapman GH. (1947). A superior culture medium for the enumeration and differentiation of coliforms. J. Bacteriol. 53:504.
- Kulp W, Mascoli C, Tavshanjian O (1953). Use of tergitol-7 triphenyl tetrazolium chloride agar as the coliform confirmatory medium in routine sanitary water analysis. Am.J.Public Health 43:1111.

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LIOFILCHEM® S.r.l.