MULLER KAUFFMANN Supplement

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
E-Mail: liofilchem@liofilchem.net

MULLER KAUFFMANN *Supplement*

Supplement for the preparation of *Muller Kauffmann* selective enrichment broth for isolation of salmonella.

DESCRIPTION

MULLER KAUFFMANN *Supplement* is a selective supplement to use for the preparation of Muller Kauffmann selective enrichment broth indicated for isolation of salmonella in water, food and clinical samples. It is made of two components: the Iodine Solution and the Brilliant Green Solution. MULLER KAUFFMANN *Supplement* is used as selective supplement of MULLER KAUFFMANN BROTH cod. 610035 or 620035.

KIT CONTENTS

Each kit contains:

- 2 bottles of Iodine Solution 50 ml.
- 1 bottle of Brilliant Green Solution 0.1% 50 ml.
- 1 instruction sheet

PRINCIPLE OF THE METHOD

MULLER KAUFFMANN *Supplement* is a selective supplement made of two components:

- Iodine Solution whose components react with Sodium Thiosulphate contained in medium MULLER KAUFFMANN BROTH cod. 610035 or 620035 forming the compound Sodium Tetrathionate that inhibits most of coliforms and intestinal bacteria except most of salmonella.
- Brilliant Green is a selective agent that inhibits many Gram-positive and Gram-negative microorganisms except most of salmonella.

COMPOSITION

MULLER KAUFFMANN <i>Supplement</i>			
Iodine Solution		Brilliant Green 0.1%	
<i>Contents / bottle</i>		<i>Contents / bottle</i>	
Iodine	10.0 g	Brilliant Green	50.0 mg
Potassium Iodide	12.5 g	Distilled water	50.0 ml
Distilled water	50.0 ml		

PROCEDURE FOR USE

1. Aseptically take 19 ml of Iodine Solution and 9.5 ml of Brilliant Green 0.1% and add them to 1 litre of medium MULLER KAUFFMANN BROTH cod. 610035 or 620035, boiled and cooled to 45-50 °C.
2. Mix with care.
3. Distribute into the final containers.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation for MULLER KAUFFMANN BROTH cod. 610035 or 620035.

QUALITY CONTROL

- Visual inspection:

the Iodine Solution is clear, without precipitate, the colour is yellow dark-brown;
the Brilliant Green Solution 0.1% is light green colour, with white precipitate.

- Microbiological control:

prepare the selective enrichment medium using as base MULLER KAUFFMANN BROTH cod. 610035 or 620035 added with MULLER KAUFFMANN *Supplement* in the quantity listed upon. The medium is inoculated with the strains indicated in the microbiological control table and it is incubated at 42-43 °C for one night. The broth culture obtained is afterwards inoculated on non selective solid medium and incubated at 36 ± 1 °C for 18-24 hours.

Microbiological control

Control strains	Growth on non selective solid medium	
<i>Escherichia coli</i>	ATCC 25922	Partial or total inhibition
<i>Proteus mirabilis</i>	ATCC 25933	Good or partial inhibition
<i>Salmonella typhimurium</i>	ATCC 14028	Good

PRECAUTIONS

The product MULLER KAUFFMANN *Supplement* is not classified as dangerous under current legislation; it is nevertheless recommended that the Safety Data Sheet be consulted on its correct use.

MULLER KAUFFMANN *Supplement* is a supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store MULLER KAUFFMANN *Supplement* at 2-8 °C in its original packaging. In such conditions MULLER KAUFFMANN *Supplement* will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

REFERENCES

- Kauffmann, F. (1935). Z. F. Hyg. **117**: 26-32.
- Kauffmann, F. (1930). Z. F. Hyg. **113**: 148-157.
- Muller, L. (1923) C.R. Soc. Biol. (Paris) **89**: 434-443.
- Jeffries, L. (1959). J. Clin. Path. **12**: 568-570.

PRESENTATION

product	REF	Σ
MULLER KAUFFMANN <i>Supplement</i>	80047	3 bottles x 50 ml

TABLE OF SYMBOLS

IVD In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

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