Legionella Growth Supplement

Supplement for presumptive identification of *Legionella* spp.

DESCRIPTION

Legionella Growth Supplement is a lyophilized supplement used for the preparation of Legionella BCYE Agar Base (ref. 610049, 620049). The complete medium is used for detection and enumeration of *Legionella* spp from clinical specimens and environmental samples, according to ISO 11731 (all parts).

KIT CONTENTS

Each kit contains:

- 10 vials of freeze-dried Legionella Growth Supplement;
- 1 instruction sheet.

PRINCIPLE OF THE METHOD

Legionella Growth Supplement is used for confirmation of presumed *Legionella* colonies by a comparative examination demonstrating the organism requirements for L-cysteine.

COMPOSITION

| | Content / vial | Content / liter of medium |
|--------------------|----------------|---------------------------|
| L-Cysteine | 0.04 g | 0.4 g |
| Iron Pyrophosphate | 0.025 g | 0.25 g |

PROCEDURE FOR USE

1. Aseptically reconstitute the content of one vial of Growth Supplement with 10 ml sterile distilled water.
2. Mix to complete dissolution and add to 90 ml Legionella BCYE Agar Base, autoclaved and cooled to 45-50°C.
3. Mix well and dispense in Petri dishes.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation of the medium being prepared.

QUALITY CONTROL

1. Control of the appearance: freeze-dried white in colour.
2. Microbiological control.
Prepare the medium per label directions. Inoculate the plates with the strains indicated below and incubate at 36 ± 2°C for 2-5 days in humidified atmosphere.

| Control strains | | Growth |
|-------------------------------|-------------|--------|
| <i>Legionella pneumophila</i> | WDCM 00107 | Good |
| <i>Escherichia coli</i> | ATCC® 25922 | Good |

PRECAUTIONS

Legionella Growth Supplement is classifiable as hazardous under current legislation; it is recommended that the Safety Data Sheet be consulted on its use. The product is a selective supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration. Once reconstituted, the product can be stored for a maximum duration of 30 days at -20°C, shielded from light.

REFERENCES

- ISO 11731-2:2004. Water quality – Detection and enumeration of Legionella. Part 2: Direct membrane filtration method for waters with low bacterial count.
- ISO 11731:1998. Water quality – Detection and enumeration of Legionella.
- Edelstein P.H. (1981) Improved semiselective medium for the isolation of *Legionella pneumoniae* from contaminated clinical and environmental specimens. J. Clin. Microbiol. 14(3):298.

PRESENTATION

| Product | Ref. | Contents |
|------------------------------|-------|----------|
| Legionella Growth Supplement | 80056 | 10 vials |

One vial is sufficient to prepare 100 ml of medium.

TABLE OF SYMBOLS

| | | | | |
|---|--|--|--|--|
|  Batch code |  In Vitro Diagnostic Medical Device |  Manufacturer |  Contains sufficient for <n> tests |  Temperature limitation |
|  Catalogue number |  Fragile, handle with care |  Use by |  Caution, consult accompanying documents |  Do not reuse |

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LIOFILCHEM® S.r.l.