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F02813
Rev.3.1 / 16.02.2016

Kit for pretreating and decontaminating the samples that are to be analyzed for mycobacteria research.

DESCRIPTION

DECONTAM-KIT is a kit for pretreating and decontaminating the samples that are to be analyzed for mycobacteria research. The configuration is available in the 90-test version.

CONTENT OF THE PACKAGES

Each package contains 9 bottles of Reagent A, 9 100-ml bottles of Reagent B, 9 bottles of Reagent C and 1 instructions sheet.

PRINCIPLE OF THE METHOD

The material to be analyzed often contains associated bacterial flora that can cover the culture medium in such a way as to hide or inhibit the growth of any mycobacteria that may be present. For most samples, it is also useful to fluidify the material to be analysed. Pre-treatment must in fact make the sample fluid and homogenous and enable elimination of the associated microbial flora to obtain a better concentration of any mycobacteria that may be present.

Adding Tween 80 before centrifugation encourages the subsequent separation by centrifugation of the mycobacteria, whose density will remain below 1.0.

The method that uses N-acetyl-L-cysteine/NaOH is based on:

1. the decontaminating capacities of sodium hydroxide, which acts as a mucolytic agent;
2. the centrifugation process that facilitates liquefaction;
3. the capacity of N-Acetyl-L-cysteine to fluidify sputum and to enable treatment with low concentrations of NaOH, thereby increasing the capacity to recover the mycobacteria.

COMPOSITION

Reagent A	Sodium citrate	1.45 g	N-acetyl-L-cysteine	0.5 g
Reagent B	Sodium hydroxide	2.0 g	Deionized water	100 ml
Reagent C	Monobasic potassium phosphate	1.452 g	Dibasic sodium phosphate	1.659 g

COLLECTION OF SAMPLES

1. The samples in which bacterial flora is not normally present (such as cerebrospinal fluid, blood or other matter obtained in sterile conditions) do not need to be decontaminated and can be sown directly onto the culture medium after centrifugation.
2. The liquid samples with a volume greater than 10 ml must be pre-concentrated before any decontamination.
3. Wash the buffers with a sterile sodium chloride solution. The washing water will be treated as a material to be analysed.
4. The textiles must be homogenised by grinding with sterile sand in a mortar or by using a homogenizer.
5. Liquefy the faeces with a sterile sodium chloride solution in the proportion 1:5 and decontaminate for no more than 25 minutes.

PREPARATION PROCEDURE

N-acetyl-L-cysteine/NaOH Solution

Pour the contents of a bottle of Reagent A into the bottle of Reagent B and shake until contents have dissolved completely. The solution that is thus made up is unstable and should be used within a day. The solution must be stored at 2-8°C and the quantity is sufficient for treating 10 samples.

Phosphate Buffer Solution

In sterile conditions, add 5 ml of sterile distilled water to a bottle of Reagent C 5. Shake until complete dissolution and add aseptically to 295 ml of sterile distilled water. Shake well. The buffer solution thereby obtained has a pH value of 6.8 (check the pH of the solution and if necessary adjust it to 6.8 ± 0.1) and remains stable for several weeks if stored at 2-8°C. Once it is prepared, sterilise the phosphate buffer at 121 °C for 20 minutes or filter it with a 0.2 micron filter before use.

TEST PROCEDURE

1. Insert up to 10 ml of the material to be analysed into a disposable centrifugal test tube with capacity from 30 to 50 ml.
2. Add an equal volume of the previously made up N-acetyl-L-cysteine/NaOH solution, close the test tube tightly and mix for 20 minutes in a laboratory stirrer.
3. Then fill the test tube up to 1 cm below the edge with the previously prepared Phosphate Buffer Solution. 1 to 2 drops of 5% Tween 80 Solution can be added.
4. Centrifuge at a minimum speed of 3000 rpm for 20 minutes.
5. Eliminate the supernatant and resuspend the sediment with 1 ml of sterile sodium chloride solution. If desired, 1 ml of Phosphate Buffer Solution can also be used.
6. Transfer the suspension to culture media.

PRECAUTIONS

The DECONTAM-KIT product contains 2% sodium hydroxide in the Reagent B configuration. At this concentration, the product is classed as hazardous by current legislation: before using it, consult the safety data sheet.

DECONTAM-KIT is a disposable device only for diagnostic *in vitro* professional use. It must be used in the laboratory by properly trained operators using approved aseptic and safety methods for pathogenic agents.

STORAGE

Store DECONTAM-KIT at 10-25°C in their original packaging. Do not store near sources of heat and do not expose to excessive temperature variations. In such conditions, DECONTAM-KIT can be used until the expiry date shown on the label. Do not use after this date. Dispose of if they show signs of deterioration.

DISPOSAL OF USED MATERIAL

After use, DECONTAM-KIT and the material that comes into contact with the sample must be decontaminated and disposed of in accordance with current laboratory techniques for the decontamination and disposal of potentially infected material.

BIBLIOGRAPHY

1. DIN 58934-3, 1996-12. Tuberkulosediagnostik, Teil 3: Kulturelle Methoden zum Nachweis von Mykobakterien.
2. Kubica G.P., DYE W.E., Cohn M.L. ed al. Sputum digestion and decontamination with N-acetyl-L-cysteine-sodium hydroxide for culture of mycobacteria. Amer.Rev.Resp.Dis.,1963, 87, 775.
3. Kubica G.P., Kaufmann A.J., DYE W.E. Comments on the use of the new mucolytic agent N-acetyl-L-cysteine, as a sputum digestant for the isolation of mycobacteria. Amer.Rev.Resp.Dis. 1964, 89, 284-288

PRESENTATION

Product	Ref.	Test
DECONTAM-KIT	80060	90 tests

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains Sufficient for <n> tests	 Caution, consult accompanying documents	 Do not reuse

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LIOFILCHEM® S.r.l.