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CLOSTRIDIUM difficile Supplement

ENGLISH

Selective supplement for the enrichment of the medium CLOSTRIDIUM difficile AGAR BASE for the isolation of *Clostridium difficile*

DESCRIPTION

CLOSTRIDIUM difficile Supplement is a selective supplement for the isolation of *Clostridium difficile*, made of a una freeze-dried D-Cicloserin and Cefoxitin mixture. CLOSTRIDIUM difficile Supplement is used for the selective enrichment of CLOSTRIDIUM difficile AGAR BASE medium code 610115 or 620115.

PACKAGE CONTENTS

Each package contains:

- 10 bottles of freeze-dried CLOSTRIDIUM difficile Supplement
- 1 instruction sheet

PRINCIPLE OF THE METHOD

The selectivity of CLOSTRIDIUM difficile Supplement, set by Levett, is due to the activity of D-Cicloserin and Cefoxitin which inhibit the growth of most of Enterobacteriaceae, such as *Enterococcus faecalis*, stafilococci, Gram-negative anaerobic non sporogens bacilli and some strains of Clostridia except for *C. difficile*.

COMPOSITION

CLOSTRIDIUM difficile Supplement		
	Content / bottle	Content / l of medium
Cicloserin	125.0 mg	250.0 mg
Cefoxitin	4.0 mg	8.0 mg

TEST PROCEDURE

1. Reconstitute aseptically the content of one bottle of CLOSTRIDIUM difficile Supplement with 2 ml of sterile distilled water. Shake until completely dissolved, avoiding foam formation.
2. Add aseptically the entire content of one bottle (2 ml) to 500 ml of CLOSTRIDIUM difficile Agar Base medium code 610115-620115 autoclaved, cooled at 45-50 °C and with the addition of 7% of defibrinated horse blood (in this case preferred to ram's blood).
3. Mix with care.
4. Distribute into Petri dishes.

TECHNIQUE AND INTERPRETATION OF RESULTS

Refer to the technical documentation for CLOSTRIDIUM difficile AGAR BASE code 610115 or 620115.

QUALITY CONTROL

1. Control of the appearance: freeze-dried product, white colour.
2. Microbiological control .

Prepare the plates using as base CLOSTRIDIUM difficile AGAR BASE code 610115 or 620115 supplemented with CLOSTRIDIUM difficile Supplement (1 bottle in 500 ml of medium) and with 7% of defibrinated horse blood. Plates are inoculated with the strains indicated in the microbiological control table.

Incubation conditions: 24 h at 36 ± 1 °C, in anaerobiosis.

Microbiological control:

	Control strains	Growth
<i>Clostridium difficile</i>	ATCC 11204	Good
<i>Escherichia coli</i>	ATCC 25922	Inhibited

PRECAUTIONS

The product CLOSTRIDIUM difficile Supplement is classifiable as hazardous under current legislation; it is recommended that the Safety Data Sheet be consulted on its use.

CLOSTRIDIUM difficile Supplement is a selective supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store CLOSTRIDIUM difficile Supplement at 2-8 °C in its original packaging. In such conditions CLOSTRIDIUM difficile Supplement will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

REFERENCES

- Levett (1985) J. Clin. Pathol. 38. 233-234.
- Hall, I. and O'Toole, E. (1935) Am. J. Dis. Child. 49. 390.

PRESENTATION

Product	REF	Σ
CLOSTRIDIUM difficile Supplement	81007	10 bottles

One bottle is sufficient to prepare 500 ml of medium.

TABLE OF SYMBOLS

IVD In Vitro Diagnostic Medical Device	Do not reuse	Manufacturer	Contains sufficient for <n> tests	Temperature limitation
REF Catalogue number	Fragile, handle with care	Use by	Caution, consult accompanying documents	LOT Batch code

LIOFILCHEM Bacteriology Products
Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

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