# HAEMOPHILUS Supplement

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Rev.1 / 21.01.2008

# HAEMOPHILUS Supplement

Selective supplement for the preparation of G.C. Medium for *Haemophilus spp.* cultivation.

## DESCRIPTION

HAEMOPHILUS Supplement is a freeze-dried selective supplement constituted by Bacitracin, Vancomycin and Clindamycin to use as a supplement of the culture medium G.C. Medium cod. 610022 or 620022 for *Haemophilus spp.* cultivation.

## KIT CONTENTS

Each kit contains:

- 10 vials of HAEMOPHILUS Supplement freeze-dried.
- 1 instructions sheet

## PRINCIPLE OF THE METHOD

HAEMOPHILUS Supplement supports *Haemophilus spp.* growth thanks to the selective activity of the antibiotics: Bacitracin, Vancomycin and Clindamycin inhibit Gram-positive bacteria growth interfering with cell wall synthesis.

## COMPOSITION

HAEMOPHILUS Supplement		
	Content / vial	Content / l of medium
Bacitracin	150.0 mg	300.0 mg
Vancomycin	5.0 mg	10.0 mg
Clindamycin	0.5 mg	1.0 mg

## PROCEDURE FOR USE

1. Aseptically reconstitute the content of one vial of HAEMOPHILUS Supplement with 5 mL of sterile distilled water .
2. Mix to complete dissolution and add to 500 mL of G.C. Medium cod. 610022 or 620022, autoclaved, cooled to 45-50°C and added with 5% of defibrinated horse blood cooked as chocolate and with 1 vial of Vitalex growth Supplement reconstituted with 5 mL of Vitalex restoring solution (cod. 81023).
3. Mix well and dispense in Petri dishes.

## TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation for G.C. Medium cod. 610022 or 620022.

## QUALITY CONTROL

1. Control of the appearance: freeze-dried pale yellow in colour.
2. Microbiological control.

Prepare the plates using as base G.C. Medium cod. 610022 or 620022 added with the content of 1 vial of HAEMOPHILUS Supplement (1 vial in 500 mL of medium). The plates are inoculated with the strains indicated in the microbiological control table. Incubation conditions: 36 ± 1 °C for 48-72 hours in atmosphere containing 5-10% CO<sub>2</sub>. Microbiological control:

	Control strains	Growth
<i>Haemophilus influenzae</i>	ATCC 10211	Good
<i>Staphylococcus aureus</i>	ATCC25923	Inhibited

## PRECAUTIONS

HAEMOPHILUS Supplement contains dangerous substances according to directives 1999/45/CE and 2001/60/CE. HAEMOPHILUS Supplement is a selective supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

## STORAGE

Store HAEMOPHILUS Supplement at 2-8 °C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration. Once reconstituted , the product can be stored for a maximum duration of 30 days at -20°C, shielded from light.

## REFERENCES

- Chapin C.K., G.V. Doern (1983). J. Clin. Microbiol. 17: 1163-1165.
- Martin, J.E., Armstrong J.H., Smith, P.B. (1974). Appl. Microbiol. 27: 802-805.
- NCCLS document M22-A2, 1996. Quality Assurance for Commercially Prepared Microbiological Culture Media- Second Ed. Approved Standard.

## PRESENTATION

Product	REF	Σ
HAEMOPHILUS Supplement	81014	10 vials

One vial is sufficient to prepare 500 mL of medium

## TABLE OF SYMBOLS

 In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	 Batch code

**LIOFILCHEM Bacteriology Products**

Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

Tel. +390858930745 Fax +390858930330 Website: [www.liofilchem.net](http://www.liofilchem.net) E-Mail: [liofilchem@liofilchem.net](mailto:liofilchem@liofilchem.net)

Rev.1 / 21.01.2008