VITALEX growth Supplement

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Rev.2 /03.07.2013

VITALEX growth Supplement

Growth supplement for the isolation of fastidious microorganisms.

DESCRIPTION

VITALEX growth Supplement is a growth supplement for the isolation of fastidious microorganisms and it is made of a freeze-dried mixture of aminoacids, vitamins and other growth substrata. VITALEX growth Supplement is used to enrich the G.C. Medium ref. 610022 or 620022.

CONTENT OF THE PACKAGES

Each package contains:

- 10 bottles of freeze-dried VITALEX growth Supplement + 10 bottles of VITALEX growth diluent.
- 1 instruction sheet.

PRINCIPLE OF THE METHOD

VITALEX growth Supplement provides vitamins, amino acids, coenzymes, glucose, ferric ions and other factors that specifically promote neisseria growth.

COMPOSITION

	VITALEX growth Supplement		VITALEX growth diluent	
	Content / bottle	Content / L of medium		Content / bottle
Glutamine	100.00 mg	200.00 mg	Glucose	0.5 g
Adenine	10.00 mg	20.00 mg	Distilled water	5.0 ml
Guanine	0.30 mg	0.60 mg		
Aminobenzoic acid	0.13 mg	0.26 mg		
NAD	2.50 mg	5.00 mg		
Coccarboxylase	1.00 mg	2.00 mg		
Ferric nitrate	0.20 mg	0.40 mg		
Thiamine	0.03 mg	0.06 mg		
Vitamin B12	0.10 mg	0.20 mg		

TEST PROCEDURE

1. Aseptically reconstitute the content of one bottle of VITALEX growth Supplement with 5 mL of VITALEX growth diluent (10% glucose solution). Shake until completely dissolved, avoiding foam formation.
2. Aseptically add the entire content of one vial (5 mL) to 500 mL of G.C. Medium ref. 610022 or 620022, autoclaved and cooled at 45-50 °C.
3. Mix with care.
4. Distribute into Petri dishes.

TECHNIQUE AND INTERPRETATION OF RESULTS

Refer to the technical sheet of the culture medium G.C. Medium ref. 610022 or 620022 to be enriched with VITALEX growth Supplement.

QUALITY CONTROL FOR THE USER

Appearance: pink colour freeze-dried medium.

Microbiological control:

Prepare plates of G.C. Medium ref. 610022 enriched with VITALEX growth Supplement (2 bottles per 1000 mL of medium). The plates are seeded with the strains indicated in the microbiological control table.

Incubation conditions: 24-48 h at 36 ± 1 °C.

Microbiological control.

Control strains		Growth
<i>Escherichia coli</i>	ATCC® 25922	good
<i>Haemophilus parainfluenzae</i>	ATCC® 7901	good
<i>Neisseria meningitidis</i>	ATCC® 13090	good

PRECAUTIONS

The product VITALEX growth Supplement is not classified as dangerous under current legislation; it is nevertheless recommended that the Safety Data Sheet be consulted on its correct use. VITALEX growth Supplement is a supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store VITALEX growth Supplement at 2-8 °C in its original package. In such conditions the product is valid until the expiry date shown on the label. Do not use them beyond that date. Dispose of them if they show signs of deterioration.

DISPOSAL OF USED MATERIAL

After the use, VITALEX growth Supplement and the material that has come into contact with the sample must be decontaminated and disposed of in accordance with the laboratory procedures for the decontamination and disposal of potentially infected material.

REFERENCES

1. Piccolomini, R., Di Bonaventura, G., Festi, D., et al. (1997) J. Clin. Microbiol. 35: 1541-1544.

PRESENTATION

Product	REF	Σ
VITALEX growth Supplement	81023	10 bottles + 10 bottles

One bottle is sufficient to prepare 500 mL of medium.

TABLE OF SYMBOLS

IVD In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

LIOFILCHEM® S.r.l.

Via Scozia, Zona Ind.le - 64026, Roseto degli Abruzzi (TE) - ITALY
Tel +39 0858930745 Fax +39 0858930330 Website: www.liofilchem.net E-mail: liofilchem@liofilchem.net

Rev.2 /03.07.2013