V.C.N.T. Supplement

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Rev.2 / 26.04.2012

V.C.N.T. Supplement

Selective supplement for the preparation of G.C. Medium for *Neisseria* spp. cultivation

DESCRIPTION

V.C.N.T. *Supplement* is a freeze-dried selective supplement constituted by Vancomycin, Colistin, Nystatin and Trimethoprim to use as a supplement of the culture medium G.C. Medium ref. 610022 or 620022 for *Neisseria* spp. cultivation.

KIT CONTENTS

Each kit contains:

- 10 vials of V.C.N.T. *Supplement* freeze-dried.
- 1 instructions sheet

PRINCIPLE OF THE METHOD

V.C.N.T. *Supplement* supports *Neisseria* spp growth thanks to the selective activity of the antibiotics: Vancomycin is a glycopeptide which inhibits Gram-positive bacteria growth interfering with cell wall synthesis, while Colistin inhibits Gram-negative bacteria growth; Nystatin inhibits fungi and yeasts interfering with fungal membrane permeability by binding to sterols. The antibacterial power of Trimethoprim is due to its inhibitory effect on Gram-positive and Gram-negative dihydrofolate reductase.

COMPOSITION

V.C.N.T. Supplement		
	Content / vial	Content / l of medium
Vancomycin	1.5 mg	2.0 mg
Colistin	3.75 mg	7.5 mg
Nystatin	6250 IU	12500 IU
Trimethoprim	2.5 mg	5.0 mg

PROCEDURE FOR USE

1. Aseptically reconstitute the content of one vial of V.C.N.T. *Supplement* with 5 ml of sterile distilled water.
2. Mix to complete dissolution and add to 500 ml of G.C. Medium ref. 610022 or 620022, autoclaved, cooled to 45-50°C and added with 5% of defibrinated horse blood cooked as chocolate and with 1 vial of Vitalex Growth *Supplement* reconstituted with 5 ml of Vitalex Restoring Solution ref. 81023.
3. Mix well and dispense in Petri dishes.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation for G.C. Medium ref. 610022 or 620022.

QUALITY CONTROL

1. Control of the appearance: freeze-dried pale yellow in colour.
2. Microbiological control.

Prepare the plates using as base G.C. Medium ref. 610022 or 620022 added with the content of 1 vial of V.C.N.T. *Supplement* (1 vial in 500 ml of medium). The plates are inoculated with the strains indicated in the microbiological control table. Incubation conditions: 36 ± 1°C for 48-72 hours in atmosphere containing 5-10% CO₂. Microbiological control:

	Control strains	Growth
<i>Neisseria gonorrhoeae</i>	ATCC 43069	Good
<i>Neisseria meningitidis</i>	ATCC 13090	Good
<i>Escherichia coli</i>	ATCC 25922	Inhibited
<i>Staphylococcus aureus</i>	ATCC 25923	Inhibited

PRECAUTIONS

V.C.N.T. *Supplement* contains dangerous substances according to directives 1999/45/CE and 2001/60/CE. V.C.N.T. *Supplement* is a selective supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store V.C.N.T. *Supplement* at 2-8°C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration. Once reconstituted, the product can be stored for a maximum duration of 30 days at -20°C, shielded from light.

REFERENCES

- Chapin C.K., G.V. Doern (1983). J. Clin. Microbiol. 17: 1163-1165.
- Martin, J.E., Armstrong J.H., Smith, P.B. (1974). Appl. Microbiol. 27: 802-805.
- NCCLS document M22-A2, 1996. Quality Assurance for Commercially Prepared Microbiological Culture Media- Second Ed. Approved Standard.

PRESENTATION

Product	REF	Σ
V.C.N.T. <i>Supplement</i>	81024	10 vials

One vial is sufficient to prepare 500 ml of medium.

TABLE OF SYMBOLS

IVD In Vitro Diagnostic Medical Device	 Do not reuse	Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

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LIOFILCHEM® S.r.l.