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Selective enrichment supplement for the isolation of enterococci when added to the medium KANAMYCIN AESCULIN AZIDE AGAR BASE.

DESCRIPTION

KANAMYCIN *Supplement* is a selective enrichment supplement for the isolation of enterococci from food. KANAMYCIN *Supplement* is used for selective enrichment of the medium KANAMYCIN AESCULIN AZIDE AGAR BASE, ref. 610121 or 620121.

KIT CONTENTS

Each kit contains:

- 10 bottles of lyophilised KANAMYCIN *Supplement*
- 1 instruction sheet.

PRINCIPLE OF THE METHOD

Kanamycin is an antibiotic that belongs to the family of aminoglycosides effective against aerobic gram-negative bacilli and some gram-positive bacteria.

COMPOSITION

KANAMYCIN <i>Supplement</i>		
	Content / bottle	Content / liter of medium
Kanamycin	10 mg	20 mg

PROCEDURE FOR USE

1. Aseptically reconstitute the content of a bottle of KANAMYCIN *Supplement* with 5 ml of sterile distilled water. Shake until completely dissolved avoiding foam formation.
2. Aseptically add the whole contents of one bottle (5 ml) to 500 ml of the medium KANAMYCIN AESCULIN AZIDE AGAR BASE, ref. 610121 or 620121, autoclaved and cooled to 45-50°C.
3. Mix with care.
4. Distribute into Petri dishes.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation of the medium being prepared.

QUALITY CONTROL

1. Visual inspection: clear, colourless solution.
2. Microbiological control.

Prepare the plates using KANAMYCIN AESCULIN AZIDE AGAR BASE, ref. 610121 or 620121, supplemented with KANAMYCIN *Supplement* (1 bottle per 500 ml of medium).

The plates are seeded with the strains listed below. Incubation conditions: 18-24 hours at 36±1°C, in aerobiosis.

Microorganism		Growth	Features
<i>Enterococcus faecium</i>	ATCC® 19434	good	white-gray colonies with black halo
<i>Enterococcus bovis</i>	ATCC® 27960	good	white-gray colonies with black halo
<i>Escherichia coli</i>	ATCC® 25922	inhibited	---
<i>Bacillus subtilis</i>	ATCC® 6633	inhibited	---

PRECAUTIONS

The product KANAMYCIN *Supplement* is classified as dangerous under current legislation. For its use is recommended to consult the safety data sheet. KANAMYCIN *Supplement* is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store KANAMYCIN *Supplement* at 2-8°C in its original packaging. In such conditions KANAMYCIN *Supplement* will remain valid until the expiry date indicated on the label. Do not use beyond that date. Eliminate without using if there are signs of deterioration.

REFERENCES

- Mossel, D.A.A., P.G.H. Bijker, and I. Eelederink I. (1978) Arch. Lebensmittel Hg. 29, 121-127.
- Mossel, D.A.A., I. Eelederink, H. de Vor, and E.D. Keizer. (1976) Lab. Practice 25: 393-395.
- Reuter, G. (1985). Inter. J. Food Microbiol. 2: 103-114.

PRESENTATION

Product	REF	
KANAMYCIN <i>Supplement</i>	81030	10 bottles

TABLE OF SYMBOLS

 Batch code	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	