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Rev.1 / 21.01.2008

GARDNERELLA VAGINALIS Supplement

Selective supplement for the preparation of COLUMBIA AGAR BASE for *Gardnerella vaginalis* isolation.

DESCRIPTION

GARDNERELLA VAGINALIS Supplement is a freeze-dried selective supplement constituted by Gentamycin, Nalidixic Acid and Amphotericin B, to use as a supplement of the culture medium COLUMBIA AGAR BASE cod. 610013 or 620013 for *Gardnerella vaginalis* isolation.

KIT CONTENTS

Each kit contains:

- 10 vials of GARDNERELLA VAGINALIS Supplement freeze-dried.
- 1 instructions sheet

PRINCIPLE OF THE METHOD

GARDNERELLA VAGINALIS Supplement supports *Gardnerella vaginalis* growth thanks to the selective activity of antibiotics: Gentamycin is an antibiotic of wide spectrum which inhibit Gram-positive and Gram-negative bacteria growth; Nalidixic Acid is mainly active against Gram-negative bacteria; Amphotericin B inhibits fungi and yeasts interfering with fungal membrane permeability.

COMPOSITION

GARDNERELLA VAGINALIS Supplement		
	Content / vial	Content / l of medium
Gentamycin	2.0 mg	4.0 mg
Nalidixic Acid	15.0 mg	30.0 mg
Amphotericin B	1.0 mg	2.0 mg

PROCEDURE FOR USE

1. Aseptically reconstitute the content of one vial of GARDNERELLA VAGINALIS Supplement with 5 mL of sterile distilled water/ ethyl alcohol mixture (relation 1: 1).
2. Mix to complete dissolution and add to 500 mL of COLUMBIA AGAR BASE cod. 610013 or 620013, autoclaved, cooled to 45-50°C and added with 10% of defibrinated sheep, horse or rabbit blood .
3. Mix well and dispense in Petri dishes.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation for COLUMBIA AGAR BASE cod. 610013 or 620013.

QUALITY CONTROL

1. Control of the appearance: freeze-dried yellow in colour.
2. Microbiological control.

Prepare the plates using as base COLUMBIA AGAR BASE cod. 610013 or 620013 added with the content of 1 vial of GARDNERELLA VAGINALIS Supplement (1 vial in 500 mL of medium). The plates are inoculated with the strains indicated in the microbiological control table. Incubation conditions: 36 ± 1 °C for 18-48 hours in atmosphere containing 5-10% CO₂. Microbiological control:

Control strains	Growth
<i>Enterococcus faecalis</i>	Poor
<i>Escherichia coli</i>	Inhibited
<i>Yersinia enterocolitica</i>	Good

PRECAUTIONS

GARDNERELLA VAGINALIS Supplement contains dangerous substances according to directives 1999/45/CE and 2001/60/CE. GARDNERELLA VAGINALIS Supplement is a selective supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store GARDNERELLA VAGINALIS Supplement at 2-8 °C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration. Once reconstituted , the product can be stored for a maximum duration of 30 days at -20°C, shielded from light.

REFERENCES

- Ellner, P.D., C.J. Stoessel., E. Drakeford, and F. Vasi (1966). A new culture medium for medical bacteriology. Am. J.Clin. Path. **45**, 502-504.
- Isenberg, H.D. (ed.) (1992). Clinical microbiology procedures handbook, vol. 1 American Society for Microbiology, Washington, DC.

PRESENTATION

Product	REF	
GARDNERELLA VAGINALIS Supplement	81040	10 vials

One vial is sufficient to prepare 500 mL of medium

TABLE OF SYMBOLS

 In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	 Batch code

LIOFILCHEM Bacteriology Products

Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

Tel. +390858930745 Fax +390858930330 Website: www.liofilchem.net

Rev.1 / 21.01.2008

E-Mail: liofilchem@liofilchem.net