V.C.A.T. Supplement

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Rev.1 / 21.01.2008

V.C.A.T. Supplement

Selective supplement for the preparation of G.C. Medium for *Neisseria spp.* cultivation.

DESCRIPTION

V.C.A.T. *Supplement* is a freeze-dried selective supplement constituted by Vancomycin, Colistin, Amphotericin B and Trimethoprim to use as a supplement of the culture medium G.C. Medium cod. 610022 or 620022 for *Neisseria spp.* cultivation.

KIT CONTENTS

Each kit contains:

- 10 vials of V.C.A.T. *Supplement* freeze-dried.
- 1 instructions sheet

PRINCIPLE OF THE METHOD

V.C.A.T. *Supplement* supports *Neisseria spp* growth thanks to the selective activity of the antibiotics: Vancomycin is a glycopeptide which inhibits Gram-positive bacteria growth interfering with cell wall synthesis, while Colistin inhibits Gram-negative bacteria growth; Amphotericin B inhibits fungi and yeasts interfering with fungal membrane permeability by forming channels in the membrane and causing small molecules to leak out. The antibacterial power of Trimethoprim is due to its inhibitory effect on Gram-positive and Gram-negative dihydrofolate reductase.

COMPOSITION

V.C.A.T. Supplement		
	Content / vial	Content / l of medium
Vancomycin	1.0 mg	2.0 mg
Colistin	3.75 mg	7.5 mg
Amphotericin B	0.5 mg	1.0 mg
Trimethoprim	1.5 mg	3.0 mg

PROCEDURE FOR USE

1. Aseptically reconstitute the content of one vial of V.C.A.T. *Supplement* with 4 mL of ethyl alcohol/ sterile distilled water mixture (relation 1: 1).
2. Mix to complete dissolution and add to 500 mL of G.C. Medium cod. 610022 or 620022, autoclaved, cooled to 45-50°C and added with 5% of defibrinated horse blood cooked as chocolate and with 1 vial of Vitalex growth *Supplement* reconstituted with 5 mL of Vitalex restoring solution (cod. 81023).
3. Mix well and dispense in Petri dishes.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation for G.C. Medium cod. 610022 or 620022.

QUALITY CONTROL

1. Control of the appearance: freeze-dried yellow in colour.
2. Microbiological control.

Prepare the plates using as base G.C. Medium cod. 610022 or 620022 added with the content of 1 vial of V.C.A.T. *Supplement* (1 vial in 500 mL of medium). The plates are inoculated with the strains indicated in the microbiological control table. Incubation conditions: 36 ± 1 °C for 48-72 hours in atmosphere containing 5-10% CO₂. Microbiological control:

Control strains	Growth
<i>Neisseria gonorrhoeae</i> ATCC 43069	Good
<i>Neisseria meningitidis</i> ATCC 13090	Good
<i>Escherichia coli</i> ATCC 25922	Inhibited
<i>Staphylococcus aureus</i> ATCC 25923	Inhibited

PRECAUTIONS

V.C.A.T. *Supplement* contains dangerous substances according to directives 1999/45/CE and 2001/60/CE. V.C.A.T. *Supplement* is a selective supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store V.C.A.T. *Supplement* at 2-8 °C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration. Once reconstituted, the product can be stored for a maximum duration of 30 days at -20°C, shielded from light.

REFERENCES

- Chapin C.K., G.V. Doern (1983). J. Clin. Microbiol. 17: 1163-1165.
- Martin, J.E., Armstrong J.H., Smith, P.B. (1974). Appl. Microbiol. 27: 802-805.
- NCCLS document M22-A2, 1996. Quality Assurance for Commercially Prepared Microbiological Culture Media- Second Ed. Approved Standard.

PRESENTATION

Product	REF	Σ
V.C.A.T. <i>Supplement</i>	81041	10 vials

One vial is sufficient to prepare 500 mL of medium

TABLE OF SYMBOLS

IVD In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
REF Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	LOT Batch code

LIOFILCHEM Bacteriology Products

Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

Tel. +390858930745

Fax +390858930330

Website: www.liofilchem.net

Rev.1 / 21.01.2008

E-Mail: liofilchem@liofilchem.net