CNA (Staf/ Strep) Supplement

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E-Mail: liofilchem@liofilchem.net

Rev.1 / 21.01.2008

Selective supplement for the preparation of COLUMBIA AGAR BASE
for staphylococci, streptococci e pneumococci isolation.

DESCRIPTION

CNA (Staf/ Strep) Supplement is a freeze-dried selective supplement constituted by Nalidixic Acid and Colistin, to use as a supplement of the culture medium COLUMBIA AGAR BASE cod. 610013 or 620013 for staphylococci, streptococci e pneumococci isolation.

KIT CONTENTS

Each kit contains:

- 10 vials of CNA (Staf/ Strep) Supplement freeze-dried.
- 1 instructions sheet

PRINCIPLE OF THE METHOD

CNA (Staf/ Strep) Supplement supports staphylococci, streptococci e pneumococci growth thanks to the selective activity of antibiotics, Nalidixic Acid and Colistin, which inhibit Gram-negative bacteria growth.

COMPOSITION

CNA (Staf/ Strep) Supplement		
	Content / vial	Content / l of medium
Nalidixic Acid	7.5 mg	15.0 mg
Colistin sulphate	5.0 mg	10.0 mg

PROCEDURE FOR USE

1. Aseptically reconstitute the content of one vial of CNA (Staf/ Strep) Supplement with 5 mL of sterile distilled water.
2. Mix to complete dissolution and add to 500 mL of COLUMBIA AGAR BASE cod. 610013 or 620013, autoclaved, cooled to 45-50°C and added with 5% of defibrinated sheep blood.
3. Mix well and dispense in Petri dishes.

TECHNIQUE AND INTERPRETATION OF THE RESULTS

Refer to the technical documentation for COLUMBIA AGAR BASE cod. 610013 or 620013.

QUALITY CONTROL

1. Control of the appearance: freeze-dried yellow in colour.
2. Microbiological control.

Prepare the plates using as base COLUMBIA AGAR BASE cod. 610013 or 620013 added with the content of 1 vial of CNA (Staf/ Strep) Supplement (1 vial in 500 mL of medium). The plates are inoculated with the strains indicated in the microbiological control table. Incubation conditions: 36 ± 1 °C for 18-48 hours in atmosphere containing 5-10% CO₂. Microbiological control:

Control strains	Growth	
<i>Streptococcus pyogenes</i>	ATCC 19615	Good
<i>Streptococcus pneumoniae</i>	ATCC 6303	Good
<i>Staphylococcus aureus</i>	ATCC 25923	Good
<i>Escherichia coli</i>	ATCC 25922	Inhibited

PRECAUTIONS

CNA (Staf/ Strep) Supplement contains dangerous substances according to directives 1999/45/CE and 2001/60/CE. CNA (Staf/ Strep) Supplement is a selective supplement to be used only for *in vitro* diagnostic use. It is intended for use in a professional environment and must be used in the laboratory by properly trained personnel, using approved asepsis and safety methods for handling pathogenic agents.

STORAGE

Store CNA (Staf/ Strep) Supplement at 2-8 °C in its original packaging. Keep away from sources of heat and avoid excessive changes of temperature. Use until the expiry date indicated on the label. Eliminate without using if there are signs of deterioration. Once reconstituted, the product can be stored for a maximum duration of 30 days at -20°C, shielded from light.

REFERENCES

- Ellner, P.D., C.J. Stoessel., E. Drakeford, and F. Vasi (1966). A new culture medium for medical bacteriology. Am. J.Clin. Path. **45**, 502-504.
- Isenberg, H.D. (ed.) (1992). Clinical microbiology procedures handbook, vol. 1 American Society for Microbiology, Washington, DC.

PRESENTATION

Product	REF	Σ
CNA (Staf/ Strep) Supplement	81048	10 vials

One vial is sufficient to prepare 500 mL of medium

TABLE OF SYMBOLS

 In Vitro Diagnostic Medical Device	 Do not reuse	 Manufacturer	 Contains sufficient for <n> tests	 Temperature limitation
 Catalogue number	 Fragile, handle with care	 Use by	 Caution, consult accompanying documents	 Batch code

LIOFILCHEM Bacteriology Products

Via Scozia Zona Ind.le - 64026 Roseto D.A. (TE) - Italy

Tel. +390858930745 Fax +390858930330 Website: www.liofilchem.net E-Mail: liofilchem@liofilchem.net

Rev.1 / 21.01.2008